


Catalyst E Hangers - 70 Degree Extension Tube



Item Number	Part Number	Description	Quantity	Retail
1-17	115918_ _ _	70° Cleat Hanger Assembly, 4-Way-In-Line, Catalyst E	1	\$78.07
1-17	115919_ _ _	70° Cleat Hanger Assembly, 4-Way - Out Left, Catalyst E _ _ _	1	\$78.07
1-17	115920_ _ _	70° Cleat Hanger Assembly, 4-Way - Out Right, Catalyst E _ _ _	1	\$78.07
12 - 17	003225	Assembled 4-Way In Line Hanger Release	1	\$38.61
12 - 17	003225INST	Assembled 4-Way In Line Hanger Release with Instruction	1	\$38.61
12 - 17	003643	Assembled 4-Way Out Hanger Release, Left	1	\$38.61

Item Number	Part Number	Description	Quantity	Retail
12 - 17	003643INST	Assembled 4-Way Out Hanger Release, Left with Instruction	1	\$38.61
12 - 17	003642	Assembled 4-Way Out Hanger Release, Right	1	\$38.61
12 - 17	003642INST	Assembled 4-Way Out Hanger Release, Right with Instruction	1	\$38.61
8 - 11	115845	Cleat Hanger Adjustment Hardware	1	\$11.95
2	000032	Footrest Hanger Pivot Saddle	1	\$5.51
3	100673	M8 x 45 BTN SHCS BLKZC	1	\$1.66
4	100659	M8 Nylock Nut Blk Zn	1	\$1.10
5	100804	Tube End Plug 1"	1	\$1.66
6	100669	M5 x 10 BTN SHCS BLKZC w/Patch	2	\$1.66
7	000073	M5 x 24 Threaded Barrel	1	\$3.97
8	111949	Hanger Cleat Housing	1	\$5.75
9	111950	Hanger Cleat	1	\$3.05
10	111948	Cleat Hanger Sleeve	1	\$1.15
11	114457	M8 x 8 x 1.25 Cup Point Set Screw w/ Patch	1	\$2.00
12	003228	Hanger Latch, 4-Way, Lever	1	\$3.31
13	001374	M5 x 40 BHCS P/THD BLKZC	1	\$2.25
14	100657	M5 Nylock Nut Blk Zn	1	\$0.56
15	000119	Hanger Spring	1	\$3.20
16	003227	Hanger Latch, 4-Way, Pin	1	\$20.96
17	003226	Hanger Latch, 4-Way, Bracket	1	\$5.51

Matte Finish	
Code	Color Name
CGM	Charcoal Gray
EGM	Emerald Green
GGM	Graphite Gray*
IBM	Iridescent Burgundy
INM	Iridescent Navy
AGM	Military Green**
RRM	Retro Red
RBM	Royal Blue
SBM	Shadow Black
SWM	Smooth White

* Discontinued Color

** Formerly called Army Green

Gloss Finish	
Code	Color Name
BLK	Black
BLU	Blue
SKY	Blue Skies
BRO	Bright Orange*
PNK	Bubblegum Pink
BUR	Burgundy*
CBL	Candy Blue
PUR	Candy Purple
RED	Candy Red
COP	Copper
DPR	Deep Purple*
EBL	Electric Blue*
GRN	Forest Green*
FUS	Fuchsia*
GSM	Granny Smith
MTF	Metallic Fuchsia*
MOR	Metallic Orange
ORG	Orange*
BRD	Red
SAP	Sapphire*
SIL	Silver
TEA	Teal
YLW	Yellow
FLM	Graphic Image Transfer