

Positioning Car Seats

Parents' and therapists' **first choice** for assisting in the **safe transport of individuals** with physical disabilities.

Features

- Two sizes available, accommodating users 20-130 lbs and up to 5'6" tall
- Durable, comfortable padded machine-washable cover
- Vacuum-formed ABS plastic shell, chromed steel side tubes
- Meets Fire-Retardant Code 302 for automotive upholstery
- Latex and silicone free

Benefits

- Fits in most cars, vans and school buses
- Crash tested*
- Fits most airplane seats and is compliant for air travel**
- Lightest positioning seat for use in vehicles available, weighing much less than comparable car seats
- Multiple accessories allow positioning seat to grow with the child



#2000 (20 - 102 lbs capacity)



#2500 (65 - 130 lbs capacity)

*#2000 meets or exceeds crash testing requirements in the US (FMVSS #213), Canada (NCSVA #213.3), and Europe (ECE R44.03)

*#2500 meets or exceeds crash testing requirements in the US (FMVSS #213)

**#2000 only

Available colors

- Designer Blue
- Silver Fox Grey
- Sunset Pink

Includes

- Shoulder Harness Pads
- Adjustable Head Pads
- Lateral Positioning Pads
- Easy to Follow Installation Instructions
- 5-Point Safety Harness
- Tether Strap Kit
- Crotch Strap Pad



TheraPedic™ POSITIONING CAR SEATS

MODEL NUMBER	2000	2500
Product Description	Child Positioning Seat	Small Adult Positioning Seat
Product Weight	11 lbs	18 lbs
User Weight Capacity	20 - 102 lbs	65 - 130 lbs ¹
User Height Capacity	Up to 60	54 - 66
(a) Seat-to-Top Height	29	33
(b) Seat Depth	10	12
Seat Depth with Extender ²	12½ or 15	16
(c) Inside Width at Front	15½	19½
Inside Width at Back	13	17½
Outside Width	19	23
Seat to Back Angle	88°	95°

¹ User must be at least 65 lbs to use this product

² Seat depth increases with use of the Seat Depth Extender. See page 6 for details.

Note: Dimensions are in inches unless otherwise noted

Mobility Bases for the TheraPedic Positioning Seat #2000

Columbia Medical offers two lightweight folding frames designed specifically for the #2000 TheraPedic Positioning Seat. Both present alternatives to traditional positioning wheelchairs and can be used when shopping at the mall, going to the playground or while traveling.



TheraPedic Explorer Mobility Base

- Comfortable - five angles of tilt and height-adjustable footrest
- Built to go everywhere - large wheels for a smooth ride over any surface
- Convenient - folds to carry, or store compactly
- Easy to use - seat locks easily onto the base
- Durable - strong yet lightweight frame

MODEL NUMBER	2100
Product Weight	25 lbs
Open	22 1/4W x 49L x 36 1/4H
Folded	22 1/4W x 41L x 11H
Leg Drop Range	5 - 16
Front Wheels	7 diameter
Rear Wheels	10 diameter

Note: Dimensions are in inches unless otherwise noted



Tote Bag

Attaches to the back rail of the #2100 Explorer Mobility Base.

MODEL # 2175

TheraPedic Umbrella Stroller

- Versatile - can be used with or without the TheraPedic Positioning Seat.
- Built to go everywhere - large wheels for a smooth ride over any surface
- Convenient - folds to carry, or store compactly
- Easy to use - seat locks easily onto the base
- Durable - strong yet lightweight frame

Optional Accessories

- Heel Loops #2301
- Shoulder Straps #2303
- Canopy #2304

TheraPedic™ UMBRELLA STROLLER

MODEL NUMBER	2300
Product Weight	22 lbs
Open	25W x 36 1/2L x 47 1/2H
Folded	14W x 47 1/2L x 15H
Leg Drop Range	14 - 18
Seat Width	16
Seat Depth	14
Seat Back Height	25
Wheel Diameter	7

Note: Dimensions are in inches unless otherwise noted



Available colors

-  Electric Blue
-  Purple
-  Turquoise

Positioning Seat Accessories



Custom Inserts

A unique way to provide custom molded seating for our TheraPedic Positioning Seats. Easily mold and capture a shape, and then deliver the finished product to the user in one hour, without a simulator.


MODEL # 2030 - for use with #2000 **2530** - for use with #2500



Buckle Guard

Keeps "Little Houdini" from releasing the shoulder harness. For use with all TheraPedic Positioning seats.

MODEL # 2075



Block Style Positioning Pads

Additional pads for special positioning needs. Available in four packs. For use with all TheraPedic Positioning seats.

MODEL # 2008

Block style, two pads 2½" thick & two pads 1½" thick



Seat Depth Extender

Increases the seat depth.


For Child Positioning Seat #2000

MODEL # 2025 - Increases seat depth by 2½"

2026* - Increases seat depth by 5"

For Small Adult Positioning Seat #2500

MODEL # 2525* - Increases seat depth by 4"



Closed

Open

Swing Away Abductors

Improves hip alignment and stabilizes seating posture.

Dimensions: 3½"W x 4"H x 3"D

MODEL # 2015* - For use with #2000

2018* - For use with #2000 with 2½" extender #2025

2019* - For use with #2000 with 5" extender #2026

2515* - For use with #2500

2519* - For use with #2500 with 4" extender #2525



Lateral Positioning Pads

Four additional pads to improve lateral positioning.

MODEL # 2010 - For use with #2000

2510 - For use with #2500



Extensor Thrust Wedge

Changes seat to back angle. Increases height by 2" or 3".

For Child Positioning Seat #2000

MODEL # 2050 - 2" front height to 1" rear

2055 - 3" front height to 1" rear

For Small Adults Positioning Seat #2500

MODEL # 2550 - 2" front height to 1" rear

2555 - 3" front height to 1" rear

* When ordering this item for an existing TheraPedic seat, new Seat Cover Upholstery (#2002 for the #2000 or #2502 for the #2500) must be purchased.