


Catalyst E Handrims - Superlight



Superlight Handrim with Micro Grip



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel. <p>

Not available with Heavy Duty Option<p>

Tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
2b, 3	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy and MAXX Spoke wheels</i>	1	\$16.55
2a	502170	Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Spoke</i>	1	\$16.55
1a	101161	Handrim Superlight Aluminum Anodized 24", 6 Tab	1	\$98.18
1b	505585	Handrim Superlight Aluminum Anodized 24" w/ Micro Grip	1	\$137.50
1c	101160	Handrim Superlight Aluminum Anodized 25", 6 Tab	1	\$98.18
1d	505586	Handrim Superlight Aluminum Anodized 25" w/ Micro Grip	1	\$138.68

Item Number	Part Number	Description	Quantity	Retail
2a	100536	M5 x 12 BTN SHCS ZC <i>Used with Superlight Spoke</i>	6	\$2.21
2b	100669	M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinergy and MAXX Spoke wheels</i>	6	\$1.66
3	100657	M5 Nylock Nut Blk Zn <i>Used with Spinergy and MAXX Spoke wheels</i>	6	\$0.56

Handrim Hardware Chart

Wheel		Handrim Connection				Handrim								
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim		
18" Ki Spoke ¹	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" Ki Spoke ¹	200530	6				Natural Fit: 100835	200536		200542					
22" Ki Spoke ¹	200531						100560		100576	100569	200538	200201		
24" Ki Spoke ¹	200532				200349			100577	200547	100793	200202			
25" Ki Spoke ¹	200533				200350		101870	200548	200539	200540				
26" Ki Spoke ¹	200534				200351		100578	200549	100907	101454				
24" Superlight	101159		6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460		
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963						
22" 5-Spoke X Core	101962					101898		101964						
24" 5-Spoke X Core	100960					100975		100976	200546					
24" 5-Spoke X Core	100960				101893					100768	100769			
20" Spinerger Spox	See Spinerger Spox Page	6			Screw: 100669 Nut: 100657	103125		103179						
22" Spinerger Spox						100827		100808		100889	100888			
24" Spinerger Spox						100766	101161	100615		100830	100828			
25" Spinerger Spox						100767	101160	100765		101464	101460			
26" Spinerger Spox						101477		101148		200200	100950			
22" Spinerger LX	See Spinerger LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888			
24" Spinerger LX						100766	101161	100615		100830	100828	See Spinerger LX Page		
25" Spinerger LX						100767	101160	100765		101464	101460			
26" Spinerger LX						101477		101148		200200	100950			
18" MAXX Spoke	109285	3	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" MAXX Spoke	108244	6				Natural Fit: 100835	200536		200542					
22" MAXX Spoke	105135						100560		100576	100569	200538	200201		
24" MAXX Spoke	105136				200349		101161	100577	200547	100793	200202			
25" MAXX Spoke	107436				200350		101870	200548	200539	200540				
26" MAXX Spoke	107437				200351		100578	200549	100907	101454				

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.

¹ Discontinued