


Axiom AS / ASP Solid Seat Insert



When Seat Rail Modification is present, you'll need a solid seat insert 1" shorter in width. When Notch for Back Post Modification is present, you'll need a solid seat insert that is 1" shorter in depth. When both the Seat Rail Modification and Notch for Back Post Modification is present, you'll need a solid seat insert that is 1" shorter in width and 1" shorter in depth.

Item Number	Part Number	Description	Quantity	Retail
1a, 2	XI1414A	Axiom Solid Seat Insert 14 x 14 with AB Tape	1	\$62.00
1b, 2	XI1415A	Axiom Solid Seat Insert 14 x 15 with AB Tape	1	\$62.00
1c, 2	XI1416A	Axiom Solid Seat Insert 14 x 16 with AB Tape	1	\$62.00
1d, 2	XI1417A	Axiom Solid Seat Insert 14 x 17 with AB Tape	1	\$62.00
1e, 2	XI1418A	Axiom Solid Seat Insert 14 x 18 with AB Tape	1	\$62.00
1f, 2	XI1419A	Axiom Solid Seat Insert 14 x 19 with AB Tape	1	\$62.00
1g, 2	XI1420A	Axiom Solid Seat Insert 14 x 20 with AB Tape	1	\$62.00
1h, 2	XI1514A	Axiom Solid Seat Insert 15 x 14 with AB Tape	1	\$62.00

Item Number	Part Number	Description	Quantity	Retail
1i, 2	XI1515A	Axiom Solid Seat Insert 15 x 15 with AB Tape	1	\$62.00
1j, 2	XI1516A	Axiom Solid Seat Insert 15 x 16 with AB Tape	1	\$62.00
1k, 2	XI1517A	Axiom Solid Seat Insert 15 x 17 with AB Tape	1	\$62.00
1l, 2	XI1518A	Axiom Solid Seat Insert 15 x 18 with AB Tape	1	\$62.00
1m, 2	XI1519A	Axiom Solid Seat Insert 15 x 19 with AB Tape	1	\$62.00
1n, 2	XI1520A	Axiom Solid Seat Insert 15 x 20 with AB Tape	1	\$62.00
1o, 2	XI1521A	Axiom Solid Seat Insert 15 x 21 with AB Tape	1	\$62.00
1p, 2	XI1522A	Axiom Solid Seat Insert 15 x 22 with AB Tape	1	\$62.00
1q, 2	XI1614A	Axiom Solid Seat Insert 16 x 14 with AB Tape	1	\$62.00
1r, 2	XI1615A	Axiom Solid Seat Insert 16 x 15 with AB Tape	1	\$62.00
1s, 2	XI1616A	Axiom Solid Seat Insert 16 x 16 with AB Tape	1	\$62.00
1t, 2	XI1617A	Axiom Solid Seat Insert 16 x 17 with AB Tape	1	\$62.00
1u, 2	XI1618A	Axiom Solid Seat Insert 16 x 18 with AB Tape	1	\$62.00
1v, 2	XI1619A	Axiom Solid Seat Insert 16 x 19 with AB Tape	1	\$62.00
1w, 2	XI1620A	Axiom Solid Seat Insert 16 x 20 with AB Tape	1	\$62.00
1x, 2	XI1621A	Axiom Solid Seat Insert 16 x 21 with AB Tape	1	\$62.00
1y, 2	XI1622A	Axiom Solid Seat Insert 16 x 22 with AB Tape	1	\$62.00
1z, 2	XI1715A	Axiom Solid Seat Insert 17 x 15 with AB Tape	1	\$62.00
1aa, 2	XI1716A	Axiom Solid Seat Insert 17 x 16 with AB Tape	1	\$62.00
1ab, 2	XI1717A	Axiom Solid Seat Insert 17 x 17 with AB Tape	1	\$62.00
1ac, 2	XI1718A	Axiom Solid Seat Insert 17 x 18 with AB Tape	1	\$62.00
1ad, 2	XI1719A	Axiom Solid Seat Insert 17 x 19 with AB Tape	1	\$62.00
1ae, 2	XI1720A	Axiom Solid Seat Insert 17 x 20 with AB Tape	1	\$62.00
1af, 2	XI1721A	Axiom Solid Seat Insert 17 x 21 with AB Tape	1	\$62.00
1ag, 2	XI1722A	Axiom Solid Seat Insert 17 x 22 with AB Tape	1	\$62.00
1ah, 2	XI1816A	Axiom Solid Seat Insert 18 x 16 with AB Tape	1	\$62.00
1ai, 2	XI1817A	Axiom Solid Seat Insert 18 x 17 with AB Tape	1	\$62.00
1aj, 2	XI1818A	Axiom Solid Seat Insert 18 x 18 with AB Tape	1	\$62.00

Item Number	Part Number	Description	Quantity	Retail
1ak, 2	XI1819A	Axiom Solid Seat Insert 18 x 19 with AB Tape	1	\$62.00
1al, 2	XI1820A	Axiom Solid Seat Insert 18 x 20 with AB Tape	1	\$62.00
1am, 2	XI1821A	Axiom Solid Seat Insert 18 x 21 with AB Tape	1	\$62.00
1an, 2	XI1822A	Axiom Solid Seat Insert 18 x 22 with AB Tape	1	\$62.00
1ao, 2	XI1916A	Axiom Solid Seat Insert 19 x 16 with AB Tape	1	\$62.00
1ap, 2	XI1917A	Axiom Solid Seat Insert 19 x 17 with AB Tape	1	\$62.00
1aq, 2	XI1918A	Axiom Solid Seat Insert 19 x 18 with AB Tape	1	\$62.00
1ar, 2	XI1919A	Axiom Solid Seat Insert 19 x 19 with AB Tape	1	\$62.00
1as, 2	XI1920A	Axiom Solid Seat Insert 19 x 20 with AB Tape	1	\$62.00
1at, 2	XI1921A	Axiom Solid Seat Insert 19 x 21 with AB Tape	1	\$62.00
1au, 2	XI1922A	Axiom Solid Seat Insert 19 x 22 with AB Tape	1	\$62.00
1av, 2	XI2016A	Axiom Solid Seat Insert 20 x 16 with AB Tape	1	\$62.00
1aw, 2	XI2017A	Axiom Solid Seat Insert 20 x 17 with AB Tape	1	\$62.00
1ax, 2	XI2018A	Axiom Solid Seat Insert 20 x 18 with AB Tape	1	\$62.00
1ay, 2	XI2019A	Axiom Solid Seat Insert 20 x 19 with AB Tape	1	\$62.00
1az, 2	XI2020A	Axiom Solid Seat Insert 20 x 20 with AB Tape	1	\$62.00
1ba, 2	XI2021A	Axiom Solid Seat Insert 20 x 21 with AB Tape	1	\$62.00
1bb, 2	XI2022A	Axiom Solid Seat Insert 20 x 22 with AB Tape	1	\$62.00
1bc, 2	XI2118A	Axiom Solid Seat Insert 21 x 18 with AB Tape	1	\$62.00
1bd, 2	XI2119A	Axiom Solid Seat Insert 21 x 19 with AB Tape	1	\$62.00
1be, 2	XI2120A	Axiom Solid Seat Insert 21 x 20 with AB Tape	1	\$62.00
1bf, 2	XI2121A	Axiom Solid Seat Insert 21 x 21 with AB Tape	1	\$62.00
1bg, 2	XI2122A	Axiom Solid Seat Insert 21 x 22 with AB Tape	1	\$62.00
1bh, 2	XI2216A	Axiom Solid Seat Insert 22 x 16 with AB Tape	1	\$62.00
1bi, 2	XI2217A	Axiom Solid Seat Insert 22 x 17 with AB Tape	1	\$62.00
1bj, 2	XI2218A	Axiom Solid Seat Insert 22 x 18 with AB Tape	1	\$62.00
1bk, 2	XI2219A	Axiom Solid Seat Insert 22 x 19 with AB Tape	1	\$62.00
1bl, 2	XI2220A	Axiom Solid Seat Insert 22 x 20 with AB Tape	1	\$62.00

Item Number	Part Number	Description	Quantity	Retail
1bm, 2	XI2221A	Axiom Solid Seat Insert 22 x 21 with AB Tape	1	\$62.00
1bn, 2	XI2222A	Axiom Solid Seat Insert 22 x 22 with AB Tape	1	\$62.00