


Catalyst E Anti-tips



Item Number	Part Number	Description	Quantity
1-2, 3a, 4, 5a, 6-8	115303	Anti Tip Assembly Short, Catalyst E <i>Used with 11.5" - 12.5" RSH</i>	1
1-2, 3b, 4, 5a, 6-8	115304	Anti Tip Assembly, Medium, Catalyst E <i>Used with 13" - 14" RSH</i>	1
1-2, 3c, 4, 5a, 6-8	115922	Anti Tip Assembly, Long, Catalyst E <i>Used with 14.5" - 18" RSH</i>	1
1-2, 3e, 4, 5b, 6-8	115923	Anti Tip Assembly Extra Long, Catalyst E <i>Used with 18.5" - 20" RSH</i>	1
1	111943	Anti-Tip Sleeve, Catalyst E	1
2	100533	Detent Button	1
3a	112552	Anti-Tip Tube Upper Short, CAT E <i>Used with 11.5" - 12.5" RSH</i>	1
3b	112553	Anti-Tip Tube Upper Medium, CAT E <i>Used with 13" - 14" RSH</i>	1
4	100534	Detent Button Lower	1

Item Number	Part Number	Description	Quantity
5a	000992	Anti-Tip Extension Short	1
5b	000076	Anti-Tip Extension Long, Anodized	1
6	100676	M6 x 55 CLS10.9 DIN931 HHCS ZC	1
8	100658	M6 Nylock Nut Blk Zn	1