

Catalyst 5Ti Tires - Shox



<p>

Please note tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
1a	102717	Tire 20" x 1" Shox Spoke/Mag <i>Used on 20" Spoke and Mag Wheels</i>	1	\$73.54
1b	102718	Tire 22" x 1" Shox Spoke/Mag <i>Used on 22" Spoke and Mag Wheels</i>	1	\$73.54
1c	100640	Tire 24" x 1" Shox Spoke/Mag <i>Used on 24" Spoke and Mag Wheels</i>	1	\$73.54
1d	100894	Tire 24" x 1" Shox Spinergy/Topolino/Superlight <i>Used on 24" Spinergy and Superlight Spoke Wheels</i>	1	\$73.54
1e	101453	Tire 25" x 1" Shox Spoke <i>Used on 25" Spoke Wheels</i>	1	\$73.54

Item Number	Part Number	Description	Quantity	Retail
1f	101146	Tire 25" x 1" Shox Spinergy/Topolino/Superlight <i>Used on 25" Spinergy and Superlight Spoke Wheels</i>	1	\$73.54
1g	100884	Tire 26" x 1" Shox Spoke <i>Used on 26" Spoke Wheels</i>	1	\$73.54
1h	100895	Tire 26" x 1" Shox Spinergy <i>Used on 26" Spinergy Wheels</i>	1	\$73.54

Wheel & Tire Outside Diameter Measurement Reference

Size	Tire Description	Outside Diameter
18"	Street Pneumatic 1-3/8"	18.9
18"	Full Profile Poly 1-3/8"	18.5
18"	Pneumatic w/Airless Insert 1-3/8"	18.8
18"	Low Profile Poly 1"	17.7
20"	Street Pneumatic 1-3/8"	20.7
20"	Full Profile Poly 1-3/8"	20.2
20"	Pneumatic w/Airless Insert 1-3/8"	20.5
20"	Low Profile Poly 1"	19.3
20"	Sentinel Solid 1"	19.6
20"	High Pressure 1"	20.3
20"	SHOX	19.6
22"	Street Pneumatic 1-3/8"	22.5
22"	Full Profile Poly 1-3/8"	22.3
22"	Pneumatic w/Airless Insert 1-3/8"	22.4
22"	Low Profile Poly 1"	21.4
22"	Sentinel Solid 1"	21.9
22"	High Pressure 1"	22.3
22"	SHOX	22
22"	Mountain 2"	24.1
24"	Street Pneumatic 1-3/8"	24
24"	Full Profile Poly 1-3/8"	23.6
24"	Pneumatic w/Airless Insert 1-3/8"	23.9
24"	Schwalbe HS-127 1-3/8"	24
24"	Euro Trek Knobby 1-3/8"	24.1
24"	Low Profile Poly 1"	22.8
24"	Sentinel Solid 1"	23.3
24"	High Pressure 1"	23.6
24"	Schwalbe One 1"	23.8
24"	Schwalbe Marathon Plus EVO 1"	24
24"	SHOX	23.3
24"	Schwalbe Downtown 1"	23.6
24"	Mountain 2"	25.6
25"	Sentinel Solid 1"	24
25"	High Pressure 1"	24.3
25"	Schwalbe One 1"	24.3
25"	Schwalbe Marathon Plus EVO 1	24.4
25"	SHOX	24.1
25"	Mountain 2"	26
26"	Street Pneumatic 1-3/8"	26
26"	Full Profile Poly 1-3/8"	25.8
26"	Pneumatic w/Airless Insert 1-3/8"	25.8
26"	Sentinel Solid 1"	25.1
26"	High Pressure 1"	25.5
26"	Schwalbe Marathon Plus EVO 1"	25.9
26"	SHOX	25.3

Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part Number	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke	200532	100635	100622	100633	100635 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke	200534	100636			100636 / 100505	102978	100884	102960		100865			
24" Superlight	101159	100635					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
12" Mag ¹	N/A	001424	001426		001425								
16" Mag ¹	N/A	001159	001161	003453	001160								
20" 5-Spoke X Core (Black)	101961	100911	100917	102085	100911 / 100913		102717						
22" 5-Spoke X Core (Black)	101962	100634	100621	100632	100634 / 100503		102718						
24" 5-Spoke X Core (Black)	100960	100635	100622	100633	100635 / 100504		100640		102982			102920	
20" Spinergy Spox	See Spinergy Spox Page	100911	100917	102085	100911 / 100913	102974		102957					
22" Spinergy Spox		100634	100621	100632	100634 / 100503	103165		102961					101475
24" Spinergy Spox		100635	100622	100633	100635 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy Spox						102977	101146	102959		101357	103177		100166
26" Spinergy Spox		100636			100636 / 100505	102978	100895	102960		100865			
22" Spinergy LX	See Spinergy LX Page	100634	100621	100632	100634 / 100503			102961					101475
24" Spinergy LX		100635	100622	100633	100635 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy LX						102977	101146	102959		101357	103177		100166
26" Spinergy LX		100636			100636 / 100505	102978	100895	102960		100865			
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	100635	100622	100633	100635 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			

¹ 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly