


# Liberty Tilt Mechanism - Growth

CABLE LENGTHS			
BACK POST TYPE	SEAT DEPTH RANGE		
	SHORT (14-16D)	MEDIUM (16-18D)	LONG (18-20D)
HEIGHT ADJ SHORT	003288	003289	003290
HEIGHT ADJ MEDIUM	003289	003290	003291
HEIGHT ADJ TALL	003290	003291	003292
STROLLER HANDLE SHORT	003289	003290	003291
STROLLER HANDLE TALL	003291	003292	003293



*Did your depth change? If yes, evaluate if a different size cable is necessary*

Item Number	Part Number	Description	Quantity
1-3, 4a	506801	Cable/Trigger Kit with Short Cable, Liberty	1
1-3, 4b	506802	Cable/Trigger Kit with Medium Cable, Liberty	1
1-3, 4c	506803	Cable/Trigger Kit with Long Cable, Liberty	1
1-3, 4d	506804	Cable/Trigger Kit with XL Cable, Liberty	1
1-3, 4e	506805	Cable/Trigger Kit with XXL Cable, Liberty	1
1-3, 4f	506806	Cable/Trigger Kit with XXXL Cable, Liberty	1
4a	003288	Cable, Tilt, Short	2

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
4b	003289	Cable, Tilt, Medium	2
4c	003290	Cable, Tilt, Long	2
4d	003291	Cable, Tilt, XL	2
4e	003292	Cable, Tilt, XXL	2
4f	003293	Cable, Tilt, XXXL	2