


Focus CR Base Frame - Growth



Growth Strut	Length
Short	6.6'
Medium	8.6'
Long	10.6'
Extra Long	12.6'

Has your width changed? Will your current base struts work? If not, please select appropriately. To grow the chair, a quantity of two base struts will be needed. One additional base strut is needed if the chair configuration includes an attendant wheel lock. If quoting a Left Frame for Australia or New Zealand, 122921 is required when transit is not present.

Item Number	Part Number	Description	Quantity
1b	001011	Growth Strut, Base, Small Tilt <i>Used on Small Frame Width Range</i>	1
1b	001012	Growth Strut, Base, Medium Tilt <i>Used on Medium Frame Width Range</i>	1
1c	001013	Growth Strut, Base, Large Tilt <i>Used on Large Frame Width Range</i>	1
1d	112005	Growth Strut, Base, X-Large Tilt <i>Used on X-Large Frame Width Range. This size can only be used with the Heavy Duty Base Frame.</i>	1
	122921	Not AS/NZS 3696.19 Compliant On-Chair Label with Instruction	1