LOW PROFILE ANTI-TIP ORDERING GUIDE QX STANDARD/HEMI FRAME							
AXLE HOLE POSITION	REAR WHEEL SIZE						
	20"	22"	24"	25"	26"		
5	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
4	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM		
3	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE		
2	MEDIUM	MEDIUM	LARGE	LARGE	LARGE		
1	MEDIUM	LARGE	LARGE	LARGE	LARGE		

LOW PROFILE ANTI-TIP ORDERING GUIDE QXI STANDARD FRAME										
	REAR WHEEL SIZE									
AXLE HOLE POSITION	20"	22"			24"			25"	26"	
	ALL TIRES	HP CLINCHER BLACK	FULL POLY	ALL OTHER TIRES	SCHWALBE MARATHON PLUS	HP CLINCHER BLACK	ALL OTHER TIRES	ALL TIRES	ALL TIRES	
13	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
12	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
11	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
10	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
9	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
8	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	
7	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE	
6	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE	LARGE	
5	MEDIUM	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE	LARGE	LARGE	
4	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE	LARGE	LARGE	LARGE	
3	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	
2	MEDIUM	MEDIUM	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	
1	MEDIUM	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	LARGE	

LOW PROFILE ANTI-TIP ORDERING GUIDE QXI HEMI FRAME							
AXLE HOLE POSITION	REAR WHEEL SIZE						
	20"	22"	24"	25"	26"		
13	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
12	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
11	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
10	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
9	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
8	SMALL	SMALL	MEDIUM	MEDIUM	MEDIUM		
7	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM		
6	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM		
5	SMALL	MEDIUM	MEDIUM	MEDIUM	MEDIUM		
4	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE		
3	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE		
2	MEDIUM	MEDIUM	MEDIUM	MEDIUM	LARGE		
1	MEDIUM	MEDIUM	MEDIUM	LARGE	LARGE		