


Removable Pediatric Side Guard



If retrofitting, side guard receivers may be required

Item Number	Part Number	Description	Quantity
1-4	503832	Pediatric Side Guard Assembly with Ki Logo Right Side	1
1-4	503833	Pediatric Side Guard Assembly with Ki Logo Left Side	1
1-4	502180	Pediatric Side Guard Assembly with Ki Logo PAIR	1

Item Number	Part Number	Description	Quantity
1	003101	Pediatric Side Guard Panel w/Ki Logo	1
2	000140	Side Guard Post	1
3	101184	M6 x 10 BTN SHCS BLKZC w/Patch	3
4	001342	Side Guard Post Stop	1