


(Discontinued) Spark Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

Ki Mobility Mag Wheel/ Tire/ Handrim Kit Item Number Reference

| Rear Tire Type | Handrim | Wheel | | |
|------------------------------------|--------------------------------------|------------------------|------------------------|------------------------|
| | | 20" 5-Spoke X Core Mag | 22" 5-Spoke X Core Mag | 24" 5-Spoke X Core Mag |
| 1 3/8" Street Pneumatic | Aluminum Anodized | 502090 | 502100 | 502110 |
| | Plastic Coated | 502091 | 502101 | 502111 |
| | Projection | | | 502112 |
| | Natural Fit Standard Omit Grip | | | 502113 |
| | Natural Fit Standard w/ Regular Grip | | | 502114 |
| | Natural Fit Standard w/ Super Grip | | | 502115 |
| | Natural Fit LT Omit Grip | | | 502116 |
| | Natural Fit LT w/ Regular Grip | | | 502117 |
| | Natural Fit LT w/ Super Grip | | | 502118 |
| 1 3/8" Full Profile Poly | Aluminum Anodized | 502094 | 502104 | 502127 |
| | Plastic Coated | 502095 | 502105 | 502128 |
| | Projection | | | 502129 |
| | Natural Fit Standard Omit Grip | | | 502130 |
| | Natural Fit Standard w/ Regular Grip | | | 502131 |
| | Natural Fit Standard w/ Super Grip | | | 502132 |
| | Natural Fit LT Omit Grip | | | 502133 |
| | Natural Fit LT w/ Regular Grip | | | 502134 |
| | Natural Fit LT w/ Super Grip | | | 502135 |
| 1 3/8" Pneumatic w/ Airless Insert | Aluminum Anodized | 502096 | 502106 | 502136 |
| | Plastic Coated | 502097 | 502107 | 502137 |
| | Projection | | | 502138 |
| | Natural Fit Standard Omit Grip | | | 502139 |
| | Natural Fit Standard w/ Regular Grip | | | 502140 |
| | Natural Fit Standard w/ Super Grip | | | 502141 |
| | Natural Fit LT Omit Grip | | | 502142 |
| | Natural Fit LT w/ Regular Grip | | | 502143 |
| | Natural Fit LT w/ Super Grip | | | 502144 |

Ki Mobility Mag Wheel/ Tire/ Handrim Kit Item Number Reference

| Rear Tire Type | Handrim | Wheel | | |
|-------------------------|--------------------------------------|------------------------|------------------------|------------------------|
| | | 20" 5-Spoke X Core Mag | 22" 5-Spoke X Core Mag | 24" 5-Spoke X Core Mag |
| 1 3/8" Euro Trek Knobby | Aluminum Anodized | | | 502153 |
| | Plastic Coated | | | 502154 |
| | Projection | | | 502155 |
| | Natural Fit Standard Omit Grip | | | 502156 |
| | Natural Fit Standard w/ Regular Grip | | | 502157 |
| | Natural Fit Standard w/ Super Grip | | | 502158 |
| | Natural Fit LT Omit Grip | | | 502159 |
| | Natural Fit LT w/ Regular Grip | | | 502160 |
| | Natural Fit LT w/ Super Grip | | | 502161 |
| 1" Low Profile Poly | Aluminum Anodized | 502092 | 502102 | 502119 |
| | Plastic Coated | 502093 | 502103 | 502120 |
| | Projection | | | |
| | Natural Fit Standard Omit Grip | | | 502121 |
| | Natural Fit Standard w/ Regular Grip | | | 502122 |
| | Natural Fit Standard w/ Super Grip | | | 502123 |
| | Natural Fit LT Omit Grip | | | 502124 |
| | Natural Fit LT w/ Regular Grip | | | 502125 |
| | Natural Fit LT w/ Super Grip | | | 502126 |
| 1" Sentinel Solid Tire | Aluminum Anodized | | | |
| | Plastic Coated | | | |
| | Projection | | | |
| | Natural Fit Standard Omit Grip | | | |
| | Natural Fit Standard w/ Regular Grip | | | |
| | Natural Fit Standard w/ Super Grip | | | |
| | Natural Fit LT Omit Grip | | | |
| | Natural Fit LT w/ Regular Grip | | | |
| | Natural Fit LT w/ Super Grip | | | |