


Tsunami AL Seat Upholstery - Growth



SCREW (ITEM 3) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	8
16" - 18"	10
19" - 21"	12
22"	14

*Seat upholstery width measures approximately 1" less than the chair size <p>
 Seat upholstery depth measures true to the chair size*

Item Number	Part Number	Description	Quantity	Retail
1a	101000	Seat Upholstery Rigid/LFT 1414	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1b	101005	Seat Upholstery Rigid 1415	1	\$92.70
1c	101012	Seat Upholstery Rigid/LFT 1416	1	\$92.70
1d	101019	Seat Upholstery Rigid 1417	1	\$92.70
1e	101026	Seat Upholstery Rigid/LFT 1418	1	\$92.70
1f	101033	Seat Upholstery Rigid 1419	1	\$92.70
1g	101074	Seat Upholstery Rigid/LFT 1420	1	\$92.70
1h	101001	Seat Upholstery Rigid/LFT 1514	1	\$92.70
1i	101006	Seat Upholstery Rigid 1515	1	\$92.70
1j	101013	Seat Upholstery Rigid/LFT 1516	1	\$92.70
1k	101020	Seat Upholstery Rigid 1517	1	\$92.70
1l	101027	Seat Upholstery Rigid/LFT 1518	1	\$92.70
1m	101034	Seat Upholstery Rigid 1519	1	\$92.70
1n	101040	Seat Upholstery Rigid/LFT 1520	1	\$92.70
1o	101002	Seat Upholstery Rigid/LFT 1614	1	\$92.70
1p	101007	Seat Upholstery Rigid 1615	1	\$92.70
1q	101014	Seat Upholstery Rigid/LFT 1616	1	\$92.70
1r	101021	Seat Upholstery Rigid 1617	1	\$92.70
1s	101028	Seat Upholstery Rigid/LFT 1618	1	\$92.70
1t	101035	Seat Upholstery Rigid 1619	1	\$92.70
1u	101041	Seat Upholstery Rigid/LFT 1620	1	\$92.70
1v	101542	Seat Upholstery Rigid 1621	1	\$92.70
1w	101543	Seat Upholstery Rigid 1622	1	\$92.70
1x	101003	Seat Upholstery Rigid/LFT 1714	1	\$92.70
1y	101008	Seat Upholstery Rigid 1715	1	\$92.70
1z	101015	Seat Upholstery Rigid/LFT 1716	1	\$92.70
1aa	101022	Seat Upholstery Rigid 1717	1	\$92.70
1ab	101029	Seat Upholstery Rigid/LFT 1718	1	\$92.70
1ac	101036	Seat Upholstery Rigid 1719	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1ad	101042	Seat Upholstery Rigid/LFT 1720	1	\$92.70
1ae	101544	Seat Upholstery Rigid 1721	1	\$92.70
1af	101545	Seat Upholstery Rigid 1722	1	\$92.70
1ag	101004	Seat Upholstery Rigid/LFT 1814	1	\$92.70
1ah	101009	Seat Upholstery Rigid 1815	1	\$92.70
1ai	101016	Seat Upholstery Rigid/LFT 1816	1	\$92.70
1aj	101023	Seat Upholstery Rigid 1817	1	\$92.70
1ak	101030	Seat Upholstery Rigid/LFT 1818	1	\$92.70
1al	101037	Seat Upholstery Rigid 1819	1	\$92.70
1am	101043	Seat Upholstery Rigid/LFT 1820	1	\$92.70
1an	101546	Seat Upholstery Rigid 1821	1	\$92.70
1ao	101547	Seat Upholstery Rigid 1822	1	\$92.70
1ap	101010	Seat Upholstery Rigid 1915	1	\$92.70
1aq	101017	Seat Upholstery Rigid/LFT 1916	1	\$92.70
1ar	101024	Seat Upholstery Rigid 1917	1	\$92.70
1as	101031	Seat Upholstery Rigid/LFT 1918	1	\$92.70
1at	101038	Seat Upholstery Rigid 1919	1	\$92.70
1au	101044	Seat Upholstery Rigid/LFT 1920	1	\$92.70
1av	101548	Seat Upholstery Rigid 1921	1	\$92.70
1aw	101549	Seat Upholstery Rigid 1922	1	\$92.70
1ax	101011	Seat Upholstery Rigid 2015	1	\$92.70
1ay	101018	Seat Upholstery Rigid/LFT 2016	1	\$92.70
1az	101025	Seat Upholstery Rigid 2017	1	\$92.70
1ba	101032	Seat Upholstery Rigid/LFT 2018	1	\$92.70
1bb	101039	Seat Upholstery Rigid 2019	1	\$92.70
1bc	101045	Seat Upholstery Rigid/LFT 2020	1	\$92.70
1bd	101550	Seat Upholstery Rigid 2021	1	\$92.70
1be	101551	Seat Upholstery Rigid 2022	1	\$92.70

Item Number	Part Number	Description	Quantity	Retail
1bf	101528	Seat Upholstery Rigid/LFT 2116	1	\$92.70
1bg	101529	Seat Upholstery Rigid 2117	1	\$92.70
1bh	101530	Seat Upholstery Rigid/LFT 2118	1	\$92.70
1bi	101531	Seat Upholstery Rigid 2119	1	\$92.70
1bj	101532	Seat Upholstery Rigid/LFT 2120	1	\$92.70
1bk	101533	Seat Upholstery Rigid 2121	1	\$92.70
1bl	101534	Seat Upholstery Rigid 2122	1	\$92.70
1bm	101535	Seat Upholstery Rigid/LFT 2216	1	\$92.70
1bn	101536	Seat Upholstery Rigid 2217	1	\$92.70
1bo	101537	Seat Upholstery Rigid/LFT 2218	1	\$92.70
1bp	101538	Seat Upholstery Rigid 2219	1	\$92.70
1bq	101539	Seat Upholstery Rigid/LFT 2220	1	\$92.70
1br	101540	Seat Upholstery Rigid 2221	1	\$92.70
1bs	101541	Seat Upholstery Rigid 2222	1	\$92.70
2a	000241	Seat Rail 14" Deep Rigid/LFT	2	\$15.37
2b	000242	Seat Rail 15" Deep RIGID	2	\$15.37
2c	000243	Seat Rail 16" Deep Rigid/LFT	2	\$15.37
2d	000244	Seat Rail 17" Deep RIGID	2	\$15.37
2e	000245	Seat Rail 18" Deep Rigid/LFT	2	\$15.37
2f	000246	Seat Rail 19" Deep RIGID	2	\$15.37
2g	000247	Seat Rail 20" Deep Rigid/LFT	2	\$15.37
2h	000552	Seat Rail 21" Deep RIGID	2	\$12.35
2i	000553	Seat Rail 22" Deep RIGID	2	\$13.14
3	100734	M5 x 16 PHMS ZC	14	\$1.55