


# Rogue XP Wheels - SpinerGY LX



<p>

Recommended to order the SpinerGY LX Replacement Spoke Tool Kit (110119) when replacing spokes. <p>

Replacement SpinerGY Spokes and SpinerGY Wheels with Flex Rims may have an extended lead time. <p>

If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. <p>

Please note tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
1a	100886	22" x 1/2" SpinerGY LX Black Spoke	1	\$819.47

Item Number	Part Number	Description	Quantity	Retail
1b	100939	22" x 1/2" Spinergy LX Red Spoke	1	\$830.92
1c	100940	22" x 1/2" Spinergy LX White Spoke	1	\$830.92
1d	100938	22" x 1/2" Spinergy LX Blue Spoke	1	\$830.92
1e	100935	22" x 1/2" Spinergy LX Yellow Spoke	1	\$830.92
1f	101096	22" x 1/2" Spinergy LX Pink Spoke	1	\$830.92
1g	101735	22" x 1/2" Spinergy LX Green Spoke	1	\$830.92
1h	101736	22" x 1/2" Spinergy LX Orange Spoke	1	\$830.92
1i	107394	22" x 1/2" Spinergy LX Purple Spoke	1	\$830.92
1j	100942	24" x 1/2" Spinergy LX Black Spoke	1	\$819.47
1k	100941	24" x 1/2" Spinergy LX Red Spoke	1	\$830.92
1l	100944	24" x 1/2" Spinergy LX White Spoke	1	\$830.92
1m	100943	24" x 1/2" Spinergy LX Blue Spoke	1	\$830.92
1n	100919	24" x 1/2" Spinergy LX Yellow Spoke	1	\$830.92
1o	101738	24" x 1/2" Spinergy LX Pink Spoke	1	\$830.92
1p	101737	24" x 1/2" Spinergy LX Green Spoke	1	\$830.92
1q	200212	24" x 1/2" Spinergy LX Orange Spoke	1	\$830.92
1r	107395	24" x 1/2" Spinergy LX Purple Spoke	1	\$830.92
1s	100994	25" x 1/2" Spinergy LX Black Spoke	1	\$819.47
1t	100996	25" x 1/2" Spinergy LX Red Spoke	1	\$830.92
1u	100998	25" x 1/2" Spinergy LX White Spoke	1	\$830.92
1v	100995	25" x 1/2" Spinergy LX Blue Spoke	1	\$830.92
1w	100997	25" x 1/2" Spinergy LX Yellow Spoke	1	\$830.92
1x	101739	25" x 1/2" Spinergy LX Pink Spoke	1	\$830.92
1y	101458	25" x 1/2" Spinergy LX Green Spoke	1	\$830.92
1z	200231	25" x 1/2" Spinergy LX Orange Spoke	1	\$830.92
1aa	105106	25" x 1/2" Spinergy LX Purple Spoke	1	\$830.92
1ab	200213	24" Spinergy LX w/ Flex Rim Black Spoke	1	\$1,359.51
1ac	200215	24" Spinergy LX w/ Flex Rim Red Spoke	1	\$1,359.51

Item Number	Part Number	Description	Quantity	Retail
1ad	200216	24" Spinergy LX w/ Flex Rim White Spoke	1	\$1,359.51
1ae	200214	24" Spinergy LX w/ Flex Rim Blue Spoke	1	\$1,359.51
1af	200225	24" Spinergy LX w/ Flex Rim Yellow Spoke	1	\$1,359.51
1ag	200218	24" Spinergy LX w/ Flex Rim Pink Spoke	1	\$1,359.51
1ah	200211	24" Spinergy LX w/ Flex Rim Green Spoke	1	\$1,359.51
1ai	200217	24" Spinergy LX w/ Flex Rim Orange Spoke	1	\$1,359.51
1aj	107397	24" Spinergy LX w/ Flex Rim Purple Spoke	1	\$1,359.51
1ak	200208	25" Spinergy LX w/ Flex Rim Black Spoke	1	\$1,359.51
1al	200220	25" Spinergy LX w/ Flex Rim Red Spoke	1	\$1,359.51
1am	200221	25" Spinergy LX w/ Flex Rim White Spoke	1	\$1,359.51
1an	200219	25" Spinergy LX w/ Flex Rim Blue Spoke	1	\$1,359.51
1ao	101463	25" Spinergy LX w/ Flex Rim Yellow Spoke	1	\$1,359.51
1ap	200224	25" Spinergy LX w/ Flex Rim Pink Spoke	1	\$1,359.51
1aq	200223	25" Spinergy LX w/ Flex Rim Green Spoke	1	\$1,359.51
1ar	200222	25" Spinergy LX w/ Flex Rim Orange Spoke	1	\$1,359.51
1as	107398	25" Spinergy LX w/ Flex Rim Purple Spoke	1	\$1,359.51
2	102854	Spinergy LX Spoke Wrench <i>No charge when ordered with Spinergy Wheel or replacement spoke</i>	1	
	110119	Spinergy LX Replacement Spoke Tool Kit <i>This tool kit contains: Spinergy LX Spoke Wrench and Spinergy Spoke T-Wrench</i>	1	\$7.14
3a	101929	Replacement Spoke Spinergy LX 22", Black (511291-8.75BLA)	1	\$25.50
3b	101930	Replacement Spoke Spinergy LX 22", Red (511291-8.75RED)	1	\$25.50
3c	101931	Replacement Spoke Spinergy LX 22", White (511291-8.75WHT)	1	\$25.50
3d	101932	Replacement Spoke Spinergy LX 22", Blue (511291-8.75BLU)	1	\$25.50
3e	101933	Replacement Spoke Spinergy LX 22", Yellow (511291-8.75YEL)	1	\$25.50
3f	101934	Replacement Spoke Spinergy LX 22", Pink (511291-8.75PIN)	1	\$25.50
3g	101935	Replacement Spoke Spinergy LX 22", Green (511291-8.75GRE)	1	\$25.50

Item Number	Part Number	Description	Quantity	Retail
3h	101936	Replacement Spoke Spinergy LX 22", Orange (511291-8.75ORA)	1	\$25.50
3i	101937	Replacement Spoke Spinergy LX 24", Black (511291-9.5BLA)	1	\$25.50
3j	101938	Replacement Spoke Spinergy LX 24", Red (511291-9.5RED)	1	\$25.50
3k	101939	Replacement Spoke Spinergy LX 24", White (511291-9.5WHT)	1	\$25.50
3l	101940	Replacement Spoke Spinergy LX 24", Blue (511291-9.5BLU)	1	\$25.50
3m	101941	Replacement Spoke Spinergy LX 24", Yellow (511291-9.5YEL)	1	\$25.50
3n	101942	Replacement Spoke Spinergy LX 24", Pink (511291-9.5PIN)	1	\$25.50
3o	101943	Replacement Spoke Spinergy LX 24", Green (511291-9.5GRE)	1	\$25.50
3p	101944	Replacement Spoke Spinergy LX 24", Orange (511291-9.5ORA)	1	\$25.50
3q	101945	Replacement Spoke Spinergy LX 25", Black (511291-9.84BLA)	1	\$25.50
3r	101946	Replacement Spoke Spinergy LX 25", Red (511291-9.84RED)	1	\$25.50
3s	101947	Replacement Spoke Spinergy LX 25", White (511291-9.84WHT)	1	\$25.50
3t	101948	Replacement Spoke Spinergy LX 25", Blue (511291-9.84BLU)	1	\$25.50
3u	101949	Replacement Spoke Spinergy LX 25", Yellow (511291-9.84YEL)	1	\$25.50
3v	101950	Replacement Spoke Spinergy LX 25", Pink (511291-9.84PIN)	1	\$25.50
3w	101951	Replacement Spoke Spinergy LX 25", Green (511291-9.84GRE)	1	\$25.50
3x	101952	Replacement Spoke Spinergy LX 25", Orange (511291-9.84ORA)	1	\$25.50

# Handrim Hardware Chart

Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654  Natural Fit: 100835	100206		101106					
20" Ki Spoke	200530					200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349	101161	100577	200547	100793	200202		
25" Ki Spoke	200533					200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460	
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963					
22" 5-Spoke X Core	101962						101898		101964				
24" 5-Spoke X Core	100960						100975		100976	200546			
24" 5-Spoke X Core	100960						101893				100768	100769	
20" Spinergy Spox	See Spinergy Spox Page	6			Screw: 100669 Nut: 100657	103125		103179					
22" Spinergy Spox						100827		100808		100889	100888		
24" Spinergy Spox						100766	101161	100615		100830	100828		
25" Spinergy Spox						100767	101160	100765		101464	101460		
26" Spinergy Spox						101477		101148		200200	100950		
22" Spinergy LX	See Spinergy LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888		
24" Spinergy LX						100766	101161	100615		100830	100828	200213	
25" Spinergy LX						100767	101160	100765		101464	101460	200208	
26" Spinergy LX						101477		101148		200200	100950		
18" MAXX Spoke	109285	6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654  Natural Fit: 100835  Superlight: Screw: 100669 Nut: 100657	100206		101106					
20" MAXX Spoke	108244					200536		200542					
22" MAXX Spoke	105135					100560		100576	100569	200538	200201		
24" MAXX Spoke	105136					200349	101161	100577	200547	100793	200202		
25" MAXX Spoke	107436					200350		101870	200548	200539	200540		
26" MAXX Spoke	107437					200351		100578	200549	100907	101454		

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.

## Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part No.	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Marathon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke	200529	101587	101588	101137																	
20" Ki Spoke	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke	200532	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
24" Superlight	101159	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	100635	100589	100596				102982	100589	100596					102920	100589	100596				
20" Spinergy Spox	See Spinergy Spox Page	100911	100912	101137	102957	102844	101137														
22" Spinergy Spox		100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy Spox		100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy Spox					102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Spinergy Spox		100636	100590	100608	102960	100582	100597				100865		100582	100597							
22" Spinergy LX	See Spinergy LX Page	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy LX		100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy LX					102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Spinergy LX		100636	100590	100608	102960	100582	100597				100865		100582	100597							
18" MAXX Spoke	109285	101587	101588	101137																	
20" MAXX Spoke	108244	100911	100912	101137	102957	102844	101137														
22" MAXX Spoke	105135	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" MAXX Spoke	105136	100635	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" MAXX Spoke	107436				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" MAXX Spoke	107437	100636	100590	100608	102960	100582	100597				100865		100582	100597							

## Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part Number	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke	200532	100635	100622	100633	100635 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke	200534	100636			100636 / 100505	102978	100884	102960		100865			
24" Superlight	101159	100635					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
12" Mag <sup>1</sup>	N/A	001424	001426		001425								
16" Mag <sup>1</sup>	N/A	001159	001161	003453	001160								
20" 5-Spoke X Core (Black)	101961	100911	100917	102085	100911 / 100913		102717						
22" 5-Spoke X Core (Black)	101962	100634	100621	100632	100634 / 100503		102718						
24" 5-Spoke X Core (Black)	100960	100635	100622	100633	100635 / 100504	102976	100640		102982			102920	
20" Spinergy Spox	See Spinergy Spox Page	100911	100917	102085	100911 / 100913	102974		102957					
22" Spinergy Spox		100634	100621	100632	100634 / 100503	103165		102961					101475
24" Spinergy Spox		100635	100622	100633	100635 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy Spox						102977	101146	102959		101357	103177		100166
26" Spinergy Spox		100636				100636 / 100505	102978	100895	102960		100865		
22" Spinergy LX	See Spinergy LX Page	100634	100621	100632	100634 / 100503			102961					101475
24" Spinergy LX		100635	100622	100633	100635 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy LX						102977	101146	102959		101357	103177		100166
26" Spinergy LX		100636				100636 / 100505	102978	100895	102960		100865		
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	100635	100622	100633	100635 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly