


# CR45 Pelvic Positioning Belts - 1.5" Auto Buckle and Hardware



| Item Number | Part Number | Description   | Quantity |
|-------------|-------------|---|----------|
| 1a, 2, 3, 4 | 111387      | Belt Pelvic Positioning 1.5" Auto Buckle 56", Focus<br><i>Used with 14"-20" seat widths</i> | 1        |
| 1b, 2, 3, 4 | 111392      | Belt Pelvic Positioning 1.5" Auto Buckle 60", Focus<br><i>Used with 21"-22" seat widths</i> | 1        |
| 1a          | 100595      | Belt Pelvic Positioning 1.5" Auto Buckle 56"  | 1        |
| 1b          | 100847      | Belt Pelvic Positioning 1.5" Auto Buckle 60"  | 1        |
| 2           | 101125      | M6 x 55 BTN SHCS BLKZC  | 2        |
| 3           | 100746      | M6 Flat Washer Blk Zn   | 2        |
| 4           | 100658      | M6 Nylock Nut Blk Zn  | 2        |