


## Focus MAXX Performance Spoke Wheel

---



*If changing wheel size and the Attendant Foot Lock is present, a new linkage may be required; reference the corresponding Focus CR Attendant Foot Lock page. <p>*

*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel style, new axles may be needed. <p>*

*If Drum Brake is present, see Focus MAXX Spoke for Drum Brake page*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>	<b>Retail</b>
1, 2	105136	Wheel 24" MAXX Performance Spoke, Black Hub/Rim/Spokes	1	\$147.08
2	107408	Replacement Spoke MAXX Performance Spoke 24"	1	\$4.47

# Handrim Hardware Chart

Wheel		Handrim Connection				Handrim								
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim		
18" Ki Spoke <sup>1</sup>	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" Ki Spoke <sup>1</sup>	200530	6				Natural Fit: 100835	200536		200542					
22" Ki Spoke <sup>1</sup>	200531						100560		100576	100569	200538	200201		
24" Ki Spoke <sup>1</sup>	200532				200349			100577	200547	100793	200202			
25" Ki Spoke <sup>1</sup>	200533				200350		101870	200548	200539	200540				
26" Ki Spoke <sup>1</sup>	200534				200351		100578	200549	100907	101454				
24" Superlight	101159		6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460		
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963						
22" 5-Spoke X Core	101962					101898		101964						
24" 5-Spoke X Core	100960					100975		100976	200546					
24" 5-Spoke X Core	100960				101893				100768	100769				
20" Spinerger Spox	See Spinerger Spox Page	6			Screw: 100669 Nut: 100657	103125		103179						
22" Spinerger Spox						100827		100808		100889	100888			
24" Spinerger Spox						100766	101161	100615		100830	100828			
25" Spinerger Spox						100767	101160	100765		101464	101460			
26" Spinerger Spox						101477		101148		200200	100950			
22" Spinerger LX	See Spinerger LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888			
24" Spinerger LX						100766	101161	100615		100830	100828	See Spinerger LX Page		
25" Spinerger LX						100767	101160	100765		101464	101460			
26" Spinerger LX						101477		101148		200200	100950			
18" MAXX Spoke	109285	3	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654	100206		101106						
20" MAXX Spoke	108244	6				Natural Fit: 100835	200536		200542					
22" MAXX Spoke	105135						100560		100576	100569	200538	200201		
24" MAXX Spoke	105136				200349		101161	100577	200547	100793	200202			
25" MAXX Spoke	107436				200350		101870	200548	200539	200540				
26" MAXX Spoke	107437				200351		100578	200549	100907	101454				

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.

<sup>1</sup> Discontinued

## Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part No.	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Marathon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke <sup>1</sup>	200529	101587	101588	101137																	
20" Ki Spoke <sup>1</sup>	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke <sup>1</sup>	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke <sup>1</sup>	200532	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke <sup>1</sup>	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke <sup>1</sup>	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
24" Superlight	101159	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	111818	100589	100596				102982	100589	100596					102920	100589	100596				
20" Spinergy Spox	See Spinergy Spox Page	100911	100912	101137	102957	102844	101137														
22" Spinergy Spox		100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy Spox		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy Spox					102959	101457	101452				101357	103177	101457	101452					100166	100208	101452
26" Spinergy Spox		100636	100590	100608	102960	100582	100597				100865		100582	100597							
22" Spinergy LX	See Spinergy LX Page	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy LX		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy LX					102959	101457	101452				101357	103177	101457	101452					100166	100208	101452
26" Spinergy LX		100636	100590	100608	102960	100582	100597				100865		100582	100597							
18" MAXX Spoke	109285	101587	101588	101137																	
20" MAXX Spoke	108244	100911	100912	101137	102957	102844	101137														
22" MAXX Spoke	105135	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" MAXX Spoke	105136	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" MAXX Spoke	107436				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" MAXX Spoke	107437	100636	100590	100608	102960	100582	100597				100865		100582	100597							

<sup>1</sup> Discontinued

## Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part Number	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke <sup>2</sup>	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke <sup>2</sup>	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke <sup>2</sup>	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke <sup>2</sup>	200532	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke <sup>2</sup>	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke <sup>2</sup>	200534	100636			100636 / 100505	102978	100884	102960		100865			
24" Superlight	101159	111818					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
12" Mag (non-drum brake) <sup>1</sup>	N/A	001424	001426		001425								
16" Mag (non-drum brake) <sup>1</sup>	N/A	001159	001161	003453	001160								
20" 5-Spoke X Core (Black)	101961	100911	100917	102085	100911 / 100913		102717						
22" 5-Spoke X Core (Black)	101962	100634	100621	100632	100634 / 100503		102718						
24" 5-Spoke X Core (Black)	100960	111818	100622	100633	111818 / 100504		100640		102982			102920	
20" Spinergy Spox	See Spinergy Spox Page	100911	100917	102085	100911 / 100913	102974		102957					
22" Spinergy Spox		100634	100621	100632	100634 / 100503	103165		102961					101475
24" Spinergy Spox		111818	100622	100633	111818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy Spox						102977	101146	102959		101357	103177		100166
26" Spinergy Spox		100636				100636 / 100505	102978	100895	102960		100865		
22" Spinergy LX	See Spinergy LX Page	100634	100621	100632	100634 / 100503			102961					101475
24" Spinergy LX		111818	100622	100633	111818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy LX						102977	101146	102959		101357	103177		100166
26" Spinergy LX		100636				100636 / 100505	102978	100895	102960		100865		
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly

<sup>2</sup> Discontinued