


# Catalyst E Seat Upholstery



SCREW (ITEM 4) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	10
16" - 17"	12
18" - 19"	14
20"	16

*Seat upholstery width measures approximately 2.5" less than the chair size. Seat upholstery depth measures true to the chair size.*

Item Number	Part Number	Description	Quantity	Retail
1a	100095	Seat Upholstery Cat 1416	1	\$94.55

Item Number	Part Number	Description	Quantity	Retail
1b	100097	Seat Upholstery Cat 1418	1	\$94.55
1c	100107	Seat Upholstery Cat 1616	1	\$94.55
1d	100109	Seat Upholstery Cat 1618	1	\$94.55
1e	100111	Seat Upholstery Cat 1620	1	\$94.55
1f	100121	Seat Upholstery Cat 1816	1	\$94.55
1g	100123	Seat Upholstery Cat 1818	1	\$94.55
1h	100125	Seat Upholstery Cat 1820	1	\$94.55
1i	100131	Seat Upholstery Cat 2016	1	\$94.55
1j	100133	Seat Upholstery Cat 2018	1	\$94.55
1k	100135	Seat Upholstery Cat 2020	1	\$94.55
1l	101556	Seat Upholstery Cat 2216	1	\$94.55
1m	100247	Seat Upholstery Cat 2218	1	\$105.06
1n	100261	Seat Upholstery Cat 2220	1	\$105.06
2a	100232	Seat Rail 16" Deep CAT	2	\$17.06
2b	100234	Seat Rail 18" Deep CAT	2	\$17.06
2c	100236	Seat Rail 20" Deep CAT	2	\$17.06
3a	000059	Lifting Strap 14" Wide CAT	1	\$6.39
3b	000061	Lifting Strap 16" Wide CAT	1	\$6.39
3c	000063	Lifting Strap 18" Wide CAT	1	\$6.39
3d	000058	Lifting Strap 20" Wide CAT	1	\$6.39
3e	100875	Lifting Strap 22" Wide CAT	1	\$6.40
4	100734	M5 x 16 PHMS ZC	16	\$1.58