


CR45 24" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and handrim kits will be listed. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size, new axles/receivers and handrims may be needed. If changing wheel size and Attendant Foot Lock or T-Arms are present, spacer 112503 is required which can be found on the CR45 Axle Plate page See CR45 Axles page for axle and receiver utilization.

Item Number	Part Number	Description	Quantity
	118359	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118360	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118361	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118362	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118367	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118368	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118369	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118370	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	118136	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118137	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118138	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118139	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118144	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118145	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118146	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	118147	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount	1
	117696	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117697	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117698	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117699	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117704	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117705	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117706	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	117707	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount	1