



KidWalk

Dynamic Mobility System



Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.1 Rev A

Client _____ Age _____ Height _____ Weight _____ Dx _____
Facility _____ Therapist _____ Ph. _____
Dealer _____ PO/RFQ: _____
Contact _____ Ph. _____ Date _____
Email _____

KidWalk is available in three sizes: KidWalk EI, KidWalk, and KidWalk II

Standard Features Include: Auto-centering pelvic guides, lateral weight shift, vertical weight relief with adjustable load compensating spring assembly, angle adjustable mast, anti-tippers, wheel locks, quick release wheels, padded seat, adjustable axle position, light up rear casters, two piece frame for easy transport.

KidWalk EI features a wrap around semi-rigid soft lateral system to fit the smaller client comfortably with a minimum of bulky hardware. Ultra low seat height for true EI use and natural gait with dynamic mobility.

KW EI



KW



KW II



Specifications	KidWalk EI	KidWalk	KidWalk II
Weight Capacity	65 Lbs.	65 Lbs.	90 Lbs.
Base	W 24" x L 28"	W 24"x L 28"	W 27.5" x L 32"
Seat Height	8.5" - 18"	12" - 22"	15" - 28"
Mast Height	28"	33.5"	37.5"
Mag Wheels	16"	20"	24"
Lighted Casters	4"	4"	6"


Do Not exceed weight capacity.
Units Shown with optional equipment.
All specifications subject to change without notice.

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544

KidWalk Dynamic Mobility System	☐ KidWalk EI KIDWALKEI \$3,395.00	☐ KidWalk KIDWALK \$3,595.00	☐ KidWalk II KIDWALKII \$4,195.00	Retail Price
Frame Colors: ☐ Candy Red ☐ Candy Blue ☐ Candy Purple ☐ Neon Pink ☐ Neon Orange ☐ Shocker Green ☐ Razzberry ☐ White ☐ Black				
Hip System Configurations (See Width and Depth Measurements) Standard configuration will be supplied for orders with no selection. Note: KWII Depth Range 3.5" - 6": Supplied with non-extended hip pads.	Standard Size Ranges ☐ KW00570 W 5" - 7" D 3.5" - 5.25" Optional Size Ranges ☐ KW00500 W 5.25" - 9.25" D 4.25" - 7.5" ☐ KW00580 W 5" - 7" D 4.25" - 7.5" ☐ KW00590 W 5.25" - 9.25" D 3.5" - 5.25"	Standard Size Ranges ☐ KW00500 W 5.25" - 9.25" D 4.25" - 7.5" Optional Size Ranges ☐ KW00510 W 7.5" - 11.75" D 4.25" - 9" ☐ KW00520 W 5.25" - 9.25" D 4.25" - 9"	Standard Size Ranges ☐ KW00520 W 7.5" - 11.75" D 4.25" - 9" Optional Size Ranges ☐ KW10500 W 5.25" - 9.25" D 4.25" - 7.5" ☐ KW00530 W 7.5" - 11.75" D 3.5" - 6" ☐ KW00540 W 5.25" - 9.25" D 3.5" - 6" ☐ KW00550 W 8" - 15.5" D 4.25" - 7.5" ☐ KW00560 W 8"-15.5" D 3.5" - 6"	Included
Hip Support Flip-Up Pad Brackets (Standard configuration is Straight) <i>If no selection is made, Standard will be supplied.</i>	☐ KW00010 Straight (Std) ☐ KW00030 10 Deg. Angle	☐ KW00010 Straight (Std) ☐ KW00030 10 Deg. Angle	☐ KW00020 Straight (Std) ☐ KW00030 10 Deg. Angle	Included
Extended Hip Pads (Standard on 4.25" - 7.5" Depth) (Not available on 3.5" - 6" Depth)	☐ KW00040	☐ KW00040	☐ KW00040	\$ 139.00
Chest System Configurations (Rigid) Includes lateral and lumbar support pads and swing-up anterior support pads. (See Width and Depth Range Measurements) <i>If no selection is made, Standard will be supplied.</i>	Standard Size Ranges ☐ KW00600 W 6.5" - 10.5" D 4" - 6.75"	Standard Size Ranges ☐ KW00600 W 6.5" - 10.5" D 4" - 6.75" Optional ☐ KW00610 W 8" - 15.5" D 6.5" - 9.5" ☐ KW00620 W 6.5" - 10.5" D 6.5" - 9.5" ☐ KW00630 W 8" - 15.5" D 4" - 6.75"	Standard Size Ranges ☐ KW00610 W 8" - 15.5" D 6.5" - 9.5" Optional ☐ KW00600 W 6.5" - 10.5" D 4" - 6.75" ☐ KW00630 W 8" - 15.5" D 4" - 6.75" ☐ KW00620 W 6.5" - 10.5" D 6.5" - 9.5"	\$ 597.00
Chest Support Flip-Up Pad Brackets (Standard: 10 Degree)	☐ KW00650 10 Deg. ☐ KW00640 Straight	☐ KW00650 10 Deg. ☐ KW00640 Straight	☐ KW00650 10 Deg. ☐ KW00640 Straight	Included
EI Upper Body Soft Lateral System Flexible wrap around support with side to side adjustable lateral support.	☐ KW00660 Width 4" - 7" ☐ KW00670 Width 6" - 9"	N/A	N/A	\$ 549.00
PRIME ENGINEERING A division of Axiom Industries 4202 W. Sierra Madre Ave. Fresno, CA 93722 Ph (800)827-8263 (559)271-0991 Email: info@primeengineering.com Website: www.primeengineering.com Fax (800)800-3355 (559)276-3544				



KidWalk

Dynamic Mobility System




Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.1 Rev A

	KidWalk EI	KidWalk	KidWalk II	
Upper Body Soft Support	☐ KW00680	☐ KW00680	☐ KW00680	\$ 299.00
Headrest 3" x 4" Includes mount	☐ KW00050	☐ KW00050	☐ KW00050	\$289.00
Contoured Small Headrest 4" x 6" Incl. mount	☐ KW00060	☐ KW00060	☐ KW00060	\$ 379.00
Contoured Med. Headrest 5" x 8" Incl. mount	☐ KW00070	☐ KW00070	☐ KW00070	\$ 449.00
1/2 Square Tube Mount Headrest not included	☐ KW00080	☐ KW00080	☐ KW00080	\$ 175.00
1/2 Flip Down Headrest Mount	☐ KW00090	☐ KW00090	☐ KW00090	\$ 399.00
I2I Flip Down Headrest Mount	☐ KW00100	☐ KW00100	☐ KW00100	\$ 399.00
Wide Seat Large padded seat 4" x 10"	☐ KW00110	☐ KW00110	☐ KW00120	\$ 275.00
Extended Seats Longer sides assist with controlling abduction and scissoring. <i>Select vertical seat measurement.</i>	☐ KW00130 Short 4" ☐ KW00150 Med 6" ☐ KW00170 Long 8"	☐ KW00130 Short 4" ☐ KW00150 Med 6" ☐ KW00170 Long 8"	☐ KW00140 Short 4" ☐ KW00160 Med 6" ☐ KW00180 Long 8"	\$ 385.00
Additional Standard Seat (Ordered with unit)	☐ KW00210	☐ KW00190	☐ KW00200	\$ 305.00
Dynamic Leg Wraps Assist in supplying feedback to the client while encouraging movement in the lower extremities.	☐ KW00220	☐ KW00220	☐ KW00220	\$ 149.00
Leg Straps Reduces movement of lower extremities. (Requires 6" or 8" Extended Seat)	☐ KW00230	☐ KW00230	☐ KW00240	\$ 229.00
Sliding Ankle Prompts Helps discourage client's feet from crossing. Requires Extended Seat: 6" for KW (Sz 1). 8" for KWII (Sz. 2)	☐ KW00270	☐ KW00250	☐ KW00260	\$ 356.00
Anti-Rollback Device	☐ KW00290	☐ KW00280	☐ KW00290	\$ 359.00
Directional locks Restricts unit to straight line movement.	☐ KW00310	☐ KW00300	☐ KW00310	\$ 275.00
Non-lighted caster wheels	☐ KW00320	☐ KW00320	☐ KW00330	N/C
Cruiser Wheels Allows use as an omni-directional walker for clients up to 40" tall. (Not available on KidWalk II)	☐ KW00350	☐ KW00340	N/A	\$ 415.00
Spoke Wheels Optional replacement for standard Mag Wheels	N/A	☐ KW00360 20"	☐ KW00370 24"	\$ 499.00
Spoke Protectors	N/A	☐ KW00380 20"	☐ KW00390 24"	\$ 265.00
Mag Protectors	N/A	☐ KW00400	☐ KW00410	\$ 265.00

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544



KidWalk

Dynamic Mobility System



Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.1 Rev A



Swing-Open Hip and Chest
Lateral Support
Shown in open loading position



1/2" Square
Tube Mount



Head Pad



KidWalk EI Lateral system shown open and closed



Extended Seats

Client Measurements:

A: Chest Width _____" B: Chest Depth _____" C: Hip Width _____" D: Hip Depth _____" E: Inseam _____"

F: Top of Inseam to Axilla _____" G: Client Height _____" H: Client Weight _____lbs

Note: All Measurements Assume Symmetrical Standing Posture.

Custom Requests:

Please complete the information in the client measurement area above and list any special needs in the following area to speed up the design process and/or response time.

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544