

PO# / Client:_

Raz-CAT

MSRP – USA – USD Effective Nov 3, 2023

Date: ____

Compact Attendant Tilt

	•		
ealer / Location:		T:	F:
, =====================================			

______ Email: _

Raz-CAT –Compact Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height adjustable stainless steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. User weight capacity 300 lbs.

For help with justification, please see our **Evaluation** and **Medical Justification Form**.

_ Contact: ___

Need help with this order form? Call Raz customer service at (877) 720-5678

1. Standard Raz-CAT Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]

Z30318	Raz-CAT w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$3,395
--------	--	---------

2. Custom / Special Frame Widths, Depths and Heights - Adjusta-Back must be selected for all custom widths see [Sec.7]

Z30316	Raz-CAT w/ 16"W frame width¹ (base chair price – custom upcharge below will be added)	\$3,395				
Z3CSW16	Upcharge for custom 16"W seat and base frame – will not fit over most toilets	\$595				
Z30320	Raz-CAT w/ 20"W frame width ² (base chair price – custom upcharge below will be added)	\$3,395				
Z3CSW20	Upcharge for custom 20"W seat and base frame	\$595				
Custom Frame Configuration for 18"W frame. (Call Raz for codes and pricing for 16" / 20" widths)						
Z3CLOWF5	Upcharge for custom Low Frame (16"-20" STF)	\$1,595				
¹ 18"W Visco-Foam Interface seats can be modified to fit.						
² 18"W seats can be	18"W seats can be selected but must also select Seat Support Kit (Z4SSKCT) [Sec.5]					

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16"D, are now designated 17"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Standard CAT Seat – Visco-Foam w/ Tapered Oblong Aperture (thickness ~2.5")

Standard aperture width: 3.5" at front, 4.5" at rear

ZCSC1816	18"W	x 11"-16	D (adjust	able) sta	ndard w/	18"W fra	ame	N/C
ZCSC1616	16"W	x 11"-16'	D (adjust	able) sta	ndard w/	16"W fra	ame	N/C
ZCSC2016	20"W	x 11"-16'	D (adjust	able) sta	ndard w/	20"W fra	ame	N/C
Set seat depth a	it:	11"	12"	13"	14"	15"	16"	



DOC#: RSD012r04 Rev. Date: 2023-08-25

Other Seats Available – for seat depth adjustment up to 18" – Modified to fit CAT seat frame

Molded Seats (thickness ~2.5"-3")			Visco Foam Inter	face Seats (thickness ~3.5")	
ZMSC17FRO	18"W x 17"D – Front Access	N/C	ZVFIC17FRO	18"W x 17"D – Front Access	\$425
ZMSC17BRI	18"W x 17"D – Front Bridged	N/C	ZVFIC17BRI	18"W x 17"D – Front Bridged	\$425
			ZVFIC17LTO	18"W x 17"D – Left Access	\$425
			ZVFIC17RTO	18"W x 17"D – Right Access	\$425
			ZVFIC17RRO	18"W x 17"D – Rear Access	\$425

Seat Support Kit - These are used to support the sides of seats that are narrower than seat frames on which they are mounted. The reason for a narrower seat on a wider frame is typically to select a Molded Seat and avoid the upcharge of a custom seat. A typical example would be when a 20"W chair is needed but an 18"W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.

ZSSKCT Seat Support Kit for Raz-CAT \$145

4. Seats - Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. 18" is the standard width for CAT seat frames.

See Custom Seat Assessment Form – Resources tab at razdesigninc.com

ZCSEAT Custom Dual-Visco Foam Seat (enter specs below)	' '
ZCSCOV Custom Dual-Visco Foam Shower-Only Seat (no aperture, no drain hole)	\$595

Seat Frame Width: __ 16" / **18"** / 20"

Custom Seat Width: __ 16" / 17" / 18" / 19" / 20"

Custom Seat Depth: __ 11" / 12" / 13" / 14" / 15" / 16" / 17"

Access Opening: FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged

Aperture: Tapered Oblong: __"W (at front) __"W (at rear) [standard is 3.5"W at front and 4.5"W at rear]

Keyhole: _ _"W (3.0"-8.0"W – enter "NA" if aperture is not required)

- Oblong-shaped apertures can be specified in conjunction with any of the access openings listed.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

5. Back Frames - for standard 18"W seat frames (back heights measured from top of seat to top of upholstery)

Multi-Position Adjusta-Back: The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 11"-16". Includes tension-adjustable vinyl back support, flip-up padded arm supports, back angle adjustment (90° / 97°) with cane-style back posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

Back frames fo	r standard 18"W seat frame (Z30318)	
ZCAB845	14" Back Width / 15" Back Height	N/C
ZCAB848	14" Back Width / 18" Back Height	\$95
ZCAB865	16" Back Width / 15" Back Height	N/C
ZCAB868	16" Back Width / 18" Back Height	\$95
ZCAB885	18" Back Width / 15" Back Height	N/C
ZCAB888	18" Back Width / 18" Back Height	\$95
ZDTL	Dual Triggers – left side ¹	N/C
ZDTR	Dual Triggers – right side ¹	N/C

¹ Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).

•		
ZRBC848	Recliner Back for CAT – 14"W x 18"H ^{2,3}	\$1,095
ZRBC868	Recliner Back for CAT – 16"W x 18"H ^{2,3}	\$1,095
ZRBCR888	Recliner Back for CAT – 18"W x 18"H ^{2,3}	\$1,095

DOC#: RSD012r04 Rev. Date: 2023-08-25

Reclining Backs for 18"W seat frame

Standard trigger locations for Raz-CAT with recline are: tilt triggers on the right and recline triggers on the left.
 Tilt is limited to 25°

6. Back Frames for 16"W and 20"W Seat Frames (back heights measured from top of seat to top of upholstery)

Multi-Position Adjusta-Back: The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 11"-16". Includes tension-adjustable vinyl back support, back angle adjustment (90° / 97°) and flip-up padded arm supports. Back upholstery can be positioned 2.5" higher than factory settings stated below.

			ı		
16"W Seat Fra	ames		•	s for 16"W seat frame	
ZCAB645	14" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclini	ng Backs.
ZCAB648	14" Back Width / 18" Back Height	\$95	ZRBC648	Recliner Back for CAT – 14"W x 18"H ^{2,3}	\$1,095
ZCAB665	16" Back Width / 15" Back Height	N/C	ZRBC668	Recliner Back for CAT $- 16$ "W x 18 "H 2,3	\$1,095
ZCAB668	16" Back Width / 18" Back Height	\$95			
20"W Seat Fra	ame		Reclining Backs	s for 20"W seat frame	
ZCAB065	16" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclini	ng Backs.
ZCAB068	16" Back Width / 18" Back Height	\$95	ZRBC088	Recliner Back for CAT $- 18$ "W x 18 "H 2,3	\$1,095
ZCAB085	18" Back Width / 15" Back Height	N/C			
ZCAB088	18" Back Width / 18" Back Height	\$95			
ZCAB005	20" Back Width / 15" Back Height	N/C			
ZCAB008	20" Back Width / 18" Back Height	\$95			
ZDTL	Dual Triggers – left side ¹	N/C			
ZDTR	Dual Triggers – right side ¹	N/C			
¹ Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).				ger locations for Raz-CAT with recline are tilt t cline triggers on left. to 25°	riggers on

7. Back Support Options – These back supports replace tension-adjustable upholstery on an Adjusta-Back. Symphony Backs are recommended for infection-control environments (for 18"H backs only).

ZSB14	Symphony Back – recessed solid back support – 14" back frame width	\$295
ZSB16	Symphony Back – recessed solid back support – 16" back frame width	\$295
ZSB18	Symphony Back – recessed solid back support – 18" back frame width	\$295
ZSB20	Symphony Back – recessed solid back support – 20" back frame width	\$395

DOC#: RSD012r04 Rev. Date: 2023-08-25

8. Arm Supports and Accessories

74.54.45	4	·	75444	Di .: A AA 6	44.65
ZARMS	Short Arm Supports ¹ (pr) – 10"L pad	Std	ZPAML	Pivoting Arm Mount – Left (for ZFAPLL/ZTRL)	\$165
ZARML	Long Arm Supports (pr) – 13"L pad	N/C		• •	4.55
ZALK	Arm Support Locks (pr) for Fixed-Back (ZFB1818)	\$75	ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR / ZTRR)	\$165
Z2PLT	2-Point Arm Supports ² (pr) with 13"L pads (std with reclining backs)	\$345	ZASK05	Arm Support Spacer Kit $-1/2$ " per side for all Adjusta-Backs (pr)	\$60
Z2PST	2-Point Arm Supports ² (pr) with 10"L pads (std with reclining backs)	\$345	ZASK10	Arm Support Spacer Kit -1 " per side for all Adjusta-Backs (pr)	\$115
ZFAPSL	Short Molded Flat Armpad ³ – Left	\$90	ZPHGS	Pivoting Hand Grips (pr) for 10"L arms ⁴	\$395
2.7 02	(10"L x 3"W – fits ZARM10L)	750	ZPHGL	Pivoting Hand Grips (pr) for 13"L arms ⁴	\$395
ZFAPSR	Short Molded Flat Armpad ³ – Right (10"L x 3"W – fits ZARM10R)	\$90	ZAPSB*	Anterior Posture Support Bar ^{5,6} * Enter frame width to complete code	\$225
ZFAPLL	Long Molded Flat Armpad ³ – Left (14"L x 3.25"W – fits ZARM13L)	\$90		(14" / 16" / 18" / 20")	
ZFAPLR	Long Molded Flat Armpad ³ – Right (14"L x 3.25"W – fits ZARM13R)	\$90			
ZTRL	Molded Arm Trough ³ – Left (18.5″L x 4.5″W)	\$110			
ZTRR	Molded Arm Trough ³ – Right (18.5″L x 4.5″W)	\$110			
¹ Not compati	¹ Not compatible with ZFAPL / ZFAPS / ZTR / ZPAM			combined with 2 Doint Arms requires a sustain	
'	² For chairs with same back width and frame width only			combined with 2 Point Arms requires a custon	ı up-
l <u>-</u>				with ZFAPL / ZFPS / ZTR	
depending on Trough	medial / lateral position of the Flat Arm Pad or Ar	m	l_ '	required with ZASK05 (+1") and ZASK10 (+2")	

9. Foot and Leg Supports

	Log Capports				
ZFS	Foot Supports (pr) – 16"-20.5"	Std	ZELRL	Elevating Leg Support ³ – Left	\$495
ZFS16	Foot Supports for 16"W frames (pr)	N/C		(includes ZADML)	
	– 16"-20.5"		ZELRR	Elevating Leg Support ³ – Right	\$495
ZFET2	Footrest Extension Tubes (pr)	\$50		(includes ZADMR)	
	– 2" longer ¹		ZFOR	Forward Offset Foot Support Receiver –	\$95
ZMFX8C	MFX8 Foot Supports (pr) – 8"-12"	N/C		(ea) Moves foot support forward 2.5" to	
ZMFX12C	MFX12 Foot Supports (pr) – 12"-16.5"	N/C		accomodate deeper seats and reduced knee flexion range. Specify 1 or 2	
ZADM	Adj Angle/Depth Molded Footplates ² (pr)	\$120	ZRLSL	Universal Residual Limb Support – Left	\$195
ZADABS16	Angle/Depth Adj ABS Footplates – Left	\$195	ZRLSR	Universal Residual Limb Support – Leit	\$195
	For 16"W frames		ZNLSK	Oniversal Residual Limb Support – Right	\$195
ZVS*	V-Style Foot Support (1-piece) 16"-20.5"	\$295			
	* Enter frame width to complete code (16" / 18")				
ZVS20	V-Style Foot Support (1-piece) 16"-20.5" For 20"W frames	\$345			
7507		Ċ 4 E E			
ZFBT	Flip-Back Foot Support (1-piece) 16"-20.5"	\$455			
1 -	¹ Not compatible with ZMFX8C / ZADM / ZVS / ZVSW / ZFBT			patible with frame widths ≥ 18" and adjusted se	at
² Not compatible	e with Z30316 or ZMFX8C		depths ≤ 18" c	or ≤ 21" w/ 3" longer seat frames (Z3LNGSF3)	

10. Foot and Leg Support Accessories

	<u> </u>				
ZFSLK	Foot Support Locks (pr)	\$95	ZNC	Neoprene Cover for Standard	\$90
ZCS* Calf Strap, Snap-On – Bodypoint Aeromesh – For 16" / 18" / 20"W frames		\$105		Footplate (pr)	
			ZNCA	Neoprene Cover for Angle/Depth	\$90
ZCSIC*	Calf Strap, Snap-On – Infection Control	\$95		Adjustable Footplates	
	For 16" / 18" / 20"W base frames		ZNC16	Neoprene Covers for Narrow Molded	\$90
ZCPL Calf Panel – Bodypoint Aeromesh ¹		\$225		Footplates – 16"W frame (pr)	
	Fits frame widths (16" / 18" / 20")		ZNCABS16	Neoprene Cover for Narrow Molded ABS	\$90
ZHS1820				Angle/Depth Adjustable Footplates	
– hook and loop attachment		,	ZNCVS	Neoprene Cover for V-Style Foot Support	\$90
ZAHL Adjustable Heel Loop (pr) \$ (not available for ZVS / ZFBT)		\$60	ZNCFB	Neoprene Cover for Flip-Back Foot	\$90
				Support	
		ZLMOR	Lateral / Medial Offset Foot Support	\$95	
				Receiver (ea) – moves foot support 2"	,,,,
				laterally or medially. Specify 1 or 2	
		ZFOR	Forward Offset Foot Support Receiver –	\$95	
				(ea) Moves foot support forward 2.5" to accomodate deeper seats and reduced	·
¹ Not compatible with ZMFX8C / ZMFX12C / ZVS / ZFBT				knee flexion range. Specify 1 or 2	

11. Positioning Pads – Hip Adductors are not compatible with Reclining Backrest unless ordered on-chair.

Head Suppo	rts			
ZHSP1	Head Support Pad – Small (10" x 4")	Standard		
ZHSP2	Head Support Pad – Large (12" x 6")	\$50		
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$75		
	acic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to ZLAEP extension plate.	4.5" medial		
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	\$225		
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	\$225		
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	\$225		
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	\$225		
Leg Adducto	ors			
ZLGT1L	Leg Adductor w/ Small Pad – Left	\$225		
ZLGT1R	Leg Adductor w/ Small Pad – Right	\$225		
ZLGT2L	Leg Adductor w/ Large Pad – Left	\$225		
ZLGT2R	Leg Adductor w/ Large Pad – Right	\$225		
Hip Adductors	for 4"narrower backs			
ZHP4C1L	Hip Adductor for 4" narrower back w/ Small Pad– Left	\$225		
ZHP4C1R	Hip Adductor for 4" narrower back w/ Small Pad– Right	\$225		
ZHP4C2L	Hip Adductor for 4" narrower back Large Pad– Left	\$225		
ZHP4C2R	Hip Adductor for 4" narrower back w/ Large Pad– Right	\$225		
Hip Adductors	for backs that are same width as seat frame or 2" narrower			
ZHP02C1L	Hip Adductor for 0-2" narrower back w/ Small Pad – Left	\$225		
ZHP02C1R	Hip Adductor for 0-2" narrower back w/ Small Pad – Right			
ZHP02C2L	· ·			
ZHP02C2R	Hip Adductor for 0-2" narrower back w/ Large Pad – Right	\$225		
Hip Adductors	For Reclining Backs			
ZHPRC1L	Hip Adductor for Reclining Back w/ Small Pad – Left	\$225		
ZHPRC1R	Hip Adductor for Reclining Back w/ Small Pad – Right	\$225		
ZHPRC2L	Hip Adductor for Reclining Back w/ Large Pad – Left	\$225		
ZHPRC2R	Hip Adductor for Reclining Back w/ Large Pad – Right	\$225		
	es for Lateral Thoracic Supports and Hip / Leg Adductors llow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors			
ZLAEP	Lateral/Adductor Extension Plate (ea)	\$28		
Lateral Thorac	ic Support Left Right			
Leg Adductor	Left Right			
Hip Adductor	Left Right			
Total Quantity	_			
Medial Thigh S	Supports (Pommel)			
ZMTST1	Medial Thigh Support w/ Small Pads	\$495		
ZMTST2	Medial Thigh Support w/ Large Pads	\$495		

12. Pelvic and Chest Belts - A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts			
ZPBS	Pelvic Belt – Small (frame widths 14-16", 2-piece w/ side-release buckle)	\$50	ZCBS	Chest Belt ² – Small (I 1-piece w/ hook and	•	\$45
ZPBM	Pelvic Belt – Medium (frame widths 16"- 18", 2-piece w/ side-release buckle)	\$50	ZCBM	Chest Belt ² – Mediur 18", 1-piece w/ hook	•	\$45
ZPBL	Pelvic Belt – Large (frame widths 18"-20", 2-piece w/ side-release buckle)	\$50	ZCBL		back widths 18"-20",	\$45
ZPBICS	Pelvic Belt ¹ – Infection Control, Small 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 14"-16")	\$75	ZCBICS	Chest Belt ^{1,2} – Infect 1-piece w/ side-relea widths 14-16")	ion Control, Small	\$75
ZPBICM	Pelvic Belt ¹ – Infection Control, Medium 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 16"-18")	\$75	ZCBICM	Chest Belt ^{1,2} – Infect Medium 1-piece w/ (back widths 16"-18'	side-release buckle	\$75
ZPBICL	Pelvic Belt ¹ – Infection Control, Large 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 18"-20")	\$85	ZCBICL	Chest Belt ^{1,2} – Infect 1-piece w/ side-relea widths 18"-20")	ion Control, Large	\$85
ZPBBPS	Pelvic Belt – Bodypoint Aeromesh, Small	\$245	ZCBBPS	Chest Belt – Bodypo	int Aeromesh Small	\$275
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh, Med	\$215	ZCBBPM	Chest Belt – Bodypo		\$245
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh, Large	\$215	ZCBBPL	Chest Belt – Bodypo		\$245
Recommende	Pelvic belt ZPBBPS N/A		ZPBBPM 16"V	N base frame	ZPBBPL 18"/20"W bas	e frame
sizes for Bodypoint bel	Chest belt ZCBBPS 14"W back frame		ZCBBPM 16"V	N back frame	ZCBBPL 18"/20"W bac	k frame

 $^{^{1}\,\}mbox{Non-absorbent}$ belts — recommended for infection-control environments

13. Urine Management

	•			
ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67		
ZSS	SplaShield – Molded urine deflector (must be held in place while using)	\$125		
ZSC17	Splash Curtain for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat $1/2$ ").	\$325		
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$345		
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345		
ZBUDBH	Bard Urine Drainage Bag Holder	\$155		
¹ Bridged Front	¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for order codes and pricing.			

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

14. Other

ZDCLT	Directional Caster Lock – fits on one	\$165	ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$200
	rear caster (facilitates steering)		ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$215
ZATTRR	Anti-Tippers – Rear (pr)	\$345	ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$230
ZATTFR	Anti-Tippers – Front (pr)	\$345		·	
ZCUT4	4" Caster Upgrade (4x)	\$120	ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$245
ZCUT6	6" Caster Upgrade (4x)	\$350	7SPKC22	CAT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	\$1,425
ZTLS	ZTLS Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for	N/C	ZSPKC24	CAT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	\$1,425
additional tilt limiting. Quantity				n Spacer Kits provide extra clearance for self-ca creases height required to clear toilet.	are with

15. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907	Omit Foot Supports	- \$45	Z912	Omit Complete Head Support Assembly	- \$50

Some assembly of Raz products is required.

16. Notes

Specifications Raz-CAT (Compact Attendant Tilt)

Weight Capacity	300 lbs (136 kg)		
Overall Width	Seat width plus 4.5" (with Standard Arm Support Pads)		
Overall Length (with MFX12 foot supports and headrest adjusted to middle positions)			
at 0° tiltat 40°	36.75" (93.5cm) (rear casters trailing) 44" (112cm)		
Seat Depth	11"-16" (28-41cm) depending on seat/cushion length		
Clearance Height (under seat) • with pan brackets • without pan brackets	17.25"-21.25" (44-54cm) (1" increments) 17.75"-21.75" (45-55cm) (1" increments)		
Clearance Width (to fit around toilet)	15.75"(40cm) for standard 18"W (46cm) seat frame. For other widths, seat frame width minus 2.25" (5.75cm)		
Arm Support Height - above molded seat (1/2" increments) (1.25cm)			
15"H Adjusta-Back18"H Adjusta-Back18"H Reclining Adjusta-Back	5.5"-8" (14-20cm) 6.5"-9" (16.5-23cm) 6.5"-9" (16.5-23cm)		
Seat-to-Floor Height (with molded seat)	21"-25" (53.5-63.5cm)		
Product Weight	55 lbs (25 kg)		
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)		