


Clik Push and Pull to Lock Wheel Locks



Part #	Bar	Length
003415	Standard	3.5"
001321	Extended	6.0"
002574	Bent	6.0"
003106	Extended Bent	8.0"

Please confirm wheel lock bar length and style for correct replacement parts

Item Number	Part Number	Description	Quantity
1a, 4b, 5	503206	Wheel Lock Assembly, Push Right w/ Instruction <i>This kit contains 200251 wheel lock along with an instruction for reversing the wheel lock</i>	1
1b, 4b, 5	503205	Wheel Lock Assembly, Push Left w/ Instruction <i>This kit contains 200252 wheel lock along with an instruction for reversing the wheel lock</i>	1

Item Number	Part Number	Description	Quantity
2a, 4b, 5	503207	Wheel Lock Assembly, Pull Right w/ Instruction <i>This kit contains 200253 wheel lock along with an instruction for reversing the wheel lock</i>	1
2b, 4b, 5	503208	Wheel Lock Assembly, Pull Left w/ Instruction <i>This kit contains 200254 wheel lock along with an instruction for reversing the wheel lock</i>	1
1a, 3, 4b, 5	115311	Wheel Lock Assembly, Push w/ Extension, Right with Instruction <i>This kit contains 112041 which is temporarily replacing the 200410 wheel lock. It contains an instruction for reversing the wheel lock.</i>	1
1b, 3, 4b, 5	115312	Wheel Lock Assembly, Push w/ Extension, Left with Instruction <i>This kit contains 112042 which is temporarily replacing the 200411 wheel lock. It contains an instruction for reversing the wheel lock.</i>	1
2a, 3, 4a, 5	503211	Wheel Lock Assembly, Pull Extension Right w/ Instruction <i>This kit contains 200412 wheel lock along with an instruction for reversing the wheel lock</i>	1
2b, 3, 4a, 5	503212	Wheel Lock Assembly, Pull Extension Left w/ Instruction <i>This kit contains 200413 wheel lock along with an instruction for reversing the wheel lock</i>	1
1a, 4b, 6, 9	503063	Wheel Lock Push w/ Extended Bar Right w/ Instruction <i>This kit contains 200251 wheel lock along with an instruction for reversing the wheel lock</i>	1
1b, 4b, 6, 9	503064	Wheel Lock Push w/ Extended Bar Left w/ Instruction <i>This kit contains 200252 wheel lock along with an instruction for reversing the wheel lock</i>	1
2a, 4b, 6, 9	503065	Wheel Lock Pull w/ Extended Bar Right w/ Instruction <i>This kit contains 200253 wheel lock along with an instruction for reversing the wheel lock</i>	1
2b, 4b, 6, 9	503066	Wheel Lock Pull w/ Extended Bar Left w/ Instruction <i>This kit contains 200254 wheel lock along with an instruction for reversing the wheel lock</i>	1
1a, 3, 4a, 6, 9	503067	Wheel Lock Push w/ Extension and Extended Bar Right w/ Instruction <i>This kit contains 200410 wheel lock along with an instruction for reversing the wheel lock</i>	1
1b, 3, 4a, 6, 9	503068	Wheel Lock Push w/ Extension and Extended Bar Left w/ Instruction <i>This kit contains 200411 wheel lock along with an instruction for reversing the wheel lock</i>	1
2a, 3, 4a, 6, 9	503069	Wheel Lock Pull w/ Extension and Extended Bar Right w/ Instruction <i>This kit contains 200412 wheel lock along with an instruction for reversing the wheel lock</i>	1
2b, 3, 4a, 6, 9	503070	Wheel Lock Pull w/ Extension and Extended Bar Left w/ Instruction <i>This kit contains 200413 wheel lock along with an instruction for reversing the wheel lock</i>	1

Item Number	Part Number	Description	Quantity
1a, 4b, 7	503071	Wheel Lock Push w/ Bent Bar Right w/ Instruction <i>This kit contains 200251, used with Offset Frames</i>	1
1b, 4b, 7	503072	Wheel Lock Push w/ Bent Bar Left w/ Instruction <i>This kit contains 200252, used with Offset Frames</i>	1
2a, 4b, 7	503073	Wheel Lock Pull w/ Bent Bar Right w/ Instruction <i>This kit contains 200253, used with Offset Frames</i>	1
2b, 4b, 7	503074	Wheel Lock Pull w/ Bent Bar Left w/ Instruction <i>This kit contains 200254, used with Offset Frames</i>	1
1a, 3, 4a, 7	503075	Wheel Lock Push w/ Extension and Bent Bar Right w/ Instruction <i>This kit contains 200410, used with Offset Frames</i>	1
1b, 3, 4a, 7	503076	Wheel Lock Push w/ Extension and Bent Bar Left w/ Instruction <i>This kit contains 200411, used with Offset Frames</i>	1
2a, 3, 4a, 7	503077	Wheel Lock Pull w/ Extension and Bent Bar Right w/ Instruction <i>This kit contains 200412, used with Offset Frames</i>	1
2b, 3, 4a, 7	503078	Wheel Lock Pull w/ Extension and Bent Bar Left w/ Instruction <i>This kit contains 200413, used with Offset Frames</i>	1
1a, 4a, 8, 9	122520	Wheel Lock Push w/ Extended Bent Bar Right w/ Instruction <i>This kit contains 200251, used as needed with Offset Frames</i>	1
1b, 4a, 8, 9	122521	Wheel Lock Push w/ Extended Bent Bar Left w/ Instruction <i>This kit contains 200252, used as needed with Offset Frames</i>	1
2a, 4a, 8, 9	122522	Wheel Lock Pull w/ Extended Bent Bar Right w/ Instruction <i>This kit contains 200253, used as needed with Offset Frames</i>	1
2b, 4a, 8, 9	122523	Wheel Lock Pull w/ Extended Bent Bar Left w/ Instruction <i>This kit contains 200254, used as needed with Offset Frames</i>	1
1a, 3, 4a, 8, 9	122524	Wheel Lock Push w/Extension and Extended Bent Bar Left w/ Instruction <i>This kit contains 200411, used as needed with Offset Frames</i>	1
1b, 3, 4a, 8, 9	122525	Wheel Lock Push w/Extension and Extended Bent Bar Right w/ Instruction <i>This kit contains 200410, used as needed with Offset Frames</i>	1
2a, 3, 4a, 8, 9	122526	Wheel Lock Pull w/Extension and Extended Bent Bar Left w/ Instruction <i>This kit contains 200413, used as needed with Offset Frames</i>	1
2b, 3, 4a, 8, 9	122527	Wheel Lock Pull w/Extension and Extended Bent Bar Right w/ Instruction <i>This kit contains 200412, used as needed with Offset Frames</i>	1

Item Number	Part Number	Description	Quantity
3	200369	Wheel Lock Extension <i>6" long extension</i>	1
4a	102336	Wheel Lock Cover Plastic <i>Used with Kits 503211 - 503212, 503067 - 503070, and 503075 - 503078, and 122520 - 122527.</i>	1
4b	111247	Wheel Lock Cover Plastic <i>Used with Kits 503205 - 503208, 115311 - 115312, 503063 - 503066, 503071 - 503074.</i>	1
5	003415	Standard Wheel Lock Bar	1
6	001321	Wheel Lock Mount Bar Extended <i>Used as needed</i>	1
7	002574	Wheel Lock Bar Bent <i>Used with Offset Frame</i>	1
8	003106	Extended Wheel Lock Bar Bent <i>Used as needed with Offset Frame</i>	1
9	002835	End Cap Plastic 3/4" <i>Used as needed with Extended Wheel Lock Bars</i>	1