

Front Seat Height in centimeters

		Fork	Frame	43	44	45	46	47	48	49	50	51	52	53	54	55	
Front Caster	3"	Multi Pos.	low	0	0												
		Single Side		0	0												
		Frog Leg				0	0	0	0								
		Multi Pos.	high			0	0	0	0	0	0						
		Single Side				0	0	0	0	0	0						
		Frog Leg								0	0	0	0				
	4"	Multi Pos.	low	0	0	0	0										
		Single Side		0	0	0	0										
		Frog Leg					0	0	0								
		Multi Pos.	high					0	0	0	0	0	0				
		Single Side						0	0	0	0	0	0				
		Frog Leg								0	0	0	0	0			
5"	Multi Pos.	low			0	0	0	0									
	Single Side				0	0	0	0									
	Frog Leg						0	0	0	0							
	Multi Pos.	high							0	0	0	0	0				
	Single Side								0	0	0	0	0				
	Frog Leg										0	0	0	0			
6"	Multi Pos.	low					0	0	0	0							
	Single Side						0	0	0	0							
	Frog Leg									0	0						
	Multi Pos.	high									0	0	0	0			
	Single Side										0	0	0	0			
	Frog Leg													0	0		