


Clik / Rogue XP Fixed Height Back Post



For item 1, see Little Wave Clik Adjustable Height Rigidizer Bar page or Rogue XP Adjustable Height Rigidizer Bar page. <p>
If changing from another style backrest and the Drum Brake is present, a new Drum Brake Cable Kit may be required

Item Number	Part Number	Description	Quantity	Retail
2a, 3	110587	Fixed Height Back Post 7.438" with Tube End Plugs <i>Used with 9" back height with omit push handles</i>	1	\$39.96
2b, 3	110588	Fixed Height Back Post 8.438" with Tube End Plugs <i>Used with 10" back height with omit push handles</i>	1	\$39.96
2c, 3	110589	Fixed Height Back Post 9.438" with Tube End Plugs <i>Used with 11" back height with omit push handles</i>	1	\$39.96
2d, 3	110590	Fixed Height Back Post 10.438" with Tube End Plugs <i>Used with 12" back height with omit push handles</i>	1	\$39.96
2e, 3	110591	Fixed Height Back Post 11.438" with Tube End Plugs <i>Used with 13" back height with omit push handles</i>	1	\$39.96

Item Number	Part Number	Description	Quantity	Retail
2f, 3	110592	Fixed Height Back Post 12.438" with Tube End Plugs <i>Used with 14" back height with omit push handles</i>	1	\$39.96
2g, 3	110593	Fixed Height Back Post 13.438" with Tube End Plugs <i>Used with 15" back height with omit push handles</i>	1	\$39.96
2h, 3	110594	Fixed Height Back Post 14.438" with Tube End Plugs <i>Used with 16" back height with omit push handles</i>	1	\$39.96
2i, 3	110595	Fixed Height Back Post 15.438" with Tube End Plugs <i>Used with 17" back height with omit push handles</i>	1	\$39.96
2j, 3	110596	Fixed Height Back Post 16.438" with Tube End Plugs <i>Used with 18" back height with omit push handles</i>	1	\$39.96
2k, 3, 9	110597	Fixed Height Back Post 5.438" with Tube End Plug and Sleeve <i>Used with 10" back height with standard or low profile push handles</i>	1	\$42.41
2l, 3, 9	110598	Fixed Height Back Post 6.438" with Tube End Plug and Sleeve <i>Used with 11" back height with standard or low profile push handles or 12" - 13" back heights with fold down push handles</i>	1	\$42.41
2a, 3, 9	110599	Fixed Height Back Post 7.438" with Tube End Plug and Sleeve <i>Used with 12" back height with standard or low profile push handles or 14" back height with fold down push handles</i>	1	\$42.41
2b, 3, 9	110600	Fixed Height Back Post 8.438" with Tube End Plug and Sleeve <i>Used with 13" back height with standard or low profile push handles or 15" back height with fold down push handles</i>	1	\$42.41
2c, 3, 9	110601	Fixed Height Back Post 9.438" with Tube End Plug and Sleeve <i>Used with 14" back height with standard or low profile push handles or 16" back height with fold down push handles</i>	1	\$42.41
2d, 3, 9	110602	Fixed Height Back Post 10.438" with Tube End Plug and Sleeve <i>Used with 17" back height with fold down push handles</i>	1	\$42.41
2e, 3, 9	110603	Fixed Height Back Post 11.438" with Tube End Plug and Sleeve <i>Used with 15" back height with standard or low profile push handles or 18" back height with fold down push handles</i>	1	\$42.41
2f, 3, 9	110604	Fixed Height Back Post 12.438" with Tube End Plug and Sleeve <i>Used with 16" back height with standard or low profile push handles</i>	1	\$42.41
2g, 3, 9	110605	Fixed Height Back Post 13.438" with Tube End Plug and Sleeve <i>Used with 17" back height with standard or low profile push handles</i>	1	\$42.41
2h, 3, 9	110606	Fixed Height Back Post 14.438" with Tube End Plug and Sleeve <i>Used with 18" back height with standard or low profile push handles</i>	1	\$42.41
3	100804	Tube End Plug 1"	4	\$1.58

Item Number	Part Number	Description	Quantity	Retail
4	000073	M5 x 24 Threaded Barrel	4	\$3.78
5	100669	M5 x 10 BTN SHCS BLKZC w/Patch	4	\$1.58
6a, 7a, 8	110584	Standard Push Handle Assembly, 10" Back Height <i>Used with 10" back height and drum brake is not present</i>	1	\$34.92
6b, 7a, 8	502602	Standard Push Handle Assembly, Short <i>Used with 11" - 18" back heights and drum brake is not present</i>	1	\$34.92
6e, 7a, 8	113221	Backtube w/ Push Handle Short with Backcane Extension, Drum Brake, Rigid <i>Used when drum brake is present. To view the backcane extension included in this kit, see Rogue XP / Klik / Little Wave XP Drum Brake page.</i>	1	\$52.29
6c, 7b, 8	110585	Low Profile Push Handle Assembly, 10" Back Height <i>Used with 10" back height</i>	1	\$36.22
6d, 7b, 8	107205	Low Profile Push Handle Assembly, Short <i>Used with 11" - 18" back heights</i>	1	\$36.22
7a	109432	Ergo Push Handle Grip 3/4" Standard with Instruction	1	\$7.14
7b	109430	Ergo Push Handle Grip 3/4" Low Profile with Instruction	1	\$7.14
8	100638	Tube End Plug 3/4"	2	\$2.10
9	000034	Sleeve 3/4" x 1 1/2"	2	\$3.66
10	100734	M5 x 16 PHMS ZC <i>Used to mount back upholstery and/or Fold Down Push Handle</i>	2	\$1.58
10, 11a, 12, 14	110586	Fixed Height Fold Down Push Handle Assembly, 12" Back Height <i>Used with 12" back height</i>	1	\$114.82
10, 11b, 12, 14	502603	Fixed Height Fold Down Push Handle Assembly, 13" - 18" Back Height <i>Used with 13" - 18" back heights</i>	1	\$114.82
13a	001748	Tension Strap Folding Handle 9-11 Wide	1	\$18.32
13b	001749	Tension Strap Folding Handle 12-13 Wide	1	\$19.52
13c	001750	Tension Strap Folding Handle 14-15 Wide	1	\$20.81
13d	001751	Tension Strap Folding Handle 16-17 Wide	1	\$21.98
13e	001752	Tension Strap Folding Handle 18-19 Wide	1	\$23.28
14	103524	M5 x 10 BTN SHCS BLK ZN	4	\$0.53