



OPEN FRAME.
OPEN ROAD.



The open frame of the aluminum Aero Z calls out to the open road. The Aero Z is the only TiFit, mono-tube aluminum frame on the

market. In a class by itself, the Aero Z is lightweight and fully adjustable and it won't break the bank!









## TECHNICAL SPECIFICATIONS.

Frame Features Aluminum Mono-Tube Rigid

TiFit

Full Adjustability

Frame Options Titanium Frame Upgrade

SuperLite Edition
Armed Forces Edition
Power Adaptable Frame
Reinforced Frame

Transport Weight\* 12.4 lbs.

Weight Limit 265 lbs.

Total Chair Width\* 23"

Seat Width 10" - 20"

Front Seat Width (Seat Taper) 3 Options

Seat Depth 10" - 20"

Front Seat Height 16" - 21"
Rear Seat Height 12.5" - 20"

12.5 - 20

Camber 0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°

Front Frame Angle 75°, 80°, 85° or 90°

Footrest Width Options Standard Taper or V-Front End

Center of Gravity Adjustment 6"

Front Wheels 23 Options
Rear Wheels 67 Options
Handrims 15 Options

Convenience Transit Tie-Down Option

Folding Angle Adjustable Seat Back Depth Adjustable Back/Frame

Aesthetics 5 Aluminum Frame Finishes

19 Frame Colors

6 Color Anodize Packages7 Ultrasuede Colors8 Spinergy Wheel Colors

14000

HCPCS Code K0005

Warranty Lifetime on Frame

Limited Lifetime on Components

\*Weights and widths based on 16" x 16" chair with  $0^{\circ}$  camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lbs.

## TIFIT. PERSONAL ERGONOMICS.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is measure. The optimized fit of a TiFit chair brings unparallele

made to measure. The optimized fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

## CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!



TiLite | The Ultimate Ride™ 800.545.2266 | tilite.com | © 2014 TiLite, LLC. TiLite wheelchairs are shown with available options.