


Catalyst E Transit



Notice: Chair configuration review required prior to parts being quoted or ordered. Not available with Heavy Duty Option

Item Number	Part Number	Description	Quantity
1-9	115998	Transit Bracket Retrofit Kit Catalyst E <i>Not available with Heavy Duty Option. This kit is a complete retrofit, containing a pair of front and rear transit brackets and hardware.</i>	1
1	100558	M6 Nyl Ins Domed Cap L/N Zc	4
2	100746	M6 Flat Washer Blk Zn	4
3	000056	Saddle 1 1/8"	10
5	100510	Eccentric Washer	4
6	100677	M6 x 60 CLS10.9 DIN931 HHCS ZC	4
7	100658	M6 Nylock Nut Blk Zn	6
9	101826	M6 x 50 HHCS Blk Zn	6