


# Catalyst E Tires - Inner Tubes



See Tire Component Chart for utilization.

Item Number	Part Number	Description	Quantity	Retail
1a	100912	Inner Tube 20" x 1 3/8"	1	\$12.61
1b	100588	Inner Tube 22" x 1 3/8"	1	\$12.61
1c	100589	Inner Tube 24" x 1 3/8"	1	\$12.61
1d	100590	Inner Tube 26" x 1 3/8"	1	\$12.61
1e	102844	Inner Tube 20" x 1"	1	\$12.61
1f	100609	Inner Tube 22" x 1"	1	\$12.61
1g	100580	Inner Tube 24" x 1"	1	\$12.61

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>	<b>Retail</b>
1h	101457	Inner Tube 25" x 1"	1	\$12.61
1i	100582	Inner Tube 26" x 1"	1	\$12.61

## Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part No.	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Marathon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke <sup>1</sup>	200529	101587	101588	101137																	
20" Ki Spoke <sup>1</sup>	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke <sup>1</sup>	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke <sup>1</sup>	200532	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke <sup>1</sup>	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke <sup>1</sup>	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
24" Superlight	101159	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	111818	100589	100596				102982	100589	100596					102920	100589	100596				
20" Spinergy Spox	See Spinergy Spox Page	100911	100912	101137	102957	102844	101137														
22" Spinergy Spox		100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy Spox		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy Spox					102959	101457	101452				101357	103177	101457	101452					100166	100208	101452
26" Spinergy Spox		100636	100590	100608	102960	100582	100597				100865		100582	100597							
22" Spinergy LX	See Spinergy LX Page	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinergy LX		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinergy LX					102959	101457	101452				101357	103177	101457	101452					100166	100208	101452
26" Spinergy LX		100636	100590	100608	102960	100582	100597				100865		100582	100597							
18" MAXX Spoke	109285	101587	101588	101137																	
20" MAXX Spoke	108244	100911	100912	101137	102957	102844	101137														
22" MAXX Spoke	105135	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" MAXX Spoke	105136	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" MAXX Spoke	107436				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" MAXX Spoke	107437	100636	100590	100608	102960	100582	100597				100865		100582	100597							

<sup>1</sup> Discontinued