


# Catalyst Seat Upholstery



| SCREW (ITEM 4) QTY PER CHAIR |           |
|------------------------------|-----------|
| SEAT DEPTH                   | SCREW QTY |
| 14" - 15"                    | 10        |
| 16" - 17"                    | 12        |
| 18" - 19"                    | 14        |
| 20"                          | 16        |

Seat upholstery width measures approximately 2.5" less than the chair size. Seat upholstery depth measures true to the chair size.

| Item Number | Part Number | Description              | Quantity | Retail  |
|-------------|-------------|--------------------------|----------|---------|
| 1a          | 100093      | Seat Upholstery Cat 1414 | 1        | \$94.55 |

| Item Number | Part Number | Description              | Quantity | Retail  |
|-------------|-------------|--------------------------|----------|---------|
| 1b          | 100094      | Seat Upholstery Cat 1415 | 1        | \$94.55 |
| 1c          | 100095      | Seat Upholstery Cat 1416 | 1        | \$94.55 |
| 1d          | 100096      | Seat Upholstery Cat 1417 | 1        | \$94.55 |
| 1e          | 100097      | Seat Upholstery Cat 1418 | 1        | \$94.55 |
| 1f          | 100896      | Seat Upholstery Cat 1419 | 1        | \$94.55 |
| 1g          | 100136      | Seat Upholstery Cat 1420 | 1        | \$94.55 |
| 1h          | 100098      | Seat Upholstery Cat 1514 | 1        | \$94.55 |
| 1i          | 100099      | Seat Upholstery Cat 1515 | 1        | \$94.55 |
| 1j          | 100100      | Seat Upholstery Cat 1516 | 1        | \$94.55 |
| 1k          | 100101      | Seat Upholstery Cat 1517 | 1        | \$94.55 |
| 1l          | 100102      | Seat Upholstery Cat 1518 | 1        | \$94.55 |
| 1m          | 100103      | Seat Upholstery Cat 1519 | 1        | \$94.55 |
| 1n          | 100104      | Seat Upholstery Cat 1520 | 1        | \$94.55 |
| 1o          | 100105      | Seat Upholstery Cat 1614 | 1        | \$94.55 |
| 1p          | 100106      | Seat Upholstery Cat 1615 | 1        | \$94.55 |
| 1q          | 100107      | Seat Upholstery Cat 1616 | 1        | \$94.55 |
| 1r          | 100108      | Seat Upholstery Cat 1617 | 1        | \$94.55 |
| 1s          | 100109      | Seat Upholstery Cat 1618 | 1        | \$94.55 |
| 1t          | 100110      | Seat Upholstery Cat 1619 | 1        | \$94.55 |
| 1u          | 100111      | Seat Upholstery Cat 1620 | 1        | \$94.55 |
| 1v          | 100112      | Seat Upholstery Cat 1714 | 1        | \$94.55 |
| 1w          | 100113      | Seat Upholstery Cat 1715 | 1        | \$94.55 |
| 1x          | 100114      | Seat Upholstery Cat 1716 | 1        | \$94.55 |
| 1y          | 100115      | Seat Upholstery Cat 1717 | 1        | \$94.55 |
| 1z          | 100116      | Seat Upholstery Cat 1718 | 1        | \$94.55 |
| 1aa         | 100117      | Seat Upholstery Cat 1719 | 1        | \$94.55 |
| 1ab         | 100118      | Seat Upholstery Cat 1720 | 1        | \$94.55 |
| 1ac         | 100119      | Seat Upholstery Cat 1814 | 1        | \$94.55 |

| Item Number | Part Number | Description              | Quantity | Retail   |
|-------------|-------------|--------------------------|----------|----------|
| 1ad         | 100120      | Seat Upholstery Cat 1815 | 1        | \$94.55  |
| 1ae         | 100121      | Seat Upholstery Cat 1816 | 1        | \$94.55  |
| 1af         | 100122      | Seat Upholstery Cat 1817 | 1        | \$94.55  |
| 1ag         | 100123      | Seat Upholstery Cat 1818 | 1        | \$94.55  |
| 1ah         | 100124      | Seat Upholstery Cat 1819 | 1        | \$94.55  |
| 1ai         | 100125      | Seat Upholstery Cat 1820 | 1        | \$94.55  |
| 1aj         | 101558      | Seat Upholstery Cat 1914 | 1        | \$94.55  |
| 1ak         | 100906      | Seat Upholstery Cat 1915 | 1        | \$94.55  |
| 1al         | 100126      | Seat Upholstery Cat 1916 | 1        | \$94.55  |
| 1am         | 100127      | Seat Upholstery Cat 1917 | 1        | \$94.55  |
| 1an         | 100128      | Seat Upholstery Cat 1918 | 1        | \$94.55  |
| 1ao         | 100129      | Seat Upholstery Cat 1919 | 1        | \$94.55  |
| 1ap         | 100130      | Seat Upholstery Cat 1920 | 1        | \$94.55  |
| 1aq         | 101559      | Seat Upholstery Cat 2014 | 1        | \$94.55  |
| 1ar         | 101557      | Seat Upholstery Cat 2015 | 1        | \$94.55  |
| 1as         | 100131      | Seat Upholstery Cat 2016 | 1        | \$94.55  |
| 1at         | 100132      | Seat Upholstery Cat 2017 | 1        | \$94.55  |
| 1au         | 100133      | Seat Upholstery Cat 2018 | 1        | \$94.55  |
| 1av         | 100134      | Seat Upholstery Cat 2019 | 1        | \$94.55  |
| 1aw         | 100135      | Seat Upholstery Cat 2020 | 1        | \$94.55  |
| 1ax         | 102526      | Seat Upholstery Cat 2114 | 1        | \$94.55  |
| 1ay         | 102527      | Seat Upholstery Cat 2115 | 1        | \$94.55  |
| 1az         | 101450      | Seat Upholstery Cat 2116 | 1        | \$94.55  |
| 1ba         | 101554      | Seat Upholstery Cat 2117 | 1        | \$94.55  |
| 1bb         | 100855      | Seat Upholstery Cat 2118 | 1        | \$105.06 |
| 1bc         | 100856      | Seat Upholstery Cat 2119 | 1        | \$105.06 |
| 1bd         | 100857      | Seat Upholstery Cat 2120 | 1        | \$105.06 |
| 1be         | 200352      | Seat Upholstery Cat 2214 | 1        | \$94.55  |

| Item Number | Part Number | Description                | Quantity | Retail   |
|-------------|-------------|----------------------------|----------|----------|
| 1bf         | 200353      | Seat Upholstery Cat 2215   | 1        | \$94.55  |
| 1bg         | 101556      | Seat Upholstery Cat 2216   | 1        | \$94.55  |
| 1bh         | 101555      | Seat Upholstery Cat 2217   | 1        | \$94.55  |
| 1bi         | 100247      | Seat Upholstery Cat 2218   | 1        | \$105.06 |
| 1bj         | 100858      | Seat Upholstery Cat 2219   | 1        | \$105.06 |
| 1bk         | 100261      | Seat Upholstery Cat 2220   | 1        | \$105.06 |
| 2a          | 100230      | Seat Rail 14" Deep CAT     | 2        | \$17.06  |
| 2b          | 100231      | Seat Rail 15" Deep CAT     | 2        | \$17.06  |
| 2c          | 100232      | Seat Rail 16" Deep CAT     | 2        | \$17.06  |
| 2d          | 100233      | Seat Rail 17" Deep CAT     | 2        | \$17.06  |
| 2e          | 100234      | Seat Rail 18" Deep CAT     | 2        | \$17.06  |
| 2f          | 100235      | Seat Rail 19" Deep CAT     | 2        | \$17.06  |
| 2g          | 100236      | Seat Rail 20" Deep CAT     | 2        | \$17.06  |
| 3a          | 000059      | Lifting Strap 14" Wide CAT | 1        | \$6.39   |
| 3b          | 000060      | Lifting Strap 15" Wide CAT | 1        | \$6.39   |
| 3c          | 000061      | Lifting Strap 16" Wide CAT | 1        | \$6.39   |
| 3d          | 000062      | Lifting Strap 17" Wide CAT | 1        | \$6.39   |
| 3e          | 000063      | Lifting Strap 18" Wide CAT | 1        | \$6.39   |
| 3f          | 000064      | Lifting Strap 19" Wide CAT | 1        | \$6.39   |
| 3g          | 000058      | Lifting Strap 20" Wide CAT | 1        | \$6.39   |
| 3h          | 100248      | Lifting Strap 21" Wide CAT | 1        | \$6.40   |
| 3i          | 100875      | Lifting Strap 22" Wide CAT | 1        | \$6.40   |
| 4           | 100734      | M5 x 16 PHMS ZC            | 16       | \$1.58   |