


Catalyst Seat Upholstery



| SCREW (ITEM 4) QTY PER CHAIR | |
|------------------------------|-----------|
| SEAT DEPTH | SCREW QTY |
| 14" - 15" | 10 |
| 16" - 17" | 12 |
| 18" - 19" | 14 |
| 20" | 16 |

Seat upholstery width measures approximately 2.5" less than the chair size. Seat upholstery depth measures true to the chair size. If adding seat upholstery, rails and screws are needed.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--------------------------|----------|
| 1a | 100093 | Seat Upholstery Cat 1414 | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|--------------------|--------------------------|-----------------|
| 1b | 100094 | Seat Upholstery Cat 1415 | 1 |
| 1c | 100095 | Seat Upholstery Cat 1416 | 1 |
| 1d | 100096 | Seat Upholstery Cat 1417 | 1 |
| 1e | 100097 | Seat Upholstery Cat 1418 | 1 |
| 1f | 100896 | Seat Upholstery Cat 1419 | 1 |
| 1g | 100136 | Seat Upholstery Cat 1420 | 1 |
| 1h | 100098 | Seat Upholstery Cat 1514 | 1 |
| 1i | 100099 | Seat Upholstery Cat 1515 | 1 |
| 1j | 100100 | Seat Upholstery Cat 1516 | 1 |
| 1k | 100101 | Seat Upholstery Cat 1517 | 1 |
| 1l | 100102 | Seat Upholstery Cat 1518 | 1 |
| 1m | 100103 | Seat Upholstery Cat 1519 | 1 |
| 1n | 100104 | Seat Upholstery Cat 1520 | 1 |
| 1o | 100105 | Seat Upholstery Cat 1614 | 1 |
| 1p | 100106 | Seat Upholstery Cat 1615 | 1 |
| 1q | 100107 | Seat Upholstery Cat 1616 | 1 |
| 1r | 100108 | Seat Upholstery Cat 1617 | 1 |
| 1s | 100109 | Seat Upholstery Cat 1618 | 1 |
| 1t | 100110 | Seat Upholstery Cat 1619 | 1 |
| 1u | 100111 | Seat Upholstery Cat 1620 | 1 |
| 1v | 100112 | Seat Upholstery Cat 1714 | 1 |
| 1w | 100113 | Seat Upholstery Cat 1715 | 1 |
| 1x | 100114 | Seat Upholstery Cat 1716 | 1 |
| 1y | 100115 | Seat Upholstery Cat 1717 | 1 |
| 1z | 100116 | Seat Upholstery Cat 1718 | 1 |
| 1aa | 100117 | Seat Upholstery Cat 1719 | 1 |
| 1ab | 100118 | Seat Upholstery Cat 1720 | 1 |
| 1ac | 100119 | Seat Upholstery Cat 1814 | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|--------------------|--------------------------|-----------------|
| 1ad | 100120 | Seat Upholstery Cat 1815 | 1 |
| 1ae | 100121 | Seat Upholstery Cat 1816 | 1 |
| 1af | 100122 | Seat Upholstery Cat 1817 | 1 |
| 1ag | 100123 | Seat Upholstery Cat 1818 | 1 |
| 1ah | 100124 | Seat Upholstery Cat 1819 | 1 |
| 1ai | 100125 | Seat Upholstery Cat 1820 | 1 |
| 1aj | 101558 | Seat Upholstery Cat 1914 | 1 |
| 1ak | 100906 | Seat Upholstery Cat 1915 | 1 |
| 1al | 100126 | Seat Upholstery Cat 1916 | 1 |
| 1am | 100127 | Seat Upholstery Cat 1917 | 1 |
| 1an | 100128 | Seat Upholstery Cat 1918 | 1 |
| 1ao | 100129 | Seat Upholstery Cat 1919 | 1 |
| 1ap | 100130 | Seat Upholstery Cat 1920 | 1 |
| 1aq | 101559 | Seat Upholstery Cat 2014 | 1 |
| 1ar | 101557 | Seat Upholstery Cat 2015 | 1 |
| 1as | 100131 | Seat Upholstery Cat 2016 | 1 |
| 1at | 100132 | Seat Upholstery Cat 2017 | 1 |
| 1au | 100133 | Seat Upholstery Cat 2018 | 1 |
| 1av | 100134 | Seat Upholstery Cat 2019 | 1 |
| 1aw | 100135 | Seat Upholstery Cat 2020 | 1 |
| 1ax | 102526 | Seat Upholstery Cat 2114 | 1 |
| 1ay | 102527 | Seat Upholstery Cat 2115 | 1 |
| 1az | 101450 | Seat Upholstery Cat 2116 | 1 |
| 1ba | 101554 | Seat Upholstery Cat 2117 | 1 |
| 1bb | 100855 | Seat Upholstery Cat 2118 | 1 |
| 1bc | 100856 | Seat Upholstery Cat 2119 | 1 |
| 1bd | 100857 | Seat Upholstery Cat 2120 | 1 |
| 1be | 200352 | Seat Upholstery Cat 2214 | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|----------------------------|----------|
| 1bf | 200353 | Seat Upholstery Cat 2215 | 1 |
| 1bg | 101556 | Seat Upholstery Cat 2216 | 1 |
| 1bh | 101555 | Seat Upholstery Cat 2217 | 1 |
| 1bi | 100247 | Seat Upholstery Cat 2218 | 1 |
| 1bj | 100858 | Seat Upholstery Cat 2219 | 1 |
| 1bk | 100261 | Seat Upholstery Cat 2220 | 1 |
| 2a | 100230 | Seat Rail 14" Deep CAT | 2 |
| 2b | 100231 | Seat Rail 15" Deep CAT | 2 |
| 2c | 100232 | Seat Rail 16" Deep CAT | 2 |
| 2d | 100233 | Seat Rail 17" Deep CAT | 2 |
| 2e | 100234 | Seat Rail 18" Deep CAT | 2 |
| 2f | 100235 | Seat Rail 19" Deep CAT | 2 |
| 2g | 100236 | Seat Rail 20" Deep CAT | 2 |
| 3a | 000059 | Lifting Strap 14" Wide CAT | 1 |
| 3b | 000060 | Lifting Strap 15" Wide CAT | 1 |
| 3c | 000061 | Lifting Strap 16" Wide CAT | 1 |
| 3d | 000062 | Lifting Strap 17" Wide CAT | 1 |
| 3e | 000063 | Lifting Strap 18" Wide CAT | 1 |
| 3f | 000064 | Lifting Strap 19" Wide CAT | 1 |
| 3g | 000058 | Lifting Strap 20" Wide CAT | 1 |
| 3h | 100248 | Lifting Strap 21" Wide CAT | 1 |
| 3i | 100875 | Lifting Strap 22" Wide CAT | 1 |
| 4 | 100734 | M5 x 16 PHMS ZC | 16 |