

First Class School Chair

Designed to assure good seating posture, the First Class school chair is comfortable, attractive and adjustable. It is ideal for physically challenged children with mild to moderate involvement for use in the classroom, clinic or home. Its multiple adjustments and accessories meet the needs of each child as they grow and develop.

So easy to adjust!! Simplifies use for multiple children in schools and clinics.

Basic model includes pelvic strap and armrests.

FC 2000 Small, 1/cs

FC 4000 Large, 1/cs

FEATURES:

- Seat - adjustable in height and depth with a 15 degree anterior or posterior tilt.
- Armrest - adjustable in height and width.
- Pelvic Belt - adjustable seat belt helps maintain proper hip position.
- Two year limited warranty.



ACCESSORIES:

Support Kit

FC 8025 Support Kit for FC 2000 & FC 4000, 1/bx

Includes:

- Harness, trunk harness in an H configuration for anterior chest support
- Abductor, adjustable in depth.
- Lateral Supports, individually adjustable in height and width.

Clear Tray

FC 2026 Clear Tray for FC 2000, 1/bx

FC 4026 Clear Tray for FC 4000, 1/bx

Clear Horizontal Positional Device.

Tray

FC 2024 Tray for FC 2000, 1/bx

FC 4024 Tray for FC 4000, 1/bx

Horizontal Positional Device (Tray).

Footrest

FC 2028 Footrest for FC 2000, 1/bx

FC 4028 Footrest for FC 4000, 1/bx

Two-piece flip up footplate with heel loop.



Headrest

FC 8000 Adjustable Headrest for
FC 2000 & FC 4000, 1/bx

Height, depth and angle adjustable
headrest rotates and can be placed off center.

Headrest Wedges

WD 8000 Wedges, 1 pr/bx

Wedges can be placed anywhere on the
headrest for additional positioning.

Hip Guides

FC 8027 Hip Guides for FC 2000 & FC 4000, 1pr/bx

Hip Guides are width adjustable.

Mobility Base

FC 2029 Mobility Base for FC 2000, 1/bx

FC 4029 Mobility Base for FC 4000, 1/bx

Four wheeled mobility base safely transports child.

Leg Extensions

FC 8050 Leg extensions for FC 2000 & FC 4000, 1pr/bx

Set of four leg extensions for custom height.
Specify height range required.

Anti-tipper

FC 2030 Anti-tipper for FC 2000, 1 pr/bx

FC 4030 Anti-tipper for FC 4000, 1 pr/bx

Anti-tipper provides additional stability.



positioning
seating accessories

Wenzel[®] Re/hab

To order call toll free:
New York: 877.224.0946 • California: 800.371.2266

ITEM #	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT TO FLOOR	SEAT HEIGHT TO FOOTPLATE	BACK HEIGHT	LENGTH OF ARMREST	SEAT TO ARMREST HEIGHT	WEIGHT CAPACITY
FC 2000	12"	9" - 14"	9" - 13"	6" - 12"	15"	10"	1.5" - 9.5"	100 lbs.
FC 4000	16"	13" - 17"	14" - 23"	6" - 22"	17"	10"	1.5" - 11.5"	200 lbs.

ITEM #	SPACE BETWEEN ARMRESTS	HIP GUIDES	SPACE BETWEEN HIP GUIDES	LATERALS LATERAL HEIGHT	SEAT TO LATERALS	SPACE BETWEEN	HEADREST
FC 2000	9.5" - 14.5"	4" H X 9" D	8" - 13"	4" H X 5" D	1.5" - 6"	6" - 13"	ADDS 4.5" - 8" HEIGHT
FC 4000	13.5" - 18.5"	4" H X 9" D	12" - 17"	4" H X 5" D	1.5" - 9"	10" - 17"	ADDS 4.5" - 8" HEIGHT