



Raz-AT

MSRP – USA – USD
Effective Jan 1, 2022

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AT – Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height-adjustable stainless-steel frame, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan.
User weight capacity 400 lbs. See our [Evaluation and Medical Justification Form](#).

1. Chair Codes by size – Must Choose One

Z300S18	Raz-AT w/ 18"W seat and base frame – Standard size, fits over most toilets ¹	\$2,995
Z300S__*	Non-Standard Width Raz-AT – Must add upcharge for custom width: see below * add custom width to complete code (16" / 20" / 22" / 24") – Adjusta-Back must be selected for all custom widths [Sec.5]	\$2,995
Z3SW16	Upcharge for custom 16"W seat and base frame ² – will not fit over most toilets	\$595
Z3SW20	Upcharge for custom 20"W seat and base frame ³	\$595
Z3SW22	Upcharge for custom 22"W seat and base frame ⁴	\$695
Z3SW24	Upcharge for custom 24"W seat and base frame ⁴	\$795

¹ Fits all 18"W seats

² Visco-Foam Interface seats can be modified to fit.

³ 18"W seats can be selected but must also select Seat Support Adapter Kit Z4SSKT [Sec.2]

⁴ Z4V23N Visco Foam Interface HD seat can be modified to fit. Seat Support Kit (Z4SSKT) required for 24"W seat frame.

Custom Seat / Base Frame Length and Seat Height

Z3SFL3	Upcharge for custom 3" Longer Seat Frame: required for seats ≥ 21"D	\$595
Z3BFL3	Upcharge for custom 3" Longer Base Frame: required for seats ≥ 22"D (typically used in conjunction with Z3SFL3)	\$595
Z3SL5	Upcharge for custom Super-Low Frame (5" lower) – height adjustment range (molded seat): 16"-20" ¹	\$1,595

¹ To lower height adjustment range by 1" (15"-19"), select 4" caster upgrade option Z64T (will not fit over most toilets)

2. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
Z405	18"W x 16"D – Front Access	Standard	Z436BRI	18"W x 16"D – Front Bridged	\$425
Z410	18"W x 16"D – Front Bridged	N/C	Z436FRO	18"W x 16"D – Front Access	\$425
Z415	18"W x 19"D – Front Access	\$135	Z436LTO	18"W x 16"D – Left Access	\$425
Z420	18"W x 19"D – Front Bridged	\$135	Z436RTO	18"W x 16"D – Right Access	\$425
Specialty Seats			Z436RRO	18"W x 16"D – Rear Access	\$425
Z4SOM	Shower-Only Seat – Molded 18"W x 15.25"-20.25"D	\$395	Z439BRI	18"W x 19"D – Front Bridged	\$495
Z4EJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$295	Z439FRO	18"W x 19"D – Front Access	\$495
Z4EJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	Z439LTO	18"W x 19"D – Left Access	\$495
Z4V23N	23"W x 18"D HD – Visco Foam – Front Access Notched at rear to maintain fore/aft adjustment range (IPAS) – 22"W frame	\$695	Z439RTO	18"W x 19"D – Right Access	\$495
			Z439RRO	18"W x 19"D – Rear Access	\$495

3. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

Z4CS	Custom Seat (enter specs below) Select if both width and depth are ≤ 21"	\$795
Z4CEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$995
Z4CSHD	Heavy Duty Custom Seats (enter specs below) Select if width and/or depth is ≥ 22"	\$1,095

Custom Specifications (for above custom seats)

Seat Frame Width:	-- 16" / 18" / 20" / 22" / 24"
Custom Seat Width:	-- 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than seat frame width)
Custom Seat Depth:	-- 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Longer Frame (Z3SFL3)]
Access Opening:	-- FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
Access Opening Width:	-- 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
Aperture Width:	-- 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

4. Seat Options

Z4SSKT	Seat Support Kit for Raz-AT/CAT (to support seats that are narrower than seat frames)	\$145	ZSCS	Splash Curtain for Bridged Front Molded and VFI Seats – helps direct urine into pan or toilet. Note: raises seat 1/2"	\$295
Z4GSB_*	Custom Seat Bracket Spacers – for adding dump angle (pr) Specify additional height at front of seat to complete code (0.5"/1.0"/1.5")	\$295	ZSCC	Splash Curtain for Bridged Front Custom Seats – helps direct urine into pan or toilet. Note: raises seat 1/2"	\$345

4. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded, vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (85°/90°/95°) with cane-style posts. Back height is measured from top of molded seat to top of back upholstery. Hip Adductors [Z7HP1/2] are not compatible with Fixed Back.

Z50818	Fixed Back – 18"W x 18"H ^{1,2}	Standard
Z5BU1818	Tension-Adjustable Back Upholstery for Fixed Back (Z508)	\$185

5. Adjusta-Back Frames




Includes height- and angle-adjustable arms with Armrest Lock Kit (Z710), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

Standard Width Adjust-Backs (18"W)			Non-Standard Width Adjusta-Backs		
Z50918	Adjusta-Back – 15"H	\$225	<i>* Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22"/24") Note: 18" is the standard width.</i>		
Z51018	Adjusta-Back – 17"H	\$225	Z509__*	Adjusta-Back – 15"H non-std width	\$225
Z51518	Adjusta-Back – 20"H	\$295	Z510__*	Adjusta-Back – 17"H non-std width	\$225
Z52018	Adjusta-Back – 24"H	\$345	Z515__*	Adjusta-Back – 20"H non-std width	\$295
Z5DTL	Dual Triggers – both on left side ¹	N/C	Z520__*	Adjusta-Back – 24"H non-std width	\$345
Z5DTR	Dual Triggers – both on right side ¹	N/C	Z515ATR__*	Reclining Adjusta-Back ² 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$895
Z515ATR18	Reclining Adjusta-Back ² 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$895	Z520ATR__*	Reclining Adjusta-Back ² 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$995
Z520ATR18	Reclining Adjusta-Back ² 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	\$995			

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths.

² On reclining backs (Z515ATR/Z520ATR), Dual Triggers for recline are mounted on the left side and for tilt on the right side.

6. Custom Back Frame Widths Different From Seat Frames – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)

Z5NB____*	Custom Back – 2" Narrower than Seat Frame	\$345	Seat Frame Width	16" ¹	18" ²	20" ³	22" ⁴	24" ⁵	
			Back Frame Width	14"	16"	18"	20"	22"	
<i>* Enter seat and back frame widths to complete code. eg. Z5NB2018 specifies a 20"W seat frame with an 18"W back frame</i>									
Z5NB____*	Custom Back – 4" Narrower than Seat Frame – requires custom seat	\$445	Seat Frame Width	N/A	18" ¹	20" ³	22" ⁴	24" ⁵	
			Back Frame Width	N/A	14"	16"	18"	20"	
<i>* Enter seat and back frame widths to complete code. eg. Z5NB2016 specifies a 20"W seat frame with a 16"W back frame</i>									
Z5WB____*	Custom Back – 2" Wider than Seat Frame	\$345	Seat Frame Width	16"	18"	20" ³	22" ⁶	N/A	
			Back Frame Width	18"	20"	22"	24"	N/A	
<i>* Enter seat and back frame widths to complete code. eg. Z5WB1820 specifies an 18"W seat frame with a 20"W back frame</i>									
<i>Consider also using Arm Spacer Kit (Z718/Z719) as another solution for providing more width between the arms.</i>									

¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats, requires Seat Support Kit Z4SSKT [Sec.2]

⁴ Custom Seat recommended. Also compatible with modified Z4V23N (18" seat depth with limited adjustability).

⁵ Custom Seat recommended. Compatible with modified Z4V23N (18" seat depth with limited adjustability, requires Seat Support Kit (Z4SSKT) [Sec.2] and modified Z4V26N (18" seat depth with limited adjustability).

⁶ Custom Seat recommended. Also compatible with Z4V23N (18" seat depth).

7. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.

Z5S14	14"W Symphony Back ¹ – recessed solid back	\$395	Z749	Harmony Back ² (contoured solid back – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back frames only	\$245
Z5S16	16"W Symphony Back ¹ – recessed solid back	\$395			
Z5S18	18"W Symphony Back ¹ – recessed solid back	\$395			
Z5S20	20"W Symphony Back ¹ – recessed solid back	\$495			
Z5S22	22"W Symphony Back ¹ – recessed solid back	\$595			
Z5S24	24"W Symphony Back ¹ – recessed solid back	\$695			
Z5S26	26"W Symphony Back ¹ – recessed solid back	\$795			

¹ Not compatible with 15"H back – all Symphony Back panels are 13"H
² Not recommended for reclining backs.

8. Arm Supports and Accessories

Z5BAR13	Flip-Up Arm Supports – 13"L pad (pr)	Std	Z7PAML	Pivoting Arm Mount – Left (for Z7FL/Z7TR)	\$165
Z710	Arm Support Lock Kit (for Z505/Z506) (pr)	\$75	Z7PAMR	Pivoting Arm Mount – Right (for Z7FL/Z7TR)	\$165
Z52PAT	2-Point Flip-up Arm Supports ¹ (pr) (standard with reclining backs)	\$345	Z718	Arm Support Spacer Kit – 1/2" per side (pr)(for Z509/Z510/Z515/Z520)	\$60
Z7FLL	Long Molded Flat Armpad ² – Left (14"L x 3.25"W)	\$95	Z719	Arm Support Spacer Kit – 1" per side (pr)(for Z509/Z510/Z515/Z520)	\$115
Z7FLR	Long Molded Flat Armpad ² – Right (14"L x 3.25"W)	\$95	Z7PHG	Pivoting Hand Grips (pr)	\$395
Z7TRL	Molded Arm Trough ² – Left (18.5"L x 4.5"W)	\$110	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (14"/16"/18") ^{3,4}	\$195
Z7TRR	Molded Arm Trough ² – Right (18.5"L x 4.5"W)	\$110	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (20"/22"/24") ^{3,4}	\$225

¹ For chairs with same back width and frame width only
² May require Z7PAM for arm support to flip up and/or for lateral supports to swing away
³ Not available with Z7FL / Z7TR
⁴ Special length required with Z718 (+1") and Z719 (+2").

9. Foot and Leg Supports – Standard Foot Supports (Z705N) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support. e.g. By selecting Elevating Leg Support – Left (Z7ELRL), the chair will be supplied with one Foot Support (Z705N) on the right side and a Z7ELRL on the left side.

Z705N	Foot Supports 16"-20.5" (pr)	Std	Z7ELRL	Elevating Leg Support ³ – Left (includes Z7ADM)	\$475
Z705N16	Foot Supports for 16"W frames 16"-20.5" (pr)	N/C	Z7ELRR	Elevating Leg Support ³ – Right (includes Z7ADM)	\$475
Z7FET2	Footrest Extension Tubes – 2" longer ¹ (pr)	\$50	Z7ELROFL	Forward Offset Elevating Leg Support ⁴ – Left (includes Z7ADM)	\$525
Z7MFX8	MFX8 Foot Supports 8"-12" (pr)	\$75	Z7ELROFR	Forward Offset Elevating Leg Support ⁴ – Right (includes Z7ADM)	\$525
Z7MFX12	MFX12 Foot Supports 12"-16.5" (pr)	\$175	Z7RLSL	Universal Residual Limb Support – Left	\$395
Z7ADM	Angle/Depth Adj Molded Footplates ² (pr)	\$190	Z7RLSR	Universal Residual Limb Support – Right	\$395
Z7FET2A	Footrest Extension Tubes for Z7ADM (pr) (2" longer)	\$50			
Z7ADABS	Angle/Depth Adj ABS Footplates (pr) For 16"W frames	\$190			
Z7VS18	V-Style Foot Support (1-piece) – for 18"W base frames	\$195			
Z7VS_*	V-Style Foot Support (1-piece) * Enter frame width to complete code (16" / 20" / 22" / 24")	\$295			
Z7FBT	Flip-Back Foot Support (1-piece) with calf pads	\$395			
¹ Not compatible with Z7MFX8 / Z7ADM / Z7VS / Z7FBT ² Not compatible with Z300S16 or Z7MFX8			³ Z7ELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 3" longer seat frames (Z3SFL3) ⁴ Z7ELROFL/R compatible with frame widths ≥ 18" for adjusted seat depths 19"-20" on standard frame or 22"-23" w/ 3" longer seat frames (Z3SFL3)		

10. Foot and Leg Support Accessories

Z7FSLK	Foot Support Lock (pr)	\$95	Z708ADJ	Adjustable Heel Loops (pr) – not available for Z7VS / Z7FBT	\$60
Z7CS_*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter seat frame width to complete code (16"/18"/20"/22"/24")	\$95	Z7NCS	Neoprene Covers for Standard Footplates (pr)	\$90
Z7CSIC_*	Calf Strap, Snap-On – Infection Control * Enter seat frame width to complete code (16" / 18" / 20")	\$155	Z7NC16	Neoprene Covers for narrow molded footplates (pr)	\$90
Z7CSIC22	Snap-On Infection Control Calf Strap – 22"W frames	\$175	Z7NCA	Neoprene Covers for Angle/Depth Adj Molded Footplates (pr)	\$90
Z737CPL	Calf Panel – Bodypoint ¹ Fits frame widths 16" / 18" / 20"	\$175	Z7NCABS	Neoprene Covers for ABS Adj Angle/Depth Footplates Z7ADABS (pr)	\$90
Z706	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$90	Z7OFLM	Lateral Offset Foot Support Receiver – moves foot support 2" laterally or medially Quantity (1 or 2) __	\$160
Z712	H-Strap – for 22" / 24"W frames – hook and loop attachment	\$110			
¹ Not compatible with Z7MFX8 / Z7FBT / Z7VS					

11. Positioning Pads

Z7HSTS	Head Support Pad – Small (10" x 4")	Std	Z7LGT1L	Leg Adductor w/ Small Pad ¹ – Left	\$320
Z7HSTL	Head Support Pad – Large (12" x 6")	\$50	Z7LGT1R	Leg Adductor w/ Small Pad ¹ – Right	\$320
Z7OHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$95	Z7LGT2L	Leg Adductor w/ Large Pad ¹ – Left	\$320
Z7LTS1L	Lateral Thoracic Support w/ Small Pad ^{1,2} – Left	\$320	Z7LGT2R	Leg Adductor w/ Large Pad ¹ – Right	\$320
Z7LTS1R	Lateral Thoracic Support w/ Small Pad ^{1,2} – Right	\$320	Z7HP1L	Hip Adductor w/ Small Pad ^{1,3} – Left	\$320
Z7LTS2L	Lateral Thoracic Support w/ Large Pad ^{1,2} – Left	\$320	Z7HP1R	Hip Adductor w/ Small Pad ^{1,3} – Right	\$320
Z7LTS2R	Lateral Thoracic Support w/ Large Pad ^{1,2} – Right	\$320	Z7HP2L	Hip Adductor w/ Large Pad ^{1,3} – Left	\$320
Z742L	Lateral/Adductor Extension Plate (5.5" x 3.5") – Left	\$25	Z7HP2R	Hip Adductor w/ Large Pad ^{1,3} – Right	\$320
Z742R	Lateral/Adductor Extension Plate (5.5" x 3.5") – Right	\$25	Z7MTST1	Medial Thigh Support – Small	\$445
			Z7MTST2	Medial Thigh Support – Large	\$445

¹ Adjustable Pad Sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3".
² Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For additional medial offset, specify Z742 extension plate.
³ Fits Adjusta-Back only

12. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBS	Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$35	ZCBS	Chest Belt ² – Small (back widths 14"-16", 1-piece w/ hook and loop closure)	\$35
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$35	ZCBM	Chest Belt ² – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)	\$35
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$45	ZCBL	Chest Belt ² – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	\$45
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$55	ZCBXL	Chest Belt – XLarge ² (back frame widths 22"-30", 1-piece w/ hook and loop closure)	\$55
ZPBICS	Infection Control Pelvic Belt ¹ – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$75	ZCBICS	Infection Control Chest Belt ^{1,2} – Small (back widths 14"-16", 1-piece w/ side-release buckle)	\$75
ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$85	ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	\$75
ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$95	ZCBICL	Infection Control Chest Belt ^{1,2} – Large (back widths 18"-22", 1-piece w/ side-release buckle)	\$85
ZPBICXL	Infection Control Pelvic Belt ¹ – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$105	ZCBICXL	Infection Control Chest Belt ^{1,2} – XLarge (back frame widths 22"-30", 1-piece w/ side-release buckle)	\$95
ZPBBPS	Pelvic Belt – Bodypoint Aeromesh – Small	\$245	ZCBBPS	Chest Belt – Bodypoint Aeromesh – Small	\$255
ZPBBP _ *	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	\$195	ZCBBP _ *	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	\$205
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$165			

Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPS	14"/16"W seat frame	ZPBBPM	16"/18"W seat frame	ZPBBPL	18"-22"W seat frame	ZPBBPXL	22"-26"W seat frame
	Chest belt	ZCBBPS	14"/16"W back frame	ZCBBPM	16"/18"W back frame	ZCBBPL	18"/20"W back frame	ZCBBPXL	22"-26"W back frame

¹ Non-absorbent belts – recommended for infection-control environments

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

13. Other

Z675	Directional Caster Lock – fits on one rear caster (facilitates steering)	\$145	ZIHW	Whizard – winged urine deflector ²	\$67
Z66T4	Rear anti-Tippers for 4" casters (pr)	\$345	ZBUDBH_ _ *	Bard Urine Drainage Bag Holder * Enter size to complete code (16"/18"/20"/22"/24")	\$155
Z66T5	Rear anti-Tippers for 5" casters (pr)	\$345			
Z66T6	Rear anti-Tippers for 6" casters (pr)	\$345	ZSPKT22	AT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	\$1395
Z64T	4" Caster Upgrade (4x)	\$280	ZSPKT24	AT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	\$1395
Z66T	6" Caster Upgrade (4x)	\$350			
ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$200			
ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$215			
ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$230			
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$245			
ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ___	N/C		¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Reduces toilet bowl clearance height. ² Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	

14. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907	Omit Foot Supports - if no other Foot / Leg Support Option is selected	-\$45	Z912	Omit Complete Head Support Assembly	-\$50

Some assembly of Raz products is required.

15. Notes:

--

Specifications

Raz-AT (Attendant Tilt)

Weight Capacity	400 lbs
Overall Width	Seat width plus 4.5" with Standard Arm Pads
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	43"
- at 40°	54"
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat)	
- with pan brackets	17"-21" (1" increments)
- without pan brackets	17.75"-21.75" (1" increments)
Toilet Clearance Width	15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25"
Arm Support Height - above molded seat (1/2" increments)	
- 18"H Fixed Back - standard	8.5"
- 15"H Adjusta-Back	5.5"-8.5" (1/2" increments)
- 17"H Adjusta-Back	6.5"-9.5" (1/2" increments)
- 20"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 20"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	21"-25" in 1" increments
Product Weight (with 18"W seat)	61 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)