



ICON™ MID BACK



Icon™ Mid Back Size (in.)	12	14	15	16	17	18	20	22	24
A. Back Height (in.) ¹	9.7	11.2	12	12.7	13.5	14.25	15.75	15.75	15.75
B. Bottom of back to beginning of vertical slot (in.)	1.3	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5
C. Bottom of back to top of vertical slot (in.)	5.3	6.5	7.0	7.5	8.0	8.5	8.5	8.5	8.5
D. Range of Icon™ hardware adjustment in vertical slot (in.)	0.5	1.5	2.0	2.5	3.0	3.5	3.5	3.5	3.5
E. Outside width of shell (in.)	9.5	11.4	12.4	13.4	14.4	15.4	17.4	19.4	21.4
F. Inside width of support surface (in.)	8.7	10.6	11.6	12.6	13.6	14.6	16.6	18.6	20.6
Will fit wheelchair widths (in.)	12-13	13-15	14-16	15-17	16-18	17-19	19-21	21-23	23-25

1. Shell only. Complete back system height may vary by 1.5 in. due to the cushion rollover.