


Discontinued Catalyst Half Folding Back Tension Adjustable Back Upholstery



BACK HEIGHT STRAIGHT BACK POSTS	TOTAL STRAP QTY ITEM # 3	D-RING QTY ITEM #4
12.5" - 13.5"	4	4
14.5" - 15.5"	5	6
16.5" - 17.5"	6	8
18.5" - 19.5"	7	10
BACK HEIGHT 8" BEND BACK POST	TOTAL STRAP QTY ITEM # 3	D-RING QTY ITEM #4
13.5" - 14.5"	5	6
15.5" - 17.5"	6	8
18.8" - 19.5"	7	10

Item Number	Part Number	Description	Quantity	Retail
3a	003483	Back Upholstery Strap 7.3"	2	\$10.51
3b	003484	Back Upholstery Strap 9.3"	2	\$10.51
3c	003485	Back Upholstery Strap 11.3"	2	\$10.51
3d	003486	Back Upholstery Strap 13.3"	2	\$10.51
3e	003487	Back Upholstery Strap 15.3"	2	\$12.61
3f	002532	Back Upholstery Strap 17.3"	2	\$12.61
3g	002533	Back Upholstery Strap 19.3"	2	\$12.61

Item Number	Part Number	Description	Quantity	Retail
3h	002534	Back Upholstery Strap 21.3"	2	\$12.61
3i	002535	Back Upholstery Strap 23.3"	2	\$12.61
3j	002536	Back Upholstery Strap 25.3"	2	\$15.76
3k	002537	Back Upholstery Strap 27.3"	2	\$15.76
4	100973	Back Upholstery D-Ring	4	\$10.51
5a	002498	Backrest Cover 13" Rigid / 14" Cat x 11-14"	1	\$213.08
5b	002499	Backrest Cover 13" Rigid / 14" Cat x 13-16"	1	\$213.08
5c	002500	Backrest Cover 13" Rigid - 14" CAT/LFT x 15-18"	1	\$213.08
5d	003043	Backrest Cover 13" Rigid - 14" CAT/LFT x 17-20"	1	\$234.09
5e	002502	Backrest Cover 14" Rigid / 15" Cat x 11-14"	1	\$223.58
5f	002503	Backrest Cover 14" Rigid / 15" Cat x 13-16"	1	\$223.58
5g	002504	Backrest Cover 14" Rigid - 15" CAT/LFT x 15-18"	1	\$223.58
5h	003044	Backrest Cover 14" Rigid - 15" CAT/LFT x 17-20"	1	\$234.09
5i	002506	Backrest Cover 15" Rigid / 16" Cat x 11-14"	1	\$223.58
5j	002507	Backrest Cover 15" Rigid / 16" Cat x 13-16"	1	\$223.58
5k	002508	Backrest Cover 15" Rigid - 16" CAT/LFT x 15-18"	1	\$223.58
5l	003045	Backrest Cover 15" Rigid - 16" CAT/LFT x 17-20"	1	\$234.09
5m	002510	Backrest Cover 16" Rigid / 17" Cat x 11-14"	1	\$223.58
5n	002511	Backrest Cover 16" Rigid / 17" Cat x 13-16"	1	\$223.58
5o	002512	Backrest Cover 16" Rigid - 17" CAT/LFT x 15-18"	1	\$223.58
5p	003046	Backrest Cover 16" Rigid - 17" CAT/LFT x 17-20"	1	\$234.09
5q	003049	Backrest Cover 17" Rigid / 18" Cat x 11-14"	1	\$234.09
5r	003050	Backrest Cover 17" Rigid / 18" Cat x 13-16"	1	\$234.09
5s	003051	Backrest Cover 17" Rigid - 18" CAT/LFT x 15-18"	1	\$234.09
5t	003052	Backrest Cover 17" Rigid - 18" CAT/LFT x 17-20"	1	\$234.09
5u	003055	Backrest Cover 18" Rigid / 19" Cat x 11-14"	1	\$234.09
5v	003056	Backrest Cover 18" Rigid / 19" Cat x 13-16"	1	\$234.09
5w	003057	Backrest Cover 18" Rigid - 19" CAT/LFT x 15-18"	1	\$234.09

Item Number	Part Number	Description	Quantity	Retail
5x	003058	Backrest Cover 18" Rigid - 19" CAT/LFT x 17-20"	1	\$234.09
5y	003061	Backrest Cover 19" Rigid / 20" Cat x 11-14"	1	\$244.60
5z	003062	Backrest Cover 19" Rigid / 20" Cat x 13-16"	1	\$244.60
5aa	003063	Backrest Cover 19" Rigid - 20" CAT/LFT x 15-18"	1	\$244.60
5ab	003064	Backrest Cover 19" Rigid - 20" CAT/LFT x 17-20"	1	\$244.60
5ac	003067	Backrest Cover 20" Rigid / 21" Cat x 11-14"	1	\$244.60
5ad	003068	Backrest Cover 20" Rigid / 21" Cat x 13-16"	1	\$244.60
5ae	003069	Backrest Cover 20" Rigid - 21" CAT/LFT x 15-18"	1	\$244.60
5af	003070	Backrest Cover 20" Rigid - 21" CAT/LFT x 17-20"	1	\$244.60
5ag	003072	Backrest Cover 21" Rigid / 22" Cat x 13-16"	1	\$260.36
5ah	003073	Backrest Cover 21" Rigid - 22" CAT/LFT x 15-18"	1	\$260.36
5ai	003074	Backrest Cover 21" Rigid - 22" CAT/LFT x 17-20"	1	\$260.36