


Focus CR Transit



Notice: Chair configuration review required prior to parts being quoted or ordered.

| Item Number | Part Number | Description | Quantity |
|----------------|-------------|--|----------|
| 2a - 10 | 501096 | Transit Retrofit Kit Focus <i>This kit is a complete retrofit, containing a pair of front transit brackets with attaching hardware and transit backplates.</i> | 1 |
| 2b - 9 | 503318 | Transit Retrofit Kit Reclining Back Focus <i>This kit is a complete retrofit, containing a pair of front transit brackets with attaching hardware and transit backplates.</i> | 1 |
| 1, 2a, 5, 7-11 | 501095 | Backrest Hardware Kit Transit Focus (PAIR) | 1 |
| 3 - 7 | 115745 | Rotary Frame Hardware Kit Transit Focus (PAIR) <i>This kit contains a pair of front transit brackets with attaching hardware</i> | 1 |
| 1 | 001071 | Back Plate Std TILT | 2 |
| 2a | 001421 | Back Plate Transit Inner TILT | 2 |
| 2b | 001263 | Back Plate, Recline, Transit Tilt | 2 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| | | <i>For full hardware used with the Reclining Backrest, please see the Focus CR Reclining Backrest page</i> | |
| 4 | 101811 | M6 x 35 BHCS Blk Zn | 4 |
| 5 | 100658 | M6 Nylock Nut Blk Zn | 12 |
| 6 | 100746 | M6 Flat Washer Blk Zn | 4 |
| 9 | 002594 | Saddle, Dual, 1" OD X 1.75 SP | 4 |
| 10 | 101781 | M6 x 50mm Partial Thread BHCS Blk ZC | 8 |
| 11 | 002738 | Spacer .313 OD x .243 ID x 1.68 long | 2 |