


## 20" MAXX Mag Wheel / Tire / Handrim Kits

---



*Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed.*

Item Number	Part Number	Description	Quantity
	118168	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118169	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118170	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118171	20" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118172	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118173	20" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118174	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118175	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118176	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118177	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118178	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118179	20" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	117676	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117677	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117678	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117679	20" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117680	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117681	20" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117714	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117715	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117716	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117717	20" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117718	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117719	20" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
	117748	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117749	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117750	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117751	20" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117752	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117753	20" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount	1