


Catalyst 5VX Step Tube



New rear frames required with serial numbers prior to C5053765, C5HD003546, or CVX0091137 Not available with Heavy Duty Option

Item Number	Part Number	Description	Quantity
1 - 6	505420	Catalyst 5VX Step Tube Assembly Kit	1
1, 2	003535	Step Tube Assembly, CAT5VX	1
3	100860	M6 x 35 HHCS P/THD ZC	2
4	100658	M6 Nylock Nut Blk Zn	2
5	100746	M6 Flat Washer Blk Zn	2
6	003501	Step Tube/Transit Wheel Sleeve, CAT5VX	1