



Quickie Xenon²

January 2022

EIX2

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
1. FRAME									

Base Model

For sizes 14" or less please use HCPCS E1236

☐	107WM96	K0005	Xenon ² - Fixed Front Weight capacity 240 lbs	\$3,146.00
☐	107WM98	K0005	Xenon ² - Swing Away Weight capacity 265 lbs	\$3,146.00
☐	107WM97	K0005	Xenon ² - Hybrid Dual Tube Weight capacity 300 lbs	\$3,146.00

Carbon Package (No substitutions, discount given on separate line)

Proton wheel weight limit of 220 lbs, Auto fold footplate not available.

☐	107QD03		Proton rear wheel	\$215.00
			Platform CF Footplate / leg strap	\$188.00
			Carbon Side Guard w/ Large fender	\$558.00
			Retail Value	\$989.00
			Package Price	\$979.00

Ultra Lite Package (No substitutions, discount given on separate line)

Proton wheel weight limit of 220 lbs

☐	107QD04		Proton rear wheel	\$215.00
			Schwalbe Speed Run Tire	\$151.00
			Compact Wheel Lock	\$151.00
			Retail Value	\$531.00
			Package Price	\$466.00

Frame Angle

Fixed Front & Dual tube only. See **FIGURE 2**

☐	107FA8	80°	NC
☐	107FA16	88°	NC

Frame Inset Per Side

Fixed Front version only. See **FIGURE 1**. N/A with 32cm width.

☐	107FF1	0	NC
☐	107FF12	1cm	NC

Frame - Direct Mount Wheel Lock (Optional)

Only available with Fixed Front Frame and with Compact Wheel Lock, 0° - 2° Camber. No mounting clamp hardware, N/A with Aluminum or carbon fiber sideguards. N/A with mountain tires. N/A in 34-40cm seat depths. N/A with height adjustable armrests. N/A with carbon package. Not available with Desk Sideguard

☐	107WH13	Direct Mount Wheel Lock	NC
--------------------------	---------	-------------------------	----



XENON²
Fixed Front



XENON²
Swing Away



XENON² Hybrid Dual
Tube

2. SIZE / SEATING

Seat Width

See **FIGURE 1** Measured from outside of seat tube to the outside of opposite seat tube. 48-50cm N/A with Fixed Front frame. 1cm inset N/A with 32cm.

☐	107S53	32cm	12.6"	☐	107S58	42cm	16.5"
☐	107S54	34cm	13.4"	☐	107S59	44cm	17.3"
☐	107S55	36cm	14.2"	☐	107S60	46cm	18.1"
☐	107S56	38cm	15.0"	☐	107S61	48cm	18.9"
☐	107S57	40cm	15.7"	☐	107S62	50cm	19.7"

Seat Depth

See **FIGURE 2** Measured from back canes to front edge of upholstery.

☐	107SD81	34cm	13.4"	☐	107SD86	44cm	17.3"
☐	107SD82	36cm	14.2"	☐	107SD87	46cm	18.1"
☐	107SD83	38cm	15.0"	☐	107SD88	48cm	18.9"
☐	107SD84	40cm	15.7"	☐	107SD89	50cm	19.7"
☐	107SD85	42cm	16.5"				

Figure 1





Figure 2



3. CASTERS / SEAT HEIGHT

Front Seat Height

Dependence of front casters, forks and frame version, see Matrix on Page 7 and measure from the top of the seat sling to the floor like in Figure 3.


☐	107SH83	43cm	16.9"	☐	107SH90	50cm	19.7"
☐	107SH84	44cm	17.3"	☐	107SH91	51cm	20.1"
☐	107SH85	45cm	17.7"	☐	107SH92	52cm	20.5"
☐	107SH86	46cm	18.1"	☐	107SH93	53cm	20.9"
☐	107SH87	47cm	18.5"	☐	107SH94	54cm	21.3"
☐	107SH88	48cm	18.9"	☐	107SH95	55cm	21.7"
☐	107SH89	49cm	19.3"				

Caster Fork and Suspension

Multiple & Single Sided forks are std black finish unless anodize package chosen.

☐	107CF65		Multiple Position	NC
☐	107CF69		Multiple Position Single Sided Fork	\$182.00
☐	107CF53	E1015	Frog Legs 0 weight to 100 lbs.	\$467.00
☐	107CF54	E1015	Frog Legs 101 lbs. weight to 190 lbs.	\$467.00
☐	107CF55	E1015	Frog Legs 191 lbs. weight to 265 lbs.	\$467.00

Figure 3



Seating Systems

☐	107OP13		Seat sling, tension adjustable	NC
☐	107OP12		Seat sling, tension adjustable, 1 utility bag	\$135.00
☐	107OP04		Omit Seat Sling	NC

Cushion

Cushion size increments are in inches.

☐	107CU1	K0669	2" Cushion	\$161.00
--------------------------	--------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list on Page 8.

Cushion name or part number must be filled out with desired width and depth

*JAY Cushions are designed to imperial values and listed metric equivalents are not exact measurements.

☐	107JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width:	Depth:	

Caster Wheels

☐	107CW1		3" Micro	NC
☐	107CW3		4" Low Profile Polyurethane	NC
☐	107CW31		4" Micro Lighted	\$75.00
☐	107CW72		4" x 1" Aluminum Poly	\$118.00
☐	107CW44	E2219	4" x 1.25" Semi Pneumatic	\$59.00
☐	107CW5		5" Low Profile Polyurethane	NC
☐	107CW51		5" Micro Lighted	\$92.00
☐	107CW6	E2219	5" Soft Roll (Grey) N/A w/Frog Legs or Single sided fork	\$97.00
☐	107CW65		5" x 1" Aluminum Poly	\$118.00
☐	107CW7	E2219	5" x 1.5" Semi Pneumatic	\$92.00
☐	107CW34		6" Poly	NC
☐	107CW33	E2214 E2215	6" Pneumatic	\$81.00
☐	107CW35	E2219	6" x 1.5" Semi Pneumatic	\$92.00

Caster Color Wheels option

Refer to Anodize package section on Page 6 for Anodize color selection. Cannot choose Match frame color and Anodize package at same time. Casters are silver color unless you choose Match frame color or Anodize package.

☐	107CW64		4" x 1.5" Aluminum Soft Roll*	\$199.00
☐	107CW66		5" x 1.5" Aluminum Soft Roll*	\$209.00
☐	107CW68		6" x 1.5" Aluminum Soft Roll*	\$220.00
	Colors:	☐	Match frame color	\$54.00

Rear Seat Height							
49-50cm RSH only with lightweight adapter. Refer to Seat Squeeze matrix.							
☐	107RH61	37cm	14.6"	☐	107RH68	44cm	17.3"
☐	107RH62	38cm	15.0"	☐	107RH69	45cm	17.7"
☐	107RH63	39cm	15.4"	☐	107RH70	46cm	18.1"
☐	107RH64	40cm	15.7"	☐	107RH71	47cm	18.5"
☐	107RH65	41cm	16.1"	☐	107RH72	48cm	18.9"
☐	107RH66	42cm	16.5"	☐	107RH73	49cm	19.3"
☐	107RH67	43cm	16.9"	☐	107RH74	50cm	19.7"

Rear Wheel			
Some Rear wheel options may limit the chair weight capacity. Proton wheels weight limit is 220 lbs. 24" wheels N/A w/50cm RSH. 25" wheels N/A w/37cm RSH.			
☐	107RW102	24" Mag Wheel 5 Spoke	NC
☐	107RW92	24" Lite Spoke	NC
☐	107RW93	25" Lite Spoke	NC
☐	107RW10	24" Ultra Lightweight Spoke Wheel	\$268.00
☐	107RW16	25" Ultra Lightweight Spoke Wheel	\$268.00
☐	107RW95	24" Spinergy SPOX Black Spokes	\$907.00
☐	107RW96	25" Spinergy SPOX Black Spokes	\$907.00
☐	107RW98	24" Spinergy LX Black Spokes	\$1,138.00
☐	107RW99	25" Spinergy LX Black Spokes	\$1,138.00
☐	107RW50	24" Mountain	\$268.00
☐	107RW148	24" Proton Wheel	\$215.00
☐	107RW149	25" Proton Wheel	\$215.00
☐	107RW150	24" Drum brake wheel - 36 spokes	NC

Rear Wheel Tire			
☐	107RT1	Pneumatic	NC
☐	107RT2	E2211 E2213 Pneumatic Airless Insert	\$92.00
☐	107RT3	Full Poly Tire	NC
☐	107RT38	E2218 SHOX Performance Poly	\$102.00
☐	107RT39	E2211 E2212 High-Pressure Clincher Black	\$102.00
☐	107RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$151.00
☐	107RT28	E2211 E2212 Schwalbe Speed Run / One Tire	\$151.00
☐	107RT40	E2211 E2212 Schwalbe Down Town Two Grip <i>Only available in 24"</i>	\$135.00
☐	107RT41	E2211 E2212 Schwalbe Right Run <i>Only available in 24"</i>	\$92.00
☐	107RT27	E2211 E2212 Knobby V-TRAK by Kenda	\$102.00
☐	107RT12	E2211 E2212 Mountain <i>Only available in 24"</i>	\$102.00

Handrims	Wheels						
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain	Proton	Drum Brake
Aluminum Anodized / Color	24	24, 25	24, 25	24, 25	24	24, 25	24
Plastic Coat	24	24		24, 25	24		24
Natural Fit - Standard Grip	24	24, 25		24, 25			
Natural Fit - No Thumb Grip	24	24, 25		24, 25			
Natural Fit - Super Grip	24	24, 25		24, 25			
Natural Fit LT - Standard Grip	24	24, 25		24, 25			
Natural Fit LT - Super Grip	24	24, 25		24, 25			
Natural Fit LT - No Thumb Grip	24	24, 25		24, 25			
Surge - Oval w/Gription Strip		24, 25		24, 25			
Surge LT - Oval w/Gription Strip		24, 25		24, 25			
Q Grip Neoprene Coated		24, 25		24, 25			
Spinergy Flex Rim				24, 25			
Omit Handrim	24			24, 25	24		

Seat squeeze matrix	
Allowable squeeze (may vary +/- 1.5cm / 0.59")	
Seat depth	Max squeeze
34-36cm	13.4-14.2" 8cm / 3.1"
38-40cm	15.0-15.7" 9cm / 3.5"
42-44cm	16.5-17.3" 10cm / 3.9"
46-50cm	18.1-19.7" 11cm / 4.3"

Max 11cm (4.3") of seat squeeze

Handrim			
See Page 6 for color handrim option			
☐	107HR1	Aluminum Anodized	NC
☐	107HR29	Aluminum Powder coated colored	\$75.00
☐	107HR4	Plastic Coated <i>N/A w/ 25" Lite Spoke Wheels</i>	\$118.00
☐	107HR22	Natural Fit - Standard Grip	\$478.00
☐	107HR23	Natural Fit - Super Grip	\$532.00
☐	107HR25	Natural Fit - No Thumb Grip	\$414.00
☐	107HR30	Natural Fit LT - Standard Grip	\$478.00
☐	107HR31	Natural Fit LT - Super Grip	\$532.00
☐	107HR32	Natural Fit LT - No Thumb Grip	\$414.00
☐	107HR26	The Surge ® - Oval w/Gription Strip	\$429.00
☐	107HR27	The Surge LT ® - Oval w/Gription Strip	\$462.00
☐	107HR28	Q-Grip Neoprene Coated	\$376.00
☐	107HR24	Spinergy Flex Rim <i>Only available on Spinergy LX wheels 24" & 25"</i>	\$751.00
☐	107HR7	Omit Handrim <i>N/A w/ Ultra Lightweight Wheels</i>	NC

Axle			
☐	107AX1	Quick Release Steel	NC
☐	107AX13	Quad Quick Release Steel <i>Only available with Mountain Wheels</i>	\$135.00

Axle Plate			
☐	107AL36	Axle plate Standard	NC

Camber			
Light weight adapter available only 0°			
☐	107CI1	0°	NC
☐	107CI13	2°	NC
☐	107CI3	4°	NC

Center of Gravity Preset			
Lightweight adapter N/A with Single post or Desk side guard armrests.			
☐	107LQ01	Standard Adapter. Choose your setting:	
☐	1cm (0.4")	☐ 2cm (0.8")	☐ 3cm (1.2") ☐ 4cm (1.6")
☐	5cm (2.0")	☐ 6cm (2.4")	☐ 7cm (2.8") ☐ 8cm (3.1")
☐	9cm (3.5")		
OR			
☐	107LQ02	Lightweight Adapter. Choose your setting:	
☐	4cm (1.6")	☐ 5cm (2.0")	☐ 6cm (2.4") ☐ 7cm (4.3")
☐	9cm (3.5")	☐ 10cm (3.9")	

Tire	Wheels					
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Proton	Drum Brake
Pneumatic	24	24	24	24	24	24
Pneumatic w/ Airless Insert	24	24		24		24
Full Poly	24	24		24		24
SHOX Performance Poly	24	24, 25		24, 25		24
High-Pressure Clincher Blk	24	24	24	24	24	24
Schwalbe Marathon	24	24, 25	24, 25	24, 25	24, 25	24
Schwalbe Speed Run/One	24	24, 25	24, 25	24, 25	24, 25	24
Schwalbe Down Town	24	24	24	24	24	24
Schwalbe Right Run	24	24	24	24	24	24
Knobby V Trak by Kenda	24	24	24	24	24	
24" = 540 mm 25" = 559 mm						

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
5. WHEEL LOCKS									
Wheel Lock					Attendant Brake				
Additional No Gos may exist depending on wheel size and position. Recommend to order a wheel lock option in combination with the Drum brake.					Must select Drum brake rear wheel. N/A w/ Half fold down back & Omit push handles.				
☐	107WL1		High Mount Push	NC	☐	107WL79	E2228	Drum Brake	\$537.00
☐	107WL2		High Mount Pull	NC					
☐	107WL27		Ergo Scissor	\$118.00	Wheel Lock Option				
☐	107WL30		Compact	\$151.00	☐	107WL12		Do Not Mount	NC
☐	107WL10		Omit	NC	Extension Handle				
					N/A w/ ergo scissor or compact wheel locks				
					☐	107EH6	E0961	Extension Handle (Pair)	\$81.00

6. HANGER & FOOTREST OPTIONS									
Footrest (Choose ONE Option)					Footplates (Fixed Front & Dual Tube only)				
Hangers (Swing away only)					Platform footplates				
ELR only available 38cm-50cm width.					☐	107F332		Platform Aluminum - Flip Up ☐ Black NC ☐ Std Frame color\$51	NC
☐	107H168		70° Swing In/Out Hanger	NC	☐	107F334	K0040	Platform Composite Footrest - Flip up (sideways), angle adjustable	\$123.00
☐	107H169		80° Swing In/Out Hanger	NC	☐	107F337	K0040	Platform Composite High Mount - Flip up, Lower Leg Length 22-34 cm	\$107.00
☐	107H170	E0990	Elevating Hanger (0° - 90°) (pair)	\$338.00	☐	107F335	K0040	Platform Ultra Lite Footrest - Flip up (sideways) Aluminum, angle adj.	\$123.00
☐	107H32		Omit Hanger & Footplate	-\$15.00	☐	107F336		Platform Carbon Fiber footrest - Flip Up (sideways), N/A with 32cm widths	\$145.00
Footplates (Swing Away only)					Auto-folding Platform				
Z-finity footplates & Platform Alum High Mount Max weight capacity 165 lbs.					☐	107F338		Platform Carbon Fiber footrest - Flip Up (sideways), angle adjustable	\$268.00
Divided footplates					☐	107F339	K0040	Platform Composite Footrest - Flip up (sideways), angle adjustable	\$215.00
☐	107F1		Composite - Flip up	NC	Footrest Options				
☐	107F330	K0040 E0951	Composite angle adj - Flip up (sideways) 38-50cm width	\$194.00	Platform Swingaway Mechanism only available with Platform Aluminum				
☐	107F331	K0040	Alum angle adj - Divided, Flip up ☐ Black NC ☐ Std Frame color\$51	\$204.00	☐	107DI31	E0951	Heel Loop (Swing away only)	\$43.00
☐	107F320	K0040	Multi Angle Adustable (Z-finity) ☐ Black NC ☐ Std Frame color\$51	\$199.00	☐	107DI03	K0038	Leg Strap	\$43.00
Platform footplates					☐	107DI04		Extra Leg Strap	\$43.00
☐	107F332		Platform Aluminum - Flip Up ☐ Black NC ☐ Std Frame color\$51	NC	☐	107DI05		Gel Hanger Cover	\$113.00
☐	107F333	K0040	Platform Alum AA High Mount - Flip up (Z-finity)	\$161.00	☐	107DI07		Neoprene Impact Guards	\$97.00
☐	107F334	K0040	Platform Composite Footrest - Flip up (sideways), angle adjustable	\$123.00	☐	107DI06		Platform Position Plates	\$54.00
☐	107F335	K0040	Platform Ultra Lite Footrest - Flip up (sideways) angle adjustable	\$123.00	☐	107DI32		Platform Swingaway Mechanism	\$92.00
☐	107F336		Platform Carbon Fiber footrest - Flip Up (sideways), N/A with 32cm widths	\$145.00					
☐	107F6		Omit Footplate	-\$15.00					

Lower Leg Length (LLL)							
Lower Leg length standard clearance must be minimum of 5 cm (2") less than Front Seat Height. All measurements are in cm. <u>Measured from seat sling to the top - rear of the footrest.</u> See Appendix on Page 7							
☐ 22cm (8.7")	☐ 26cm (10.2")	☐ 30cm (11.8")	☐ 34cm (13.4")	☐ 38cm (15.0")	☐ 42cm (16.5")	☐ 46cm (18.1")	☐ 50cm (19.7")
☐ 23cm (9.1")	☐ 27cm (10.6")	☐ 31cm (12.2")	☐ 35cm (13.8")	☐ 39cm (15.4")	☐ 43cm (16.9")	☐ 47cm (18.5")	
☐ 24cm (9.4")	☐ 28cm (11.0")	☐ 32cm (12.6")	☐ 36cm (14.2")	☐ 40cm (15.7")	☐ 44cm (17.3")	☐ 48cm (18.9")	
☐ 25cm (9.8")	☐ 29cm (11.4")	☐ 33cm (13.0")	☐ 37cm (14.6")	☐ 41cm (16.1")	☐ 45cm (17.7")	☐ 49cm (19.3")	

7. BACKREST									
Back Type					Back Height Setting				
Half fold down option only available with 32.5cm back height and up, Fold down push handle N/A with height adj push handles or JAY backs. Half fold down back type N/A with Std EXO Uphalf fold down with fold down push handles N/A with back heights 35cm or lower.					All Back height settings are in cm.				
☐	107BT123		Fixed (90°)	NC	☐	107B254	Setting	25 - 30cm Back Height ☐ 25 (9.8") ☐ 27.5 (10.8") ☐ 30 (11.8")	NC
☐	107BT124		Fixed Half Fold down (to the front)	\$199.00	☐	107B250	Setting	30 - 35cm Back Height ☐ 30 (11.8") ☐ 32.5 (12.8") ☐ 35 (13.8")	NC
☐	107BT125		Angle Adjustable Back (75°-103°)	NC	☐	107B251	Setting	35 - 40cm Back Height ☐ 35 (13.8") ☐ 37.5 (14.8") ☐ 40 (15.7")	NC
☐	107BT126		Angle Adjustable Back Half Fold down	\$199.00	☐	107B252	Setting	40 - 45cm Back Height ☐ 40 (15.7") ☐ 42.5 (16.7") ☐ 45 (17.7")	NC
					☐	107B253	Setting	45 - 47.5cm Back Height ☐ 45 (17.7") ☐ 47.5 (18.7")	NC

7. BACKREST (cont)				
Push Handles				
Std EXO Uph N/A with Fold down or Height Adjustable push handles				
☐	107PH20		Push handle Ergo grip long	NC
☐	107PH5		Fold Down Push Handles Available w/ 32.5-47.5cm height (12.8-18.7")	\$107.00
☐	107PH21		Height adjustable push Handles Available w/ 30-47.5cm height (11.8-18.7"). N/A w/Half fold down back	\$161.00
☐	107PH25		Height adjustable push Handles for JAY Backs N/A w/Half fold down back	\$161.00
☐	107PH22		Omit Push Handles - Straight Post	NC
Backrest Upholstery				
Std EXO Uph available with Omit & Ergo grip push handles, plastic side guards, Single post height adj & padded swing away armrests.EXO Standard and Tension Adj upholstery binding is Standard black . Std Uph N/A with Stabilizer bar.				
☐	107BU88		Standard Upholstery EXO	NC
Colors:	☐ Black	☐ Blue	☐ Orange	☐ Red ☐ Silver NC
☐	107BU89	E2611	Backrest Tension Adj EXO *Reference 186BU70 for billing	\$338.00
Colors:	☐ Black	☐ Blue	☐ Orange	☐ Red ☐ Silver NC
☐	107BU90	E2611	Backrest Tension Adj EXO Pro 3DX *Reference 186BU71 for billing	\$355.00
Colors:	☐ Black	☐ Blue	☐ Orange	☐ Red ☐ Silver NC
☐	107BU14		Omit Back Upholstery	-\$15.00

Stabilizer Bar			
Only available 38cm-50cm seat width. Swing away painted match to frame. Available between 35cm-47.5cm back height.			
☐	107SZ8	Stabilizer Bar	\$123.00
☐	107SZ9	Stabilizer swing away bar Must select height adj push handles	\$161.00

Military Service Patches			
Only available on tension adjustable upholstery and when ordered separately through JAY J3 Back OF			
☐	107VU4	US Air Force Patch	\$107.00
☐	107VU2	US Army Patch	\$107.00
☐	107VU12	US Coast Guard Patch	\$107.00
☐	107VU1	US Marines Patch	\$107.00
☐	107VU3	US Navy Patch	\$107.00
☐	107VU5	Purple Heart Patch	\$107.00
☐	107VU11	Purple Heart with Air Force Patch	\$166.00
☐	107VU9	Purple Heart with Army Patch	\$166.00
☐	107VU13	Purple Heart with Coast Guard Patch	\$166.00
☐	107VU8	Purple Heart with Marines Patch	\$166.00
☐	107VU10	Purple Heart with Navy Patch	\$166.00

JAY J3 Back			
☐	107JB01	J3 Back To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s in Page 9.	Price List

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away					
Some limitations when used with the Angle Adjustable Back.					
☐	107AA162		Padded Swing-Away Fixed height (13")	NC	
☐	107AA18	E0973	Padded Swing-Away HT Adj (11" to 14")	\$242.00	
	Setting	☐ 11"	☐ 12"	☐ 13"	☐ 14"
Please select a setting for adjustable version, will default to 11" if no selection is made					

Armrest - Desk Sideguard				
Not available with Fold Down back option or combined with Side Guard				
☐	107LA01	E0973	Desk sideguard fixed flip up, no pad	\$268.00
☐	107LA02	E0973	Desk sideguard fixed flip up	\$268.00
☐	107LA03	E0973	Desk sideguard height adj. flip up, removable	\$161.00
☐	107LG01		Quad Release - option	NC

Armrest - Single post				
Desk length 33cm armpad available only.				
☐	107AA161	E0973	Single post height adj. removable, depth adjustable	\$306.00

Armpad - For Height Adjustable				
* Desk Sideguard only. **Single Post only				
☐	107AM29		Desk Length 25cm/9.8" height *	NC
☐	107AM30		Desk Length 31cm/12.2" height *	NC
☐	107AM31		Desk Length 33cm/13" height **	NC

Side Guards				
Increases overall width of chair. Some interference may exist w/ padded SA arms at lowest settings Not available w/Single post HA or Desk Sideguard armrests.				
☐	107SG33	K0108	Plastic Standard	\$156.00
☐	107SG29	K0108	Plastic Low	\$156.00
☐	107SG35	K0108	Aluminum	\$263.00
☐	107SG36	K0108	Aluminum w/ Fender	\$296.00
☐	107SG12	K0108	Carbon Fiber	\$494.00
☐	107SG41	K0108	Carbon Fiber Long w/LG Fender	\$558.00

Side Guard / Desk Sideguard - Color Option				
Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber or plastic option. Standard colors only				
☐	107GE01		Black	NC
☐	107GE02		Matching Frame- Standard Colors	NC

9. KEY OPTIONS AND ACCESSORIES								
Anti Tips					Additional Pack Accessories			
☐	107AT1	E0971	Anti Tips (Pair)	\$161.00	☐	107A1	Backpack - Adult	\$81.00
☐	107AT81	E0971	Anti Tip swing away		☐	107A2	Backpack - Kids	\$64.00
			☐ Right Mnt ☐ Left Mnt	\$135.00	☐	1071A01	Seat Pouch Black	\$54.00
☐	107AT83		Tip assist					
			☐ Right Mnt ☐ Left Mnt	\$33.00				
Tray Tables					Positioning Belt			
Swing Away Frame only. Only available with Single Post Height Adjustable armrests.					☐	107A833	E0978 Universal 1.5" Auto Buckle Belt	\$75.00
☐	107A44	E0950	Small 33-35cm (13"-14")	\$274.00				
☐	107A45	E0950	Medium 38-43cm (15"-17")	\$274.00				
☐	107A46	E0950	Large 45-50cm (18"-20")	\$274.00				
☐	107A47	E0950	X-Large 53-55cm (21"-22")	\$274.00				
					Spoke Guards			
					Only available with 24" Lite spoke wheels.			
☐	107A401	K0065	Quickie Black Spoke Guard	\$123.00				
☐	107A402	K0065	Quickie Clear Spoke Guard	\$123.00				

9. KEY OPTIONS AND ACCESSORIES (cont)

Miscellaneous

Rubber cover only available with Aluminum Anodized & Painted handrims

☐	1078A48	E2207	Crutch Holder	\$209.00
☐	1078A21		Tool Kit	\$59.00
☐	1078A01		Gel Swing Away Pad	\$161.00
☐	1078A46		Handrim rubber cover (24" only)	\$113.00
☐	1078A03		Caddy	\$54.00
☐	1078A47		Mobile phone pocket	\$27.00
☐	107A999		Quickie Water Bottle & Holder	\$59.00

Built For Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333-9011



Positioning

Pad size is 6"x9". Swing Away frame only.

☐	1078A49	E1020	Lower Extremity S/A support	\$424.00
			☐ Right ☐ Left	

10. COLOR

Standard Frame Color

Handrim color option only available if Aluminum Coated colored option was selected.

Black not available as base w/Zebra

		Handrim
☐	107C140	Aztec Gold
☐	107C118	Bead Blast
☐	107C2	Black
☐	107C141	Black Cherry
☐	107C10	Black Opal
☐	107C38	Blue Opal
☐	107C36	Candy Blue
☐	107C14	Candy Purple
☐	107C35	Candy Red
☐	107C142	Electric Blue
☐	107C143	Evergreen
☐	107C44	Glossy White
☐	107C126	Glow
☐	107C144	Green Apple
☐	107C147	Hot Sparkle Pink
☐	107C148	Mauve Pink
☐	107C103	Rootbeer
☐	107C114	Sparkle Silver
☐	107C145	Sunrise Orange
☐	107C113	Titanium Color Paint
☐	107C8	Yellow

Matte Frame Color

Handrim color option only available if Aluminum Coated colored options was selected.

		Handrim
☐	107C116	Matte Black
☐	107C151	Matte Black Cherry
☐	107C150	Matte Electric Blue
☐	107C149	Matte Evergreen
☐	107C152	Matte Purple

Personalized Frame Color

☐	107C112	Desert Camouflage	\$424.00
☐	107C43	Mossy Oak Camouflage	\$424.00
☐	107C42	Stars & Stripes	\$424.00
☐	107C108	Zebra (Choose above base color) Sid with White unless above base paint is chosen N/A w/ Black, Matte Black, Black Opal, Titanium or Bead Blast)	\$424.00

Anodized Color Package

Forks and Aluminum casters only. N/A with Frog Legs and Alum Caster color option

☐	107EB15	Blue	NC
☐	107EB19	Red	NC
☐	107EB16	Gold	\$92.00
☐	107EB17	Lime Green	\$92.00
☐	107EB18	Purple	\$92.00
☐	107EB20	Ti Grey	\$92.00

Frame Graphics

Graphic designs on front frame tubes. N/A Personalized colors & Matte Black Cherry.

N/A with Fixed Front low height frames or Swing Away frame

☐	1078A24	Bloom	\$327.00
☐	1078A25	Dragon	\$327.00
☐	1078A26	Ghost Flames	\$327.00
☐	1078A27	Retro Squares	\$327.00
☐	1078A28	Tribal	\$327.00

Bloom design



Retro Squares design



Dragon design



Tribal design



Ghost Flames design



11. APPENDIX


Swing Away frame Front Seat Height matrix

FRONT SEAT HEIGHT SELECTION: in dependence of front casters, forks and frame versions

Front seat height in centimeters

		16.9"	17.3"	17.7"	18.1"	18.5"	18.9"	19.3"	19.7"	20.1"	20.5"	20.9"	21.3"	21.7"
Front Caster	Fork	43	44	45	46	47	48	49	50	51	52	53	54	55
	3"	Multi Post.	○		○	○	○			○	○	○		
		Single Side	○		○	○	○							
		Frog Leg												
	4"	Multi Post.	○		○	○	○		○	○	○		○	
		Single Side	○		○	○	○		○	○				
		Frog Leg				○	○		○	○				
	5"	Multi Post.				○	○	○	○		○	○	○	
		Single Side				○	○	○						
		Frog Leg					○		○	○				
	6"	Multi Post.					○	○		○	○	○	○	○
		Single Side					○	○		○	○	○	○	○
		Frog Leg								○	○			

Lower Leg Length



Fixed Front/Dual Tube frame Front Seat Height matrix

FRONT SEAT HEIGHT SELECTION: in dependence of front casters, forks and frame versions

Front seat height in centimeters

		16.9"	17.3"	17.7"	18.1"	18.5"	18.9"	19.3"	19.7"	20.1"	20.5"	20.9"	21.3"	21.7"
Front Caster	Fork	43	44	45	46	47	48	49	50	51	52	53	54	55
	3"	Multi Post.	○	○										
		Single Side	○	○										
		Frog Leg												
		Multi Post.			○	○	○	○						
		Single Side			○	○	○	○						
		Frog Leg						○	○		○	○		
	4"	Multi Post.	○	○	○	○								
		Single Side	○	○	○	○								
		Frog Leg				○	○	○						
		Multi Post.					○	○	○	○	○	○		
		Single Side					○	○	○	○	○	○		
		Frog Leg						○	○	○	○	○	○	
	5"	Multi Post.			○	○	○	○						
		Single Side			○	○	○	○						
		Frog Leg					○	○	○	○				
		Multi Post.						○	○	○	○	○	○	
		Single Side						○	○	○	○	○	○	
		Frog Leg							○	○	○	○	○	○
	6"	Multi Post.							○	○	○	○	○	
		Single Side							○	○	○	○	○	
		Frog Leg								○	○	○	○	○
		Multi Post.									○	○	○	○
		Single Side									○	○	○	○
		Frog Leg										○	○	○

Half Fold Down



Back height

		12.8"	13.8"	14.8"	15.7"	16.7"	17.7"	18.7"
		32.5	35	37.5	40	42.5	45	47.5
Seat width	12.6"	32	✓	✓	✓	✓	✓	✓
	13.4"	34	○	✓	✓	✓	✓	✓
	14.2"	36	○	○	✓	✓	✓	✓
	15.0"	38		○	○	✓	✓	✓
	15.7"	40			○	✓	✓	✓
	16.5"	42			○	○	✓	✓
	17.3"	44			○	○	○	✓
	18.1"	46				○	○	○
	18.9"	48					○	○
	19.7"	50						○



=folds 90° to the front on folded chair



=can't be folded full 90° to the front on folded chair in combination with fold-down pushhandles



=can't be folded full 90° to the front on folded chair



=will fold less than 45° to the front on folded chair, folds 90° on unfolded chair



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA

XENON²

JAY® Cushion Options

Note: The JAY cushion options shown below are the most popular with Xenon² users. For additional variations of the cushions shown below or other cushion options, please visit www.SunriseMedical.com or contact Customer Service.



JAY® X2 Cushion - HPCPS Code E2622				
✓	Part Number	Width x Depth (in)*	Width x Depth (cm)*†	Price
☐	JX21414	13.75 x 14	35 x 35.6	\$461
☐	JX21416	13.75 x 16	35 x 40.6	\$461
☐	JX21515	15.25 x 15	38.7 x 38.1	\$461
☐	JX21516	15.25 x 16	38.7 x 40.6	\$461
☐	JX21616	15.75 x 16	40 x 40.6	\$461
☐	JX21617	15.75 x 17	40 x 43.2	\$461
☐	JX21618	15.75 x 18	40 x 45.7	\$461
☐	JX21620	15.75 x 20	40 x 50.8	\$461
☐	JX21716	17.25 x 16	43.8 x 40.6	\$461
☐	JX21717	17.25 x 17	43.8 x 43.2	\$461
☐	JX21816	17.75 x 16	45.1 x 40.6	\$461
☐	JX21818	17.75 x 18	45.1 x 45.7	\$461
☐	JX21820	17.75 x 20	45.1 x 50.8	\$461
☐	JX22020	19.75 x 20	50.2 x 50.8	\$461

JAY® Fusion Reduced Profile w/ Air Insert Single Valve - HPCPS Code E2622				
✓	Part Number	Width x Depth (in)*	Width x Depth (cm)*†	Price
☐	JFUSION1414RPAS	14 x 14	35.6 x 35.6	\$532
☐	JFUSION1415RPAS	14 x 15	35.6 x 38.1	\$532
☐	JFUSION1416RPAS	14 x 16	35.6 x 40.6	\$532
☐	JFUSION1417RPAS	14 x 17	35.6 x 43.2	\$532
☐	JFUSION1418RPAS	14 x 18	35.6 x 45.7	\$532
☐	JFUSION1514RPAS	15 x 14	38.1 x 35.6	\$532
☐	JFUSION1515RPAS	15 x 15	38.1 x 38.1	\$532
☐	JFUSION1516RPAS	15 x 16	38.1 x 40.6	\$532
☐	JFUSION1517RPAS	15 x 17	38.1 x 43.2	\$532
☐	JFUSION1518RPAS	15 x 18	38.1 x 45.7	\$532
☐	JFUSION1519RPAS	15 x 19	38.1 x 48.3	\$532
☐	JFUSION1520RPAS	15 x 20	38.1 x 50.8	\$532
☐	JFUSION1614RPAS	16 x 14	40.6 x 35.6	\$532
☐	JFUSION1615RPAS	16 x 15	40.6 x 38.1	\$532
☐	JFUSION1616RPAS	16 x 16	40.6 x 40.6	\$532
☐	JFUSION1617RPAS	16 x 17	40.6 x 43.2	\$532
☐	JFUSION1618RPAS	16 x 18	40.6 x 45.7	\$532
☐	JFUSION1619RPAS	16 x 19	40.6 x 48.3	\$532
☐	JFUSION1620RPAS	16 x 20	40.6 x 50.8	\$532
☐	JFUSION1714RPAS	17 x 14	43.2 x 35.6	\$532
☐	JFUSION1715RPAS	17 x 15	43.2 x 38.1	\$532
☐	JFUSION1716RPAS	17 x 16	43.2 x 40.6	\$532
☐	JFUSION1717RPAS	17 x 17	43.2 x 43.2	\$532
☐	JFUSION1718RPAS	17 x 18	43.2 x 45.7	\$532
☐	JFUSION1719RPAS	17 x 19	43.2 x 48.3	\$532
☐	JFUSION1720RPAS	17 x 20	43.2 x 50.8	\$532
☐	JFUSION1814RPAS	18 x 14	45.7 x 35.6	\$532
☐	JFUSION1815RPAS	18 x 15	45.7 x 38.1	\$532
☐	JFUSION1816RPAS	18 x 16	45.7 x 40.6	\$532
☐	JFUSION1817RPAS	18 x 17	45.7 x 43.2	\$532
☐	JFUSION1818RPAS	18 x 18	45.7 x 45.7	\$532
☐	JFUSION1819RPAS	18 x 19	45.7 x 48.3	\$532
☐	JFUSION1820RPAS	18 x 20	45.7 x 50.8	\$532
☐	JFUSION1915RPAS	19 x 15	48.3 x 38.1	\$532
☐	JFUSION1916RPAS	19 x 16	48.3 x 40.6	\$532
☐	JFUSION1917RPAS	19 x 17	48.3 x 43.2	\$532
☐	JFUSION1918RPAS	19 x 18	48.3 x 45.7	\$532
☐	JFUSION1919RPAS	19 x 19	48.3 x 48.3	\$532
☐	JFUSION1920RPAS	19 x 20	48.3 x 50.8	\$532
☐	JFUSION2015RPAS	20 x 15	50.8 x 38.1	\$532
☐	JFUSION2016RPAS	20 x 16	50.8 x 40.6	\$532
☐	JFUSION2017RPAS	20 x 17	50.8 x 43.2	\$532
☐	JFUSION2018RPAS	20 x 18	50.8 x 45.7	\$532
☐	JFUSION2019RPAS	20 x 19	50.8 x 48.3	\$532
☐	JFUSION2020RPAS	20 x 20	50.8 x 50.8	\$532

JAY® Lite Cushion - HPCPS Code E2603				
✓	Part Number	Width x Depth (in)*	Width x Depth (cm)*†	Price
☐	JL1414	14 x 14	35.6 x 35.6	\$405
☐	JL1416	14 x 16	35.6 x 40.6	\$405
☐	JL1418	14 x 18	35.6 x 45.7	\$405
☐	JL1515	15 x 15	38.1 x 38.1	\$405
☐	JL1516	15 x 16	38.1 x 40.6	\$405
☐	JL1517	15 x 17	38.1 x 43.2	\$405
☐	JL1518	15 x 18	38.1 x 45.7	\$405
☐	JL1615	16 x 15	40.6 x 38.1	\$405
☐	JL1616	16 x 16	40.6 x 40.6	\$405
☐	JL1617	16 x 17	40.6 x 43.2	\$405
☐	JL1618	16 x 18	40.6 x 45.7	\$405
☐	JL1620	16 x 20	40.6 x 50.8	\$405
☐	JL1715	17 x 15	43.2 x 38.1	\$405
☐	JL1716	17 x 16	43.1 x 40.6	\$405
☐	JL1717	17 x 17	43.1 x 43.1	\$405
☐	JL1718	17 x 18	43.1 x 45.7	\$405
☐	JL1816	18 x 16	45.7 x 40.6	\$405
☐	JL1817	18 x 17	45.7 x 43.2	\$405
☐	JL1818	18 x 18	45.7 x 45.7	\$405
☐	JL1820	18 x 20	45.7 x 50.8	\$405
☐	JL1916	19 x 16	48.3 x 40.6	\$405
☐	JL1917	19 x 17	48.3 x 43.2	\$405
☐	JL1918	19 x 18	48.3 x 45.7	\$405
☐	JL1920	19 x 20	48.3 x 50.8	\$405
☐	JL2016	20 x 16	50.8 x 40.6	\$405
☐	JL2018	20 x 18	50.8 x 45.7	\$405
☐	JL2020	20 x 20	50.8 x 50.8	\$405

* JAY® cushions are designed to imperial values and listed metric equivalents are not exact measurements. Nominal dimensions may vary by ±0.5 in. or ±1.5 cm.

† Not all Xenon² chair sizes will have a JAY® cushion equivalent. Please contact Customer Service to modify standard cushions via JAY Your Way to chair specifications.



JAY® J3 BACKREST OPTIONS

To order a J3 Backrest, first choose a contour from one of the four options (A, B, C or D) below, then select the appropriate backrest height and width.




A) JAY® J3 Posterior (PA) Back - 2.2" Contour												
JAY J3 Back Size (in)			12	14	16	18	20					
Xenon ² Chair Width (cm)			31 - 36	36 - 41	41 - 46	46-51	51-56					
Max Trunk Width (in)			11	13	15	17	19					
Back Height	Lower Thoracic (LT)		\$546									
	6.5" (S)	J3PALTSxxSH		☐	14	☐	16	☐	18			
		J3PALTSxxRH		☐	14	☐	16	☐	18			
		J3PALTSxxRHF		☐	14	☐	16	☐	18			
	8.0" (M)	J3PALTMxxSH	☐	12	☐	14	☐	16	☐	18		
		J3PALTMxxRH	☐	12	☐	14	☐	16	☐	18		
		J3PALTMxxRHF	☐	12	☐	14	☐	16	☐	18		
	9.5" (T)	J3PALTTxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
		J3PALTTxxRH	☐	12	☐	14	☐	16	☐	18	☐	20
		J3PALTTxxRHF	☐	12	☐	14	☐	16	☐	18	☐	20
	Mid Thoracic (MT)		\$546									
	12.0" (S)	J3PAMTSxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
	13.5" (M)	J3PAMTMxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
15.0" (T)	J3PAMTTxxSH	☐	12	☐	14	☐	16	☐	18	☐	20	

B) JAY® J3 Carbon Fiber Posterior (PA) Back - 2.2" Contour												
JAY J3 Back Size (in) Xenon² Chair Width (cm) Max Trunk Width (in)					14		16		18		20	
					36 - 41		41 - 46		46-51		51-56	
					11.5		13.5		15.5		17.5	
Back Height	Lower Thoracic (LT)				\$915							
	6.5" (S)	J3CFPALTSxxRH				☐	14	☐	16	☐	18	
		J3CFPALTSxxRHF				☐	14	☐	16	☐	18	
	9.5" (T)	J3CFPALTTxxRH				☐	14	☐	16	☐	18	☐ 20
		J3CFPALTSTxxRHF				☐	14	☐	16	☐	18	☐ 20
	Mid Thoracic (MT)				\$915							
	12.0" (S)	J3CFPAMTSxxRH				☐	14	☐	16	☐	18	☐ 20
		J3CFPAMTSxxRHF				☐	14	☐	16	☐	18	☐ 20
	15.0" (T)	J3CFPAMTTxxRH				☐	14	☐	16	☐	18	☐ 20
		J3CFPAMTSTxxRHF				☐	14	☐	16	☐	18	☐ 20

C) JAY® J3 Posterior Deep Contour (PDC) Back - 6" Contour												
JAY J3 Back Size (in)				14		16		18		20		
Xenon ² Chair Width (cm)				36 - 41		41 - 46		46-51		51-56		
Max Trunk Width (in)				11.5		13.5		15.5		17.5		
Back Height	Lower Thoracic (LT)			\$765								
	8.0" (M)	J3PDCLTMxxSH		☐	14	☐	16	☐	18			
		J3PDCLTMxxRH		☐	14	☐	16	☐	18			
		J3PDCLTMxxRHF		☐	14	☐	16	☐	18			
	9.5" (T)	J3PDCLTTxxSH		☐	14	☐	16	☐	18	☐	20	
		J3PDCLTTxxRH		☐	14	☐	16	☐	18	☐	20	
		J3PDCLTTxxRHF		☐	14	☐	16	☐	18	☐	20	
	Mid Thoracic (MT)			\$765								
	12.0" (S)	J3PDCMTSxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
15.0" (T)	J3PDCMTTxxSH	☐	12	☐	14	☐	16	☐	18	☐	20	

D) JAY® J3 Posterior Lateral (PL) Back - 3" Contour												
JAY J3 Back Size (in)			12		14		16		18		20	
Xenon ² Chair Width (cm)			31 - 36		36 - 41		41 - 46		46-51		51-56	
Max Trunk Width (in)			9.5		11.5		13.5		15.5		17.5	
Back Height	Mid Thoracic (MT)		\$638									
	12.0" (S)	J3PLMTSxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
	13.5" (M)	J3PLMTMxxSH	☐	12	☐	14	☐	16	☐	18	☐	20
	15.0" (T)	J3PLMTTxxSH	☐	12	☐	14	☐	16	☐	18	☐	20

Note: The J3 Back options shown are the most popular with Xenon² users. For more contours, back heights, hardware options, and accessories, please see the J3 Back order form at www.SunriseMedical.com.

JAY® J3 Back Mounting Hardware Options			
SH		RH	
		RHF	



Sunrise Medical (US) LLC
2842 Business Park Ave. • Fresno, CA 93727 • USA