


Liberty Heel Loops



Part #	Used on Chair Widths			Strap Length
	Angle Adjustable	Composite	Composite Angle Adjustable	
100591	N/A	14"-15"	N/A	9.4"
100592	14"-15"	16"-17"	15"-16"	10.4"
100593	16"-22"	18"-22"	17"-22"	11.4"

Item Number	Part Number	Description	Quantity
1a	100591	Heel Loop Strap Short	1
1b	100592	Heel Loop Strap Medium	1
1c	100593	Heel Loop Strap Long	1
2	100551	Heel Loop Peg	1

Item Number	Part Number	Description	Quantity
3	101471	10-32 x 1-1/4 PH Pan HD T-B T/S ZC	1