

Account Information

Company Name: _____
Address: _____
City, State Zip: _____
Phone: _____
Fax: _____
Email: _____
Purchase Order: _____

Mark For: _____
Ship To: _____
Address: _____
City, State Zip: _____
Phone: _____

Order Information

To insure proper Meywalk sizing, please include required measurements below

Miniwalk I

Seat Height: 10 1/2" - 19"
Trunk Support Height: 20" - 29"
Ring Size: 8 1/4" w x 8" d
Weight Capacity: 90 lbs

106061 ☐ Miniwalk I 4,235.00 Ea
Includes: Castor Swivel Locks. Color: Turquoise

Miniwalk II

Seat Height: 14 - 25"
Trunk Support Height: 27-35"
Ring Size: 8 1/4" w x 8" d
Weight Capacity: 90 lbs

106062 ☐ Miniwalk II 4,235.00 Ea
Includes: Castor Swivel Locks. Color: Turquoise

Miniwalk Accessories

106000	☐ Nylon Hygiene Ring Cover, 70cm	95.00	Ea
106070	☐ Antitip support	358.00	Pr
106075	☐ Leg Separation Panel	726.00	Ea
104302	☐ Rear Inside Wheel Guard	660.00	Pr
106085	☐ Hip Pads	714.00	Pr
106089	☐ No-Roll Back Attachment	455.00	Pr
106096	☐ Drag brakes	358.00	Pr
11761	☐ Miniwalk Push / Pull Bar	455.00	Ea

Meywalk 2000 Small **New**

Seat Height: 19"-26"
Trunk Support Height: 30 3/4"-38 1/2"
Ring Size: 10"w x 10"d
Weight Capacity: 150lbs

104500 ☐ Meywalk 2000 Small 5,228.00 Ea
Color: Green

Meywalk 2000 Medium

Seat Height: 21"-29"
Trunk Support Height: 33"-40"
Ring Size: 12 3/4"w x 11 1/2"d
Weight Capacity: 220lbs

104175 ☐ Meywalk 2000 Medium 5,258.00 Ea
Color: Blue

Meywalk 2000 Large

Seat Height: 27" -37"
Trunk Support Height: 41" - 51"
Ring Size: 14 1/2" w x 13 3/4" d
Weight Capacity: 220lbs

104173 ☐ Meywalk 2000 Large 5,258.00 Ea
Color: Blue

Meywalk 2000 Accessories

104000	☐ Nylon Hygiene Ring Cover, 105cm	95.00	Ea
104001	☐ Nylon Hygiene Ring Cover, 120cm	95.00	Ea
104002	☐ Nylon Hygiene Ring Cover, 90cm	95.00	Ea
104023	☐ No-Roll Back Attachment	455.00	Pr
104044	☐ Leg separation plate	726.00	Ea
104301	☐ Rear Inside Wheel Guard	660.00	Pr
104110	☐ Height reducing fittings, -3"	214.00	Pr
104166	☐ Drag Brakes	358.00	Pr
104167	☐ Antitip support	358.00	Pr
104215	☐ Hip Pads	714.00	Pr
104997	☐ Front Castor Swivel Locks	276.00	Pr
104262	☐ Gel Unicycle Seat Cover	825.00	Pr
11761	☐ Push / Pull Bar (Fits Small Only)	455.00	Ea

Required Measurements:

A: Ring Height

Measure user from floor to axilla to indicate ring height.

B: Inseam Height:

Measure user from floor to inseam, with user standing. Indicates seat height.

C: Trunk Width:

Measure user across trunk, at axilla height, side to side in a straight line, to indicate proper ring size.

