

Catalyst 5 Casters - 1.5" Soft Roll Aluminum



If changing size of casters, new forks may be needed.

Item Number	Part Number	Description	Quantity
1a	116737	4 x 1.5 Aluminum Soft Roll Caster <i>Replaced part number 100810. Effective 4/2023.</i>	1
1b	116738	5 x 1.5 Aluminum Soft Roll Caster <i>Replaced part number 100812</i>	1
1c	116739	6 x 1.5 Aluminum Soft Roll Caster <i>Replaced part number 100813</i>	1