


Ethos Casters - 1.25" Soft Roll (Poly)



If changing size of casters, new forks may be needed. If changing caster size see Ethos Front Base Frame to determine if new Caster Sleeve is needed.

Item Number	Part Number	Description	Quantity
1a	100809	Caster 4 x 1.25" Soft Roll <i>Replaced part number 100517</i>	1
1b	100807	Caster 5 x 1.25" Soft Roll <i>Replaced part number 100518</i>	1
1c	100786	Caster 6 x 1.25" Soft Roll <i>Replaced part number 100521</i>	1