


Rogue XP Plastic Coated Handrim



For more information see Handrim Hardware Chart. Quantity in Handrim Hardware Kits are for one wheel.<p>
 Tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
8-10	505263	Handrim Hardware Kit, 3 Rivet Nut, EACH <i>Used with 18" Ki Spoke and 18" MAXX Spoke wheels</i>	1	\$19.31

Item Number	Part Number	Description	Quantity	Retail
11	101106	Handrim Plastic Coated 3 Rivet Nut 18" <i>Used on 18" Ki Spoke and 18" MAXX Spoke wheels</i>	1	\$74.28
8-10	502168	Handrim Hardware Kit, 6 Rivet Nut, EACH <i>Used with Ki Spoke and MAXX Spoke wheels</i>	1	\$38.61
7a	200542	Handrim Plastic Coated 6 Rivet Nut 20" <i>Used on 20" Ki Spoke and 20" MAXX Spoke wheels</i>	1	\$160.34
7b	100576	Handrim Plastic Coated 6 Rivet Nut 22" <i>Used on 22" Ki Spoke and 22" MAXX Spoke wheels</i>	1	\$143.41
7c	112821	Handrim Plastic Coated 6 Rivet Nut 24" <i>Used on 24" Ki Spoke and 24" MAXX Spoke wheels. Replaced 100577.</i>	1	\$143.41
7d	101870	Handrim Plastic Coated 6 Rivet Nut 25" <i>Used on 25" Ki Spoke and 25" MAXX Spoke wheels</i>	1	\$160.34
5, 6	502173	Handrim Hardware Kit, 5 Rivet Nut, Plastic Coated, EACH <i>Used with X Core Mag wheels</i>	1	\$33.10
4a	101963	Handrim Plastic Coated 5 Rivet Nut 20" <i>Used on 20" 5-Spoke X Core Mag wheels</i>	1	\$140.71
4b	101964	Handrim Plastic Coated 5 Rivet Nut 22" <i>Used on 22" 5-Spoke X Core Mag wheels</i>	1	\$140.71
4c	100976	Handrim Plastic Coated 5 Rivet Nut 24" <i>Used on 24" 5-Spoke X Core Mag wheels</i>	1	\$132.37
2a, 3	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinerger and Topolino wheels</i>	1	\$16.55
2b	502170	Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Spoke wheels</i>	1	\$16.55
1a	100836	Handrim Plastic Coated 6 Tab 24" <i>Used on 24" Superlight/Topolino wheels</i>	1	\$163.98
1b	101091	Handrim Plastic Coated 6 Tab 25" <i>Used on 25" Spoke/Superlight/Topolino wheels</i>	1	\$163.87
1c	103179	Handrim Plastic Coated Spinerger 6 Tab 20"	1	\$204.08
1d	100808	Handrim Plastic Coated Spinerger 6 Tab 22"	1	\$204.08
1e	112815	Handrim Plastic Coated Spinerger 6 Tab 24" <i>If chair is invoiced prior to 8/14/19, quantity 2 must be ordered due to a change in finish. Replaced 100615.</i>	1	\$194.36
1f	100765	Handrim Plastic Coated Spinerger 6 Tab 25"	1	\$204.08
2a	100669	M5 x 10 BTN SHCS BLKZC w/Patch	6	\$1.66

Item Number	Part Number	Description	Quantity	Retail
		<i>Used with Spinergy and Topolino wheels</i>		
2b	100536	M5 x 12 BTN SHCS ZC <i>Used with Superlight Spoke wheels</i>	6	\$2.21
3	100657	M5 Nylock Nut Blk Zn <i>Used with Spinergy and Topolino wheels</i>	6	\$0.56
5	101756	Handrim Spacer Plastic Coated Riv <i>Used with X Core Mag wheels</i>	5	\$0.56
6	100724	M5 x 55 PH TRS HD M/S ZC <i>Used with X Core Mag wheels</i>	5	\$1.10
8	100653	Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke wheels</i>	6	\$1.10
9	100698	Handrim Tab <i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i>	6	\$4.17
10	100654	M5 x 40 DIN7985 PH Pan HD M/S ZC <i>Used with Ki Spoke and MAXX Spoke wheels</i>	6	\$1.66

Handrim Hardware Chart

Wheel		Handrim Connection				Handrim									
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim			
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654	100206		101106							
20" Ki Spoke	200530					200536		200542							
22" Ki Spoke	200531					100560		100576	100569	200538	200201				
24" Ki Spoke	200532				6	100698	100653	Natural Fit: 100835	200349		100577	200547	100793	200202	
25" Ki Spoke	200533								200350		101870	200548	200539	200540	
26" Ki Spoke	200534								200351		100578	200549	100907	101454	
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828				
25" Superlight	101160							101160	101091		101464	101460			
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629	100724	101897		101963							
22" 5-Spoke X Core	101962						101898		101964						
24" 5-Spoke X Core	100960				Plastic Ctd: 101756	101893	100975		100976	200546					
24" 5-Spoke X Core	100960										100768	100769			
20" Spinergy Spox	See Spinergy Spox Page	6			Screw: 100669 Nut: 100657	103125		103179							
22" Spinergy Spox						100827		100808		100889	100888				
24" Spinergy Spox						100766	101161	100615		100830	100828				
25" Spinergy Spox						100767	101160	100765		101464	101460				
26" Spinergy Spox						101477		101148		200200	100950				
22" Spinergy LX	See Spinergy LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888				
24" Spinergy LX						100766	101161	100615		100830	100828	200213			
25" Spinergy LX						100767	101160	100765		101464	101460	200208			
26" Spinergy LX						101477		101148		200200	100950				
18" MAXX Spoke	109285	6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654	100206		101106							
20" MAXX Spoke	108244					200536		200542							
22" MAXX Spoke	105135					100560		100576	100569	200538	200201				
24" MAXX Spoke	105136				6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Natural Fit: 100835	200349	101161	100577	200547	100793	200202	
25" MAXX Spoke	107436								200350		101870	200548	200539	200540	
26" MAXX Spoke	107437								200351		100578	200549	100907	101454	

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.