


(Discontinued) Spinerge Backrest Release Cables



For complete cable assembly the backrest cable release springs and backrest lockout pins are required Spinerge Back Release Cable based on seat width The Spinerge Back Release Cable length will measure 1" less than the specified seat width

Item Number	Part Number	Description	Quantity
1b	102557	Spinerge Back Release Cable 8" Wide Blue	1
1c	102558	Spinerge Back Release Cable 8" Wide Red	1
1d	102559	Spinerge Back Release Cable 8" Wide Yellow	1
1e	102560	Spinerge Back Release Cable 8" Wide Green	1

Item Number	Part Number	Description	Quantity
1f	102561	Spinergy Back Release Cable 8" Wide Pink	1
1g	102562	Spinergy Back Release Cable 8" Wide Orange	1
1h	102563	Spinergy Back Release Cable 8" Wide White	1
1j	102565	Spinergy Back Release Cable 9" Wide Blue	1
1k	102566	Spinergy Back Release Cable 9" Wide Red	1
1l	102567	Spinergy Back Release Cable 9" Wide Yellow	1
1m	102568	Spinergy Back Release Cable 9" Wide Green	1
1n	102569	Spinergy Back Release Cable 9" Wide Pink	1
1o	102570	Spinergy Back Release Cable 9" Wide Orange	1
1p	102571	Spinergy Back Release Cable 9" Wide White	1
1r	102573	Spinergy Back Release Cable 10" Wide Blue	1
1s	102574	Spinergy Back Release Cable 10" Wide Red	1
1t	102575	Spinergy Back Release Cable 10" Wide Yellow	1
1u	102576	Spinergy Back Release Cable 10" Wide Green	1
1v	102577	Spinergy Back Release Cable 10" Wide Pink	1
1w	102578	Spinergy Back Release Cable 10" Wide Orange	1
1ab	102583	Spinergy Back Release Cable 11" Wide Yellow	1
1ac	102584	Spinergy Back Release Cable 11" Wide Green	1
1ad	102585	Spinergy Back Release Cable 11" Wide Pink	1
1af	102587	Spinergy Back Release Cable 11" Wide White	1
1ai	200262	Spinergy Back Release Cable 12" Wide Red	1
1aj	200263	Spinergy Back Release Cable 12" Wide Yellow	1
1an	200267	Spinergy Back Release Cable 12" Wide White	1
1ar	200271	Spinergy Back Release Cable 13" Wide Yellow	1
1au	200274	Spinergy Back Release Cable 13" Wide Orange	1
1aw	200276	Spinergy Back Release Cable 14" Wide Black	1
1bc	200282	Spinergy Back Release Cable 14" Wide Orange	1
1bh	200287	Spinergy Back Release Cable 15" Wide Yellow	1

Item Number	Part Number	Description	Quantity
1bk	200290	Spinergy Back Release Cable 15" Wide Orange	1
1bp	200295	Spinergy Back Release Cable 16" Wide Yellow	1
1bq	200296	Spinergy Back Release Cable 16" Wide Green	1
1bt	200299	Spinergy Back Release Cable 16" Wide White	1
1bw	200302	Spinergy Back Release Cable 17" Wide Red	1
1bx	200303	Spinergy Back Release Cable 17" Wide Yellow	1
1by	200304	Spinergy Back Release Cable 17" Wide Green	1
1ca	200306	Spinergy Back Release Cable 17" Wide Orange	1
1ce	200310	Spinergy Back Release Cable 18" Wide Red	1
1cf	200311	Spinergy Back Release Cable 18" Wide Yellow	1
1cg	200312	Spinergy Back Release Cable 18" Wide Green	1
1ch	200313	Spinergy Back Release Cable 18" Wide Pink	1
1ci	200314	Spinergy Back Release Cable 18" Wide Orange	1
1cj	200315	Spinergy Back Release Cable 18" Wide White	1
2	000965	Backrest Cable Release Spring RG	2
3	000909	Backrest Lockout Pin	2
4	101991	Spinergy Release Cable/SPOX Adjustment Tool <i>No charge when ordered with Spinergy Release Cable</i>	1