


# Rogue ALX Fixed Height Backrest with Non-Adjustable Rigidizer Bar - Growth (Formerly Tsunami)



|              |        |
|--------------|--------|
| Backrest 14" | 14.13" |
| Backrest 15" | 15.13" |
| Backrest 16" | 16.13" |
| Backrest 17" | 17.13" |
| Backrest 18" | 18.13" |
| Backrest 19" | 19.13" |
| Backrest 20" | 20.13" |
| Backrest 21" | 21.13" |
| Backrest 22" | 22.13" |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1a          | 108997      | Fixed Height Rigid Backrest 7.5", 14" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1b          | 108996      | Fixed Height Rigid Backrest 8.5", 14" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1c          | 108995      | Fixed Height Rigid Backrest 9.5", 14" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1d          | 108994      | Fixed Height Rigid Backrest 10.5", 14" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1e          | 108993      | Fixed Height Rigid Backrest 11.5", 14" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1f          | 108992      | Fixed Height Rigid Backrest 12.5", 14" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1g          | 108991      | Fixed Height Rigid Backrest 13.5", 14" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1h          | 108990      | Fixed Height Rigid Backrest 14.5", 14" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1i          | 108989      | Fixed Height Rigid Backrest 15.5", 14" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1j          | 108988      | Fixed Height Rigid Backrest 16.5", 14" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1k          | 109008      | Fixed Height Rigid Backrest 7.5", 15" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1l          | 109007      | Fixed Height Rigid Backrest 8.5", 15" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1m          | 109006      | Fixed Height Rigid Backrest 9.5", 15" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1n          | 109005      | Fixed Height Rigid Backrest 10.5", 15" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1o          | 109004      | Fixed Height Rigid Backrest 11.5", 15" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1p          | 109003      | Fixed Height Rigid Backrest 12.5", 15" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1q          | 109002      | Fixed Height Rigid Backrest 13.5", 15" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1r          | 109001      | Fixed Height Rigid Backrest 14.5", 15" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1s          | 109000      | Fixed Height Rigid Backrest 15.5", 15" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1t          | 108999      | Fixed Height Rigid Backrest 16.5", 15" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1u          | 109019      | Fixed Height Rigid Backrest 7.5", 16" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1v          | 109018      | Fixed Height Rigid Backrest 8.5", 16" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1w          | 109017      | Fixed Height Rigid Backrest 9.5", 16" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1x          | 109016      | Fixed Height Rigid Backrest 10.5", 16" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1y          | 109015      | Fixed Height Rigid Backrest 11.5", 16" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1z          | 109014      | Fixed Height Rigid Backrest 12.5", 16" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1aa         | 109013      | Fixed Height Rigid Backrest 13.5", 16" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1ab         | 109012      | Fixed Height Rigid Backrest 14.5", 16" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1ac         | 109011      | Fixed Height Rigid Backrest 15.5", 16" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1ad         | 109010      | Fixed Height Rigid Backrest 16.5", 16" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1ae         | 109030      | Fixed Height Rigid Backrest 7.5", 17" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1af         | 109029      | Fixed Height Rigid Backrest 8.5", 17" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1ag         | 109028      | Fixed Height Rigid Backrest 9.5", 17" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1ah         | 109027      | Fixed Height Rigid Backrest 10.5", 17" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1ai         | 109026      | Fixed Height Rigid Backrest 11.5", 17" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1aj         | 109025      | Fixed Height Rigid Backrest 12.5", 17" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1ak         | 109024      | Fixed Height Rigid Backrest 13.5", 17" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1al         | 109023      | Fixed Height Rigid Backrest 14.5", 17" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1am         | 109022      | Fixed Height Rigid Backrest 15.5", 17" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1an         | 109021      | Fixed Height Rigid Backrest 16.5", 17" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1ao         | 109041      | Fixed Height Rigid Backrest 7.5", 18" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1ap         | 109040      | Fixed Height Rigid Backrest 8.5", 18" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1aq         | 109039      | Fixed Height Rigid Backrest 9.5", 18" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1ar         | 109038      | Fixed Height Rigid Backrest 10.5", 18" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1as         | 109037      | Fixed Height Rigid Backrest 11.5", 18" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1at         | 109036      | Fixed Height Rigid Backrest 12.5", 18" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1au         | 109035      | Fixed Height Rigid Backrest 13.5", 18" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1av         | 109034      | Fixed Height Rigid Backrest 14.5", 18" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1aw         | 109033      | Fixed Height Rigid Backrest 15.5", 18" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1ax         | 109032      | Fixed Height Rigid Backrest 16.5", 18" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1ay         | 109052      | Fixed Height Rigid Backrest 7.5", 19" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1az         | 109051      | Fixed Height Rigid Backrest 8.5", 19" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1ba         | 109050      | Fixed Height Rigid Backrest 9.5", 19" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1bb         | 109049      | Fixed Height Rigid Backrest 10.5", 19" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1bc         | 109048      | Fixed Height Rigid Backrest 11.5", 19" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1bd         | 109047      | Fixed Height Rigid Backrest 12.5", 19" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1be         | 109046      | Fixed Height Rigid Backrest 13.5", 19" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1bf         | 109045      | Fixed Height Rigid Backrest 14.5", 19" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1bg         | 109044      | Fixed Height Rigid Backrest 15.5", 19" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1bh         | 109043      | Fixed Height Rigid Backrest 16.5", 19" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1bi         | 109063      | Fixed Height Rigid Backrest 7.5", 20" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1bj         | 109062      | Fixed Height Rigid Backrest 8.5", 20" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1bk         | 109061      | Fixed Height Rigid Backrest 9.5", 20" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1bl         | 109060      | Fixed Height Rigid Backrest 10.5", 20" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1bm         | 109059      | Fixed Height Rigid Backrest 11.5", 20" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1bn         | 109058      | Fixed Height Rigid Backrest 12.5", 20" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1bo         | 109057      | Fixed Height Rigid Backrest 13.5", 20" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1bp         | 109056      | Fixed Height Rigid Backrest 14.5", 20" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1bq         | 109055      | Fixed Height Rigid Backrest 15.5", 20" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1br         | 109054      | Fixed Height Rigid Backrest 16.5", 20" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 1bs         | 109074      | Fixed Height Rigid Backrest 7.5", 21" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1bt         | 109073      | Fixed Height Rigid Backrest 8.5", 21" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1bu         | 109072      | Fixed Height Rigid Backrest 9.5", 21" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1bv         | 109071      | Fixed Height Rigid Backrest 10.5", 21" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1bw         | 109070      | Fixed Height Rigid Backrest 11.5", 21" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1bx         | 109069      | Fixed Height Rigid Backrest 12.5", 21" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1by         | 109068      | Fixed Height Rigid Backrest 13.5", 21" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1bz         | 109067      | Fixed Height Rigid Backrest 14.5", 21" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1ca         | 109066      | Fixed Height Rigid Backrest 15.5", 21" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1cb         | 109065      | Fixed Height Rigid Backrest 16.5", 21" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1cc         | 109085      | Fixed Height Rigid Backrest 7.5", 22" Wide<br><i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>   | 1        |
| 1cd         | 109084      | Fixed Height Rigid Backrest 8.5", 22" Wide<br><i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>  | 1        |
| 1ce         | 109083      | Fixed Height Rigid Backrest 9.5", 22" Wide<br><i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>  | 1        |
| 1cf         | 109082      | Fixed Height Rigid Backrest 10.5", 22" Wide<br><i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i> | 1        |
| 1cg         | 109081      | Fixed Height Rigid Backrest 11.5", 22" Wide<br><i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i> | 1        |
| 1ch         | 109080      | Fixed Height Rigid Backrest 12.5", 22" Wide<br><i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i> | 1        |
| 1ci         | 109079      | Fixed Height Rigid Backrest 13.5", 22" Wide<br><i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>  | 1        |
| 1cj         | 109078      | Fixed Height Rigid Backrest 14.5", 22" Wide<br><i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>  | 1        |
| 1ck         | 109077      | Fixed Height Rigid Backrest 15.5", 22" Wide<br><i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>  | 1        |
| 1cl         | 109076      | Fixed Height Rigid Backrest 16.5", 22" Wide<br><i>Used with 18" back height with omit push handles</i>   | 1        |
| 2           | 100804      | Tube End Plug 1"<br><i>Used with omit push handles</i>   | 2        |
| 3           | 000034      | Sleeve 3/4" x 1 1/2"<br><i>Used with standard push handles, low profile push handles or fold down push handles</i>   | 2        |