

Klip



**CHOICES
PROVIDE YOU
WITH THE
BEST OPTIONS**



The **Klip** is a **lightweight posterior gait trainer for kids, teens and young adults** with special needs who require support and assistance with their day-to-day walking.

Placed behind the user, the **Klip** is open in the walking position and allows freedom of movement while it promotes walking in an upright position. It is recommended for those with Cerebral Palsy and other neurological conditions.

The **Klip** can also be used in the anterior position for all sizes other than the XS size.

It is simple and safe to use, easy to fold and can be used as a basic posterior gait trainer or as a fully supported gait trainer.

SPECIFICATIONS		Width inside Handles	Handle Height	Seat to Floor Height	Outside Dimension (Open at Base)	Product Weight	Weight Capacity	
	2 wheeled							
	KP210	13.5"	16"-19.5"		23" x 23"	7 lbs.	75 lbs.	
	KP220R	14"	21"-27"		24" x 31"	9 lbs.	85 lbs.	
	KP230	16"	24"-32"		23" x 35"	11 lbs.	150 lbs.	
	4 wheeled							
	KP410	13.5"	16"-19.5"		23" x 23"	9 lbs.	75 lbs.	
	KP420R	14"	21"-27"		24" x 31"	11 lbs.	85 lbs.	
	KP430	16"	24"-32"		23" x 35"	13 lbs.	150 lbs.	
	KP440B	17"	30"-37"		26" x 37"	13 lbs.	190 lbs.	
4 wheeled with Seat								
KP510	13.5"	16"-19.5"	9.50" - 13.50"	23" x 23"	11 lbs.	75 lbs.		
KP520R	14"	21"-27"	12" - 18"	24" x 31"	13 lbs.	85 lbs.		
KP530	16"	24"-32"	14" - 21"	23" x 35"	17 lbs.	150 lbs.		
KP540B	17"	30"-37"	17" - 24"	26" x 37"	18 lbs.	190 lbs.		
4 wheeled with 8" wheels								
KP418	13.5"	18"-21.5"		23" x 25"	12 lbs.	75 lbs.		
KP428R	14"	23"-30"		24" x 33"	15 lbs.	85 lbs.		
KP438	16"	28"X35"		25" x 37"	17 lbs.	150 lbs.		
KP448B	17"	33"X41"		26" x 40"	18 lbs.	190 lbs.		
4 wheeled with seat and 8" wheels								
KP518	13.5"	18"-21.5"	12.5"-16.5"	23" x 25"	16 lbs.	75 lbs.		
KP528R	14"	23"-30"	15" - 21"	24" x 33"	17 lbs.	85 lbs.		
KP538	16"	28"X35"	17" - 24"	25" x 37"	21 lbs.	150 lbs.		
KP548B	17"	33"X41"	20" X 28"	26" x 40"	21 lbs.	190 lbs.		

Height-adjustable to ensure proper fit and accommodate future growth.
Foot-Activated swivel wheel lock with audible “click” enables easier and greater control.* The On/Off one-directional rear wheels prevents walker from rolling backwards.
Simple to disengage as walking skills improve.**
Folds easily for compact storage and easy transport



2 Wheeled

The 2-Wheeled Klip is specifically designed for users with significant balance impairments. It's two rear crutch tips provides extra resistance and is slower moving than a 4 wheeled walker.

It is available in three sizes for toddlers, kids and young adults.



4 Wheeled

The 4-Wheeled Klip provides maximum maneuverability with front swivel wheels that makes it easy to turn.

It is available in four sizes to accommodate toddlers, children, teens and adults.



Outdoors

The 4-Wheeled Klip with 8" swivel wheels is excellent for those with an active lifestyle. The wheels are specially designed to allow easy navigation on all indoor and outdoor surfaces.

It is available in four sizes to accommodate toddlers, children, teens and adults.



Flip-up seat

The 4-Wheeled Klip with a flip-up seat is a great convenience for users who may lack the endurance to walk for long periods of time. The seat flips up for walking and standing and folds down for sitting and resting.

It is available in four sizes to accommodate toddlers, children, teens and adults.



Outdoors/Seat

Simple to use and maneuver the 8" wheels on this Klip allow for smoother rolling on both indoor and outdoor surfaces and is excellent on many fronts. The seat flips up for walking and standing and folds down for sitting and resting.

It is available in four sizes to accommodate toddlers, children, teens and adults.

**Not available on the 2 wheeled and outdoor models*

***Not available on the 2 wheeled models*

ACCESSORIES

The **Klip** is a basic gait trainer with a choice of accessories designed to improve support, posture and mobility. It offers individualized assistive support.

The accessories can be easily added or removed as walking skills and needs evolve.



FOREARM PLATFORM WITH HANDGRIP

Fully adjustable forearm platforms are height and depth adjustable. The handgrips extend forward or back, supinate, pronate, rotate and are wrist-angle adjustable.



Extends



Rotates



Wrist
Angles



Wrist
Rotation

Width and height adjustable pelvic stabilizer with a support belt.



PELVIC STABILIZER



SLING SEAT

Adjustable cushioned sling seat provides partial weight bearing and properly positions user's pelvis.