

new Meywalk Small



- Designed for the junior sized child, too large for a Miniwalk but too small for the std Med/Large
- Moves up and down during gait, allowing for a more normal heel strike, toe off and mid phase as the child moves forward
- Spring assist system to aide the child to "swing through" during gait
- Standard with castor swivel locks. Can be engaged for in line mobility
- Adjustable angle handlebar encourages more upright or forward trunk orientation
- Lever lift transfer system like the Medium and Large Meywalks

Dynamic Gait Movement

The Meywalk moves up and down with each step, allowing the user to come down into a heel strike and toe off. As the user steps forward, the spring lifts and assists with the swing through phase of gait

Wide range of accessories.....

- Anti tippers (req'd for no back and drag brakes)
- No back (reverse) brakes
- Drag Brakes
- Midline Leg Separation Panel
- Lateral Hip Pads
- Castor Swivel Locks (standard)

How to measure

Child's axilla height

Child's seat height



Sizing	Sml/Jr
Axilla Height	31-38.5"
Seat Height	17.5-25"
Max users weight	150 lb

Finally... a Gait Trainer with Dynamic Movement to promote normal gait motion!

Call / Fax / Order Online

Southwestmedical.com
 ph: 800.236.4215
 fx: 602.279.0952