


# Catalyst 4 Tires - Pneumatic



Item Number	Part Number	Description	Quantity	Retail
1a-3	110001	Tire 20" x 1 3/8" Pneumatic, Inner Tube and Rim Strip Kit	1	\$44.12
1b-3	110002	Tire 22" x 1 3/8" Pneumatic, Inner Tube and Rim Strip Kit	1	\$44.12
1c-3	110003	Tire 24" x 1 3/8" Pneumatic, Inner Tube and Rim Strip Kit	1	\$45.22
1a	100911	Tire 20" x 1 3/8" Pneumatic	1	\$26.47
1b	100634	Tire 22" x 1 3/8" Pneumatic	1	\$26.47
1c	111818	Tire 24" x 1 3/8" Pneumatic (K103) <i>Replacement for part number 100635</i>	1	\$26.27

## Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part Number	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke <sup>2</sup>	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke <sup>2</sup>	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke <sup>2</sup>	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke <sup>2</sup>	200532	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke <sup>2</sup>	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke <sup>2</sup>	200534	100636			100636 / 100505	102978	100884	102960		100865			
24" Superlight	101159	111818					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
12" Mag (non-drum brake) <sup>1</sup>	N/A	001424	001426		001425								
16" Mag (non-drum brake) <sup>1</sup>	N/A	001159	001161	003453	001160								
20" 5-Spoke X Core (Black)	101961	100911	100917		100911 / 100913		102717						
22" 5-Spoke X Core (Black)	101962	100634	100621		100634 / 100503		102718						
24" 5-Spoke X Core (Black)	100960	111818	100622		111818 / 100504		100640		102982			102920	
20" Spinergy Spox	See Spinergy Spox Page	100911	100917	102085	100911 / 100913	102974		102957					
22" Spinergy Spox		100634	100621	100632	100634 / 100503	103165		102961					101475
24" Spinergy Spox		111818	100622	100633	111818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy Spox						102977	101146	102959		101357	103177		100166
26" Spinergy Spox		100636				100636 / 100505	102978	100895	102960		100865		
22" Spinergy LX	See Spinergy LX Page	100634	100621	100632	100634 / 100503			102961					101475
24" Spinergy LX		111818	100622	100633	111818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinergy LX						102977	101146	102959		101357	103177		100166
26" Spinergy LX		100636				100636 / 100505	102978	100895	102960		100865		
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly

<sup>2</sup> Discontinued