



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.2 Rev A

Client _____ Age _____ Height _____ Weight _____ Dx _____
Facility _____ Therapist _____ Ph. _____
Dealer _____ PO/RFQ: _____
Contact _____ Ph. _____ Date _____
Email _____

Superstand HLT is available in three sizes: HLT EI, HLT 1, and HLT 2

Standard Features Include: Horizontal high-loading surface, Prone, Vertical, Supine or Multi-Position configurable; hip and chest pads: hands-free standing angle adjustments; height and depth adjustable knee system; lock out block and 3" locking casters. Section 1 of this Order Form provides for selection of individual components to create a standing system to suit your client's requirements.

Section 2 of this form offers *All-In-One* Pre-Configured Units that supply all the components needed for a complete standing program. Packages include: Hip pad with laterals, chest pad with laterals, support vest, soft hip pad, headrest, chinrest, multi-position foot system, opaque tray, hands-free standing angle adjustments, height & depth adjustable knee pads, lock out block and upgraded steel locking casters. Pre-configured Units can be accessorized to fit your client's needs with listed items in Section 2 of the order form.

Specifications	HLT EI	HLT 1	HLT 2
Weight Capacity	65 Lbs.	125 Lbs.	150 Lbs.
Base	W 19" x L 28"	W 21" x L 38"	W25 " x L 45"
Loading Height	27"	28"	29"
Client Height: Supine	20" - 38"	28" - 48"	40" - 60"
Client Height: Prone	22" - 40"	30" - 50"	40" - 62"



Do Not exceed weight capacity.
Unit shown with optional equipment.
All specifications subject to change without notice.

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.6 Rev A

Order Form Section 1: Configurable Unit Section		HLT Size EI	HLT Size 1	HLT Size 2	
Superstand HLT is configurable with the options in the below columns. Please include any custom information at the bottom of the Order Form.					
Frame Colors: ☐ Candy Red ☐ Candy Blue ☐ Candy Purple ☐ Neon Pink ☐ Neon Orange ☐ Shocker Green ☐ Razzleberry ☐ White ☐ Black					
Upholstery Colors: ☐ Black Vinyl ☐ Blue Vinyl ☐ Red Vinyl ☐ Custom: _____					
Folding Base: Includes quick release folding base, mast mounting hardware, gas cylinder and 4 plastic locking casters.	☐ HT00010 \$995.00 ☐ Supine ☐ Prone ☐ Multi-Position	☐ HT00020 \$1,295.00 ☐ Supine ☐ Prone ☐ Multi-Position	☐ HT00030 \$1,495.00 ☐ Supine ☐ Prone ☐ Multi-Position		
Mast: Includes chest and hip support pad, knee supports, and solid foot platform. EI mast fits all three sized bases, Size 1 mast fits Size 1 and Size 2 Bases, Size 2 mast only fits Size 2 base.	☐ HT00040 EI Mast 20" - 38" Client Ht.	☐ HT00040 EI Mast 20" - 38" Client Ht. ☐ HT00050 Sz1 Mast 28" - 40" Client Ht.	☐ HT00040 EI Mast 20" - 38" Client Ht. ☐ HT00050 Sz1 Mast 28" - 40" Client Ht. ☐ HT00060 Sz2 Mast 40" - 60" Client Ht.	\$1,695.00 \$1,895.00 \$2,295.00	
Lateral Positioning Supports					
Lateral Chest Supports Please choose mounting bracket configuration. See width adjustment ranges below part number. Chest Pad Dimensions Chest Lateral Pad Dimensions	☐ HT00070 90 Deg. 6" - 10" W adj. ☐ HT00080 1" offset 4" - 8" W adj. ☐ HT00090 2" offset 3" - 6" W adj. 8 1/2"W x 7"H 3 1/2"D x 3"H	☐ HT00100 90 Deg. 8" - 12" W adj. ☐ HT00110 1" offset 6" - 10" W adj. ☐ HT00120 2" offset 4" - 8" W adj. 10" W x 8 1/2" H 4 1/2" D x 4" H	☐ HT00130 90 Deg. 10" - 14" W adj. ☐ HT00140 1" offset 9" - 12" W adj. ☐ HT00150 2" offset 6" - 11" W adj. 12 1/2"W x 15"H 4 1/2" D x 6"H	\$285.00	
Sliding Adjustable Tracks for Chest Pads *sliding tracks supply maximum height and angle adjustment.	Not available on EI	☐ HT00250 Sz1: non-contoured pad only	☐ HT00250 Sz2: std pad only	\$185.00	
Lateral Hip Supports Please choose mounting bracket configuration. See width adjustment ranges below part number. Hip Pad Dimensions Hip Lateral Pad Dimensions	☐ HT00160 90 Deg. 6" - 10" W adj. ☐ HT00170 1" offset 4" - 8" W adj. ☐ HT00180 2" offset 3" - 6" W adj. 8 1/2"W x "H 3 1/2"D x 3"H	☐ HT00190 90 Deg. 8" - 12" W adj. ☐ HT00200 1" offset 6" - 10" W adj. ☐ HT00210 2" offset 4" - 8" W adj. 10"W x 8"H 4 1/2"D x 5 1/2"H	☐ HT00220 90 Deg. 10" - 14" W adj. ☐ HT00230 1" offset 9" - 12" W adj. ☐ HT00240 2" offset 6" - 11" W adj. 12 1/2"W x 9"H 4 1/2"D x 6 1/2"H	\$285.00	
Sliding Adjustable Tracks for Hip Pads Size 1: non-contoured pad Size 2: standard pad	N/A	☐ HT00260	☐ HT00260	\$185.00	
Straps and Supports Harness Color: ☐ Blk Denim ☐ Blu Denim ☐ Blk Vinyl ☐ Blu Vinyl ☐ Red Vinyl					
Chest Belt black nylon velcro	☐ HT00270	☐ HT00280	☐ HT00290	\$95.00	
Padded Chest Belt black nylon w/pad	☐ HT00300	☐ HT00310	☐ HT00320	\$155.00	
Body Support Vest Quick Release Vest	☐ HT00330	☐ HT00340	☐ HT00350	\$255.00	
Soft Hip Support Soft hip strap for Prone	☐ HT00360	☐ HT00370	☐ HT00380	\$229.00	
Single patient use body support vest, washable	☐ HT00390	N/A	N/A	\$159.00	
Single patient use prone hip support, washable	☐ HT00400	N/A	N/A	\$159.00	
Steel Caster Upgrade High performance easy-rolling steel casters that lock for swivel and roll.	☐ HT00810	☐ HT00810	☐ HT000810	\$259.00	



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.6 Rev A

	HLT EI	HLT Sz 1	HLT Sz 2	
Headrest Height and depth adjustable, padded. (Incl. 1/2" Square Tube Mount)	☐ HT00410 6 1/2"W x 3 1/2"H	☐ HT00410 6 1/2"W x 3 1/2"H	☐ HT00415 8"W x 4 1/2"H	\$325.00
1/2" Square Tube Mount Height and depth adjustable, for any headrest using 1/2" mounting bar. (mount only)	☐ HT00420	☐ HT00420	☐ HT00420	\$159.00
i2i Flip-Back Mount Flip-back mounting hardware for i2i Headrest. (mount only)	☐ HT00430	☐ HT00430	☐ HT00430	\$429.00
1/2" Square Tube Flip-Back Mount Flip-back mount for any headrest using 1/2" mounting bar. (mount only)	☐ HT00440	☐ HT00440	☐ HT00440	\$429.00
Chinrest Height adjustable.	☐ HT00450	☐ HT00450	☐ HT00450	\$155.00
Padded Head Strap Sheepskin-lined.	☐ HT00460	☐ HT00460	☐ HT00460	\$139.00
Lower Extremity Positioning				
1-Piece Foot Platform (Default standard)	☐ HT00470	☐ HT00470	☐ HT00470	N/C
Independent Footplates Individually height and depth adjustable with heel loops.	☐ HT00475	☐ HT00475	☐ HT00475	\$275.00
Multi-Position Foot System Individually adj. for plantar & dorsi flexion, rotation, and height. w/ heel loops & straps.	☐ HT00480	☐ HT00480	☐ HT00480	\$525.00
Size-Adjustable Knee Pads Width adjustment range: 2 1/2" x 4 1/2"	☐ HT00485 (pr) 2"W x 3"H	☐ HT00490 (pr) 4"W x 3 1/2"H	☐ HT00500 (pr) 4 1/4"H x 5"H	\$449.00
Contracture Knee System Supports contracted lower extremities with a fully adjustable system.	☐ HT00510	☐ HT00510	☐ HT00510	\$799.00
D-Ring Foot Straps	☐ HT00520	☐ HT00520	☐ HT00520	\$95.00
Shoe Holders Formed ABS. w/ankle & toe straps.	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	\$279.00
Height Adjustable Pommel	☐ HT00560	☐ HT00570	☐ HT00580	\$289.00
Height Adj. Pommel w/ Rigid Hip Pad Height and Depth adj. rigid hip pad.	☐ HT00590	☐ HT00600	☐ HT00600	\$549.00
Abduction Leg System W/Pommel Incl. rotating knees, fully adjustable shoe holders, and pommel with rigid pad	☐ HT00620	☐ HT00620	☐ HT00620	\$695.00

Trays				
Extremity Positioning Tray Multi-positional tray w/ height, depth & angle adj.	☐ HT00630 Opaque	☐ HT00640 Opaque	☐ HT00650 Opaque	\$599.00
	☐ HT00660 Clear	☐ HT00670 Clear	☐ HT00680 Clear	\$725.00
Tray Overlays	☐ HT00690 Opaque	☐ HT00700 Opaque	☐ HT00710 Opaque	\$239.00
	☐ HT00720 Clear	☐ HT00730 Clear	☐ HT00740 Clear	\$239.00
Tray Easel	☐ HT00750	☐ HT00750	☐ HT00750	\$289.00
Tray Edging	☐ HT00760	☐ HT00760	☐ HT00760	\$179.00
Activity Bowl Tray cut-out w 7" bowl	☐ HT00770	☐ HT00770	☐ HT00770	\$259.00
Padded Tray Covers Full Size are removable. Wrap Around are non-removable.	☐ HT00780 Full Size	☐ HT00780 Full Sized	☐ HT00780 Full Sized	\$275.00
	☐ HT00790 Wrap-arnd	☐ HT00790 Wrap-arnd	☐ HT00790 Wrap-arnd	\$279.00



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.6 Rev A



Sliding Adj. Tracks w/non-contoured pad



Multi-Adj. Foot System



Height adj. Foot System



Abduction Foot System



Size Adj. Knee Pads



Contracture Knee System



Tray Easel



Tray with Activity Bowl

Client Measurements A: Thorax Width _____" B: Hip Width _____" C: Hip to Knee _____" D: Knee to Heel _____"

E: Shoe Size _____" F: Shoulder Height _____" G: Height _____" H: Inseam _____" I: Weight _____ Lbs.


Note: All measurements assume symmetrical standing and/or sitting position.

Custom Pads All Custom sized pads must be quoted for price and manufacturing time. Describe custom shapes at the bottom.


Chest Pad: _____"W x _____"H Standard shaped pad Custom Shape Chest Laterals: _____"D x _____"H

Hip Pad: _____"W x _____"H Standard shaped pad Custom shape Hip Laterals: _____"D x _____"H

Custom Tray Surface _____"W x _____"D Standard Tray: Size 1: 22 1/2"W x 19"D Size 2: 28"W x 24"D



Include any information on
tray cutouts or submit a
paper template. Addition-
al information may be
added below



Special Requests / Additional Information

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/15/2025
Doc.# FM-10.6 Rev A

Order Form Section 2: All-in-One Pre-Configured Units

Pre-Configured Units include: High-loading platform, Folding Base and components to support prone, vertical, supine, or multi-position configurations; hip pad with laterals, chest pad with laterals, support vest, soft hip pad, headrest, chinrest, multi-position foot system, opaque tray, height & depth adjustable knee pads, and upgraded steel casters. Units can be accessorized with items listed in this section ONLY. If additional configurations or customizations are necessary, please use Section 1 for your order.

All-in-One Pre-Configured Units	☐ HLT EI HLT_EI-ALLINONE \$5,195.00	☐ HLT 1 HLTALLINONE \$5,595.00	☐ HLT 2 HLTALLINONE_2 \$6,195.00	
--	---	--	--	--

Frame Colors: ☐ Candy Red ☐ Candy Blue ☐ Candy Purple ☐ Neon Pink ☐ Neon Orange ☐ Shocker Green
☐ Razzleberry ☐ White ☐ Black

Upholstery Colors: ☐ Black Vinyl ☐ Blue Vinyl ☐ Red Vinyl ☐ Custom: _____

Harness: ☐ Black Denim ☐ Blue Denim ☐ Black Vinyl ☐ Blue Vinyl ☐ Red Vinyl

Options

Lateral Chest Supports Please choose mounting bracket configuration. See width adjustment ranges below part number. Chest Pad dimensions: Chest Pad Chest Lateral Pads	☐ HT01070 90 Deg. 6" - 10" W adj. ☐ HT01080 1" offset 4" - 8" W adj. ☐ HT01090 2" offset 3" - 6" W adj. 8 1/2"W x 7"H 3 1/2"D x 3"H	☐ HT01100 90 Deg. 8" - 12" W adj. ☐ HT01110 1" offset 6" - 10" W adj. ☐ HT01120 2" offset 4" - 8" W adj. 10"W x 8 1/2"H 4 1/2"D x 4"H	☐ HT01130 90 Deg. 10" - 14" W adj. ☐ HT01140 1" offset 9" - 12" W adj. ☐ HT01150 2" offset 6" - 11" W adj. 12 1/2"W x 15"H 4 1/2"D x 6"H	Included
Sliding Adjustable Tracks for Chest Pads *sliding tracks supply maximum height and angle adjustment. Sz 1: non-contoured pad Sz 2: standard pad	Not available on EI	☐ HT00250 Sz 1: non-contoured pad only	☐ HT00250 Sz 2: std pad only	\$185.00
Lateral Hip Supports Please choose mounting bracket configuration. See width adjustment ranges below part number. Hip Pad dimensions: Hip Pad Hip Laterals	☐ HT01160 90 Deg. 6" - 10" W adj. ☐ HT01170 1" offset 4" - 8" W adj. ☐ HT01180 2" offset 3" - 6" W adj. 8 1/2"W x 7"H 3 1/2"D x 3"H	☐ HT01190 90 Deg. 8" - 12" W adj. ☐ HT01200 1" offset 6" - 10" W adj. ☐ HT01210 2" offset 4" - 8" W adj. 10"W x 8"H 4 1/2"D x 5 1/2"H	☐ HT01220 90 Deg. 10" - 14" W adj. ☐ HT01230 1" offset 9" - 12" W adj. ☐ HT01240 2" offset 6" - 11" W adj. 12 1/2"W x 9"H 4 1/2"D x 6 1/2"H	Included
Sliding Adjustable Tracks for Hip Pads	N/A	☐ HT00260	☐ HT00260	\$185.00
Chest Belt black nylon velcro	☐ HT00270	☐ HT00280	☐ HT00290	\$95.00
Padded Chest Belt black nylon w/pad	☐ HT00300	☐ HT00310	☐ HT00320	\$155.00

Head Support Options Packages include standard Headrest. These options omit the headrest.

1/2" Square Tube Mount Only. Height and depth adjustable, for any headrest using 1/2" Horizontal mounting bar.	☐ HT01420	☐ HT01420	☐ HT01420	\$0.00
i2i Flip-Back Mount Flip-back mounting hardware for i2i Headrest. (mount only)	☐ HT01430	☐ HT01430	☐ HT01430	\$345.00
1/2" Square Tube Flip-Back Mount Flip-back mount for any headrest using 1/2" horizontal mounting bar.	☐ HT01440	☐ HT01440	☐ HT01440	\$245.00

PRIME ENGINEERING
A division of Axiom Industries

4202 W. Sierra Madre Ave. Fresno, CA 93722
Ph (800)827-8263 (559)271-0991

Email: info@primeengineering.com
Website: www.primeengineering.com
Fax (800)800-3355 (559)276-3544



Superstand HLT

Pediatric Standing System

Order Form / Price List
Effective: 6/1/2025
Doc.# FM-10.6 Rev A

All-in-One Pre-Configured Units	HLT EI	HLT Sz1	HLT Sz2	
Size-Adjustable Knee Pads Width adjustment range: 2 1/2" x 4 1/2"	☐ HT00485 (pr) 2"W x 3"H	☐ HT00490 (pr) 4"W x 3 1/2"H	☐ HT00500 (pr) 4 1/4"H x 5"H	\$449.00
Contracture Knee System Supports contracted lower extremities with a fully adjustable system.	☐ HT00510	☐ HT00510	☐ HT00510	\$799.00
D-Ring Foot Straps	☐ HT00520	☐ HT00520	☐ HT00520	\$95.00
Shoe Holders Formed ABS. w/ankle & toe straps.	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	☐ HT00530 Sm. 6"L x 4 1/4"W ☐ HT00540 Med. 8"L x 5"W ☐ HT00550 Lg. 10 1/2"L x 6"W	\$279.00
Height Adjustable Pommel	☐ HT00560	☐ HT00570	☐ HT00580	\$289.00
Height Adj. Pommel w/ Rigid Hip Pad Height and Depth adj. rigid hip pad.	☐ HT00590	☐ HT00600	☐ HT00600	\$549.00
Abduction Leg System W/Pommel Incl. size adj knees, shoe holders, multi	☐ HT00620	☐ HT00620	☐ HT00620	\$695.00

Tray Options Opaque Tray included. Clear option replaces opaque tray top.				
Extremity Positioning Tray Multi-positional tray w/ height, depth & angle adj.	☐ HT01630 Opaque ☐ HT01660 Clear	☐ HT01640 Opaque ☐ HT01670 Clear	☐ HT01650 Opaque ☐ HT01680 Clear	Included \$149.00
Tray Overlays	☐ HT00690 Opaque ☐ HT00720 Clear	☐ HT00700 Opaque ☐ HT00730 Clear	☐ HT00710 Opaque ☐ HT00740 Clear	\$239.00 \$239.00
Tray Easel	☐ HT00750	☐ HT00750	☐ HT00750	\$289.00
Tray Edging	☐ HT00760	☐ HT00760	☐ HT00760	\$179.00
Activity Bowl Tray cut-out w 7" bowl	☐ HT00770	☐ HT00770	☐ HT00770	\$259.00
Padded Tray Covers	☐ HT00780 full size removable	☐ HT00780 full size removable	☐ HT00780 full size removable	\$275.00
	☐ HT00790 Wrap around non-removable.	☐ HT00790 Wrap around non-removable.	☐ HT00790 Wrap around non-removable.	\$279.00

Special Requests / Additional Information
