


CR45 24" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and handrim kits will be listed. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size, new axles/receivers and handrims may be needed. If changing wheel size and Attendant Foot Lock or T-Arms are present, spacer 112503 is required which can be found on the CR45 Axle Plate page See CR45 Axles page for axle and receiver utilization.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118359 | 24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118360 | 24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 118361 | 24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 118362 | 24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 118367 | 24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 118368 | 24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount | 1 |
| | 118369 | 24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118370 | 24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount | 1 |
| | 118136 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 118137 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 118138 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 118139 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 118144 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 118145 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount | 1 |
| | 118146 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118147 | 24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount | 1 |
| | 117696 | 24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 117697 | 24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 117698 | 24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 117699 | 24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 117704 | 24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 117705 | 24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount | 1 |
| | 117706 | 24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount | 1 |
| | 117707 | 24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount | 1 |