


BACK UPHOLSTERY MEASURING GUIDES

MEASURING GUIDE

QUICKIE 2/M6/Rx BACK UPHOLSTERY



DESCRIPTION	"A"	"B"	"C"
BACK UPH 12 LOW	12-3/4	7-1/2	14-1/4
BACK UPH 12 MED	12-3/4	10-1/2	17-1/4
BACK UPH 12 TALL	12-3/4	13-1/2	20-1/4
BACK UPH 13 LOW	13-3/4	7-1/2	14-1/4
BACK UPH 13 MED	13-3/4	10-1/2	17-1/4
BACK UPH 13 TALL	13-3/4	13-1/2	20-1/4
BACK UPH 14 LOW	14-3/4	7-1/2	14-1/4
BACK UPH 14 MED	14-3/4	10-1/2	17-1/4
BACK UPH 14 TALL	14-3/4	13-1/2	20-1/4
BACK UPH 15 LOW	15-3/4	7-1/2	14-1/4
BACK UPH 15 MED	15-3/4	10-1/2	17-1/4
BACK UPH 15 TALL	15-3/4	13-1/2	20-1/4
BACK UPH 16 LOW	16-3/4	7-1/2	14-1/4
BACK UPH 16 MED	16-3/4	10-1/2	17-1/4
BACK UPH 16 TALL	16-3/4	13-1/2	20-1/4
BACK UPH 17 LOW	17-3/4	7-1/2	14-1/4
BACK UPH 17 MED	17-3/4	10-1/2	17-1/4
BACK UPH 17 TALL	17-3/4	13-1/2	20-1/4
BACK UPH 18 LOW	18-3/4	7-1/2	14-1/4
BACK UPH 18 MED	18-3/4	10-1/2	17-1/4
BACK UPH 18 TALL	18-3/4	13-1/2	20-1/4
BACK UPH 19 LOW	19-3/4	7-1/2	14-1/4
BACK UPH 19 MED	19-3/4	10-1/2	17-1/4
BACK UPH 19 TALL	19-3/4	13-1/2	20-1/4
BACK UPH 20 LOW	20-3/4	7-1/2	14-1/4
BACK UPH 20 MED	20-3/4	10-1/2	17-1/4
BACK UPH 20 TALL	20-3/4	13-1/2	20-1/4
BACK UPH 21 LOW	21 3/4	7 1/2	14 1/4
BACK UPH 21 MED	21 3/4	10 1/2	17 1/4
BACK UPH 21 TALL	21 3/4	13 1/2	20 1/4
BACK UPH 22 LOW	22 3/4	7 1/2	14 1/4
BACK UPH 22 MED	22 3/4	10 1/2	17 1/4
BACK UPH 22 TALL	22 3/4	13 1/2	20 1/4
BACK UPH 24 TALL	24 3/4	13 1/2	20 1/4
BACK UPH 26 TALL	26 3/4	13 1/2	20 1/4
BACK UPH 28 TALL	28 3/4	13 1/2	20 1/4
BACK UPH 30 TALL	30 3/4	13 1/2	20 1/4