


Arc Seat Frame - Growth



Is your depth changing? Will your existing seat frame achieve your new depth? If no, please select the appropriate seat frames below. Select color of parts based on Seat Frame color of chair.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1a, 4b, 6-8 | 113121_ _ _ | Seat Frame Assembly, Standard, X-Short, Arc <i>Accommodates 12-15" Depths</i> | 1 |
| 1b, 4a, 6-8 | 113122_ _ _ | Seat Frame Assembly, Standard, Short, Arc <i>Accommodates 14-17" Depths</i> | 1 |
| 1c, 4a, 6-8 | 113123_ _ _ | Seat Frame Assembly, Standard, Medium, Arc | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | | <i>Accommodates 16-19" Depths</i> | |
| 1d, 4a, 6-8 | 113124_ _ _ | Seat Frame Assembly, Standard, Long, Arc <i>Accommodates 18-20" Depths</i> | 1 |
| 2a, 4b, 6-9 | 113125_ _ _ | Seat Frame Assembly, Offset, X-Short, Right, Arc <i>Accommodates 12-15" Depths</i> | 1 |
| 2e, 4b, 6-9 | 113126_ _ _ | Seat Frame Assembly, Offset, X-Short, Left, Arc <i>Accommodates 12-15" Depths</i> | 1 |
| 2b, 4a, 6-9 | 113127_ _ _ | Seat Frame Assembly, Offset, Short, Right, Arc <i>Accommodates 14-17" Depths</i> | 1 |
| 2f, 4a, 6-9 | 113128_ _ _ | Seat Frame Assembly, Offset, Short, Left, Arc <i>Accommodates 14-17" Depths</i> | 1 |
| 2c, 4a, 6-9 | 113129_ _ _ | Seat Frame Assembly, Offset, Medium, Right, Arc <i>Accommodates 16-19" Depths</i> | 1 |
| 2g, 4a, 6-9 | 113130_ _ _ | Seat Frame Assembly, Offset, Medium, Left, Arc <i>Accommodates 16-19" Depths</i> | 1 |
| 2d, 4a, 6-9 | 113131_ _ _ | Seat Frame Assembly, Offset, Long, Right, Arc <i>Accommodates 18-20" Depths</i> | 1 |
| 2h, 4a, 6-9 | 113132_ _ _ | Seat Frame Assembly, Offset, Long, Left, Arc <i>Accommodates 18-20" Depths</i> | 1 |
| 3a, 4b, 9 | 113133_ _ _ | Seat Frame Assembly, Contracture, X-Short, Arc <i>Accommodates 12-15" Depths</i> | 1 |
| 3b, 4a, 9 | 113134_ _ _ | Seat Frame Assembly, Contracture, Short, Arc <i>Accommodates 14-17" Depths</i> | 1 |
| 3c, 4a, 9 | 113135_ _ _ | Seat Frame Assembly, Contracture, Medium, Arc <i>Accommodates 16-19" Depths</i> | 1 |
| 3d, 4a, 9 | 113136_ _ _ | Seat Frame Assembly, Contracture, Long, Arc <i>Accommodates 18-20" Depths</i> | 1 |
| 4a | 110165 | Growth Tube, Seat <i>Used on short, medium, and long seat depth ranges</i> | 1 |
| 4b | 110166 | Growth Tube, Seat, Extra Short <i>Used on x-short seat depth range</i> | 1 |
| 5a | 110206 | Mechlock, 194mm Stroke, Right <i>Used on short, medium, and long seat depth ranges</i> | 1 |
| 5b | 110207 | Mechlock, 194mm Stroke, Left <i>Used on short, medium, and long seat depth ranges</i> | 1 |
| 5c | 110208 | Mechlock, 190mm Stroke, Right | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | | <i>Used on x-short seat depth range</i> | |
| 5d | 110209 | Mechlock, 190mm Stroke, Left <i>Used on x-short seat depth range</i> | 1 |