


# Discontinued 1 Rogue XP Tension Adjustable Back Upholstery 17" - 18" Wide



| Back Height | Strap Qty<br>(Item 6) | D-Ring Qty<br>(Item 7) |
|-------------|-----------------------|------------------------|
| 9"-12"      | 3                     | 6                      |
| 11"-14"     | 3                     | 6                      |
| 13"-16"     | 4                     | 8                      |
| 15"-18"     | 5                     | 10                     |

*This Tension Adjustable Back Upholstery was effective prior to November 2016. <p>*

*Please see Rogue XP Tension Adjustable Back Upholstery 17" - 18" Wide page for current assemblies.*

| <b>Item Number</b> | <b>Part Number</b> | <b>Description</b>                      | <b>Quantity</b> | <b>Retail</b> |
|--------------------|--------------------|---|-----------------|---------------|
| 7                  | 100973             | Back Upholstery D-Ring                  | 6               | \$10.51       |
| 8                  | 100734             | M5 x 16 PHMS ZC                         | 2               | \$1.58        |
| 9a                 | 001751             | Tension Strap Folding Handle 16-17 Wide | 1               | \$21.98       |
| 9b                 | 001752             | Tension Strap Folding Handle 18-19 Wide | 1               | \$23.28       |